

QUICK REFERENCE GUIDE



SeQual[®] Eclipse 5[®]

Transportable Oxygen Concentrator with autoSAT[®]

Battery Life

Continuous Flow		Pulse Flow		
0.5 LPM – 4.4 Hours	2.0 LPM – 2.0 Hours	1.0 (16 mL) – 5.4 Hours	4.0 (64 mL) – 4.0 Hours	7.0 (128 mL) – 2.5 Hours
1.0 LPM – 3.7 Hours	3.0 LPM – 1.3 Hours	2.0 (32 mL) – 5.1 Hours	5.0 (80 mL) – 3.7 Hours	8.0 (160 mL) – 2.0 Hours
		3.0 (48 mL) – 4.9 Hours	6.0 (96 mL) – 3.5 Hours	9.0 (192 mL) – 1.7 Hours

NOTE: Pulse flow times are calculated based on 12BPM. Higher breath rates will decrease battery times on pulse mode only.

Recharge Time

1. The typical time to recharge the battery inside the Eclipse on AC or DC power is between 2–5 hours, dependent upon the flow setting.
The battery may be re-charged while the Eclipse is either ON or OFF.
2. The battery will recharge using the optional desktop charger in approximately 3 hours.
3. It is acceptable to leave the Eclipse plugged into AC power at all times while it is OFF. Note that while the unit remains plugged in and OFF, the internal fan may continue to run while the battery is charging and produce a faint humming noise.

Battery Gauge

The Eclipse 5 has a battery status gauge on the right side of the display screen. Each of the 5 horizontal gray bars represents approximately 20% of the total battery charge. When a battery is being charged, the indicator bars will blink in a waterfall-type fashion. If a battery is not installed, or is improperly installed, the battery status gauge will not appear on the screen.

On/Off

Press and hold the ON/OFF button (top left Blue Button) for two seconds to power on the Eclipse 5. Upon power up, the audible alarm should sound for 3 seconds, and all three lights on above the display screen should illuminate. It is normal for the yellow and/or red lights to remain on for up to 5 minutes during warm-up. Once warm-up is complete, only the green light should be on. Press and hold the ON/OFF button for two seconds to power off the Eclipse 5.

Changing Flow Mode

Press the flow mode button (bottom left blue button) to toggle between pulse and continuous flow. The light adjacent to the button will illuminate when in pulse mode and will blink with each breath detection and delivery. If the unit is in pulse mode but does not detect a breath, it will switch to continuous mode.

AC Power Supply

1. The Eclipse 5 will run at all flow rates and recharge a battery while connected to AC power.
2. Each AC Power supply has a removable cord. If the AC power supply is plugged into and receiving power from a wall outlet, a green LED light on the power supply brick will illuminate.
3. If the AC Power Supply is properly connected to the Eclipse 5, the external power light on the far right side of the control panel will illuminate green.
4. The Eclipse AC Power supply is rated for 100-240VAC and 50-60Hz. When traveling, the power supply will automatically adjust to the specific voltage for each country. If needed, you may use outlet adapters that can be purchased from CAIRE or from a local electronics store.

DC Power Supply

1. Each Eclipse 5 comes with a DC power supply that will allow the unit to be powered from an automotive-style power outlet.
2. The Eclipse can operate at all flow settings while using DC power.
3. The Eclipse battery will recharge while in use and connected to DC power at 2.0 LPM or lower, and at all pulse flow settings. There is not sufficient power to run the Eclipse and recharge the battery on DC power at 2.5 LPM and 3.0 LPM settings.
4. If the DC power supply is plugged into and receiving power from a DC outlet, and the power supply is properly plugged into the side of the Eclipse, a green LED on the power supply brick will illuminate.
5. If the DC Power Supply is properly connected to the Eclipse 5, the external power light on the far right side of the control panel will illuminate green. The Eclipse DC power supply is rated between 11.5-16 VDC.
6. It is suggested to always have a motor vehicle cranked and running before inserting the DC power supply into the accessory outlet.

Maintenance

The Air Intake filter must be removed and cleaned by the patient at least once weekly. The air intake filter is located on the upper rear of the Eclipse 5 and can be removed and cleaned with soap and water. Allow to dry thoroughly before re-installing the filter.

Nasal Cannula and Tubing

The Eclipse 5 uses a standard nasal cannula. Tubing length should not exceed 50 ft (15.2 m) on a continuous flow setting, or 7 ft (2.1 m) on a pulse flow setting or when using a humidifier bottle on continuous flow.

Humidifier Bottle

A humidifier bottle may be used with the Eclipse on CONTINUOUS FLOW SETTINGS ONLY, and ONLY WHEN THE UNIT IS STATIONARY. If the unit is to be transported or moved, the humidifier bottle MUST be removed.



This Quick Reference Guide is not intended to replace the detailed information available in the User Manual, so please refer to the User Manual for important critical information regarding use and operation of the product.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.