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<i>ENGLISH</i>	<i>PAGE</i>
LIQUID OXYGEN SYSTEMS	2-27

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Important Notice

Be sure to read and understand this instruction manual completely before attempting to operate your Companion® Liquid Oxygen System. Anyone who has not read these instructions must not be permitted to handle or operate this equipment.

The Companion Liquid Oxygen Systems are intended only for the delivery of medical grade oxygen as prescribed by your physician. **Oxygen supplied from this equipment is for supplemental use and is not intended to be life supporting or life-sustaining.**

This manual is provided to assist in the safe operation of the Companion Liquid Oxygen Systems and to ensure maximum benefit from its use. If you have any questions about the operation of the equipment, please consult your liquid oxygen supplier.

Table of Contents

General Safety Precautions	Page 3
Companion Liquid Oxygen Systems.....	Page 8
Companion Stationary Units	Page 10
<i>Stationary Unit Components</i>	Page 10
<i>Using the Stationary Unit</i>	Page 13
Companion Portable Units	Page 15
<i>Companion 1000/1000T Unit Components</i>	Page 15
<i>Companion 500/550 Unit Components</i>	Page 16
<i>Filling the Portable Unit</i>	Page 18
<i>Using the Portable Unit</i>	Page 26
Cleaning and Maintenance	Page 27

GENERAL SAFETY PRECAUTIONS

Throughout this manual you will see three kinds of information with emphasis in the text. Be sure to carefully read and understand these notices. Each is important and related to the adjacent text.

WARNING

Warnings are accompanied by symbols and a box which identifies a hazard. They tell of conditions that concern your personal safety and the safety of others. They are followed by the actions required to prevent injury. Ignoring these warnings can lead to injury or death.

CAUTION

Cautions inform you about possible damage to the equipment or other property. Ignoring cautions may cause damage or make the equipment inoperative.

NOTE: Notes are italicized and provide important information about using the equipment properly.

WARNING

Refer to documentation for information.

Do not smoke near this equipment – Keep cigarettes or burning tobacco away from area where equipment is operated.

Keep flammable materials away from this equipment – Oils and grease, including facial creams and petroleum jelly, ignite easily and may burn rapidly in the presence of oxygen. Never lubricate any part of this equipment.

Do not touch liquid oxygen or parts that have been exposed to liquid oxygen – Liquid oxygen is extremely cold (-297°F / -183°C). When touched, liquid oxygen or parts of the equipment that have been carrying liquid oxygen can freeze skin and body tissue.

Non flammable non-toxic gas - Oxygen is a non flammable non-toxic gas

Fire intensifying risk - High concentrations of oxygen can cause rapid burning of other substances.

Keep and use this equipment in upright position at all times – If Stationary or Portable unit is turned over, gaseous or liquid oxygen will escape. Should a liquid spill occur, ventilate area by opening doors and windows. Call your liquid oxygen supplier immediately.

MFR: Manufacturer.

Name and address of manufacturer

Authorized representative in the European Community

WARNING

Keep this equipment away from electrical appliances – Use and store Stationary unit and Portable unit at least 1.5 meters from electrical appliances that may cause heat or sparks.

Keep oxygen equipment away from open flame – Heating sources such as furnaces, water heaters, and stoves may contain open flames.

Keep equipment in well ventilated area at all times – These units periodically release small amounts of oxygen gas that must be adequately ventilated. Do not store liquid oxygen equipment in a closet, car trunk or other confined area. Do not place clothing or draperies over equipment.

Do not place the Portable unit under clothing- These units normally vent oxygen, wearing a Portable unit under clothing may saturate fabrics with oxygen and cause them to burn rapidly if ignited.

SN

Serial Number.

Fragile; Handle with care.

Keep Dry.


This device complies with the requirements of Directive 93/42/EEC concerning medical devices. It therefore bears the CE marking as shown.

This pressure vessel complies with the requirements of Directive 2010/35/EU concerning transportable pressure equipment. It therefore bears the PI marking as shown.

①	Working Pressure	SN	Serial Number
②	Test Pressure	⚠	Refer to Manual
③	Tare Mass		
④	Water Capacity		
⑤	Inspection Mark		

①	Working Pressure	SN	Serial Number
②	Test Pressure	⚠	Refer to Manual
③	Tare Mass		
④	Water Capacity		
⑤	Inspection Mark		

COMPANION LIQUID OXYGEN SYSTEMS

The Companion Liquid Oxygen System is designed to store and deliver oxygen to the patient at a prescribed rate. The system comprises of two units; the Stationary and the Portable.

Stationary Units: Each Stationary unit is a liquid oxygen reservoir that may be used to fill your Portable unit. There are two basic types of Stationary unit - a standard unit that may be used to supply breathing oxygen in the home, (including **Companion 21, Companion 31 & Companion 41**) and a *Low Loss* unit (**Companion 31LL**) that allows liquid oxygen to be held for much longer periods of time than a typical Stationary unit but **does not** provide a source of oxygen for breathing purposes.

Each Stationary unit must be periodically refilled by a liquid oxygen supplier. The frequency of refill will depend upon the rate of use.


Figure 1: Companion Stationary Units

Portable Units: The Companion Portable oxygen unit is a compact, lightweight device that provides a portable supply of oxygen for the user both in and out of the home. It is filled from the Companion Stationary unit.

Portable unit variations include - **Companion 1000** (one litre model with a maximum flow rate of 6 L/min), **Companion 1000T** (one litre model with a maximum flow rate of 15 L/min), **Companion 500** (half litre model) and **Companion 550** (half litre model with a pneumatic demand valve that increases the use time by allowing oxygen to flow only during inspiration).

(All models may not be available in all markets).


Figure 2: Companion Portable Units

COMPANION STATIONARY UNITS

Companion Stationary Unit Components


Become familiar with the various parts of the unit described below and shown in Figure 3.

- *Contents Indicator* - displays the amount of liquid oxygen remaining in the Stationary unit. It is visible through the top of the unit.
- *Pressure Gauge* - (Optional) indicates the status of the pressure inside the Stationary unit. The Stationary unit is at an acceptable pressure when the needle is in the dark blue shaded region.
- *Flow Control* - (except *Low Loss*) adjusts oxygen flow from the breathing oxygen supply. It must be set to the oxygen flow rate prescribed by your physician.

NOTE: *Set control to 0 to turn unit off.*

- *Breathing Oxygen Supply* - (except *Low Loss*) this connector is the point where a cannula or other type of breathing device is connected for direct breathing from the Stationary unit. A bubble-type humidifier may be attached to this connection before attaching the supply tubing.

- 1 — **Contents Indicator**
- 2 — **Pressure Gauge - (Optional)** *Low Loss* has Pressure Gauge located where the Flow Control is.
- 3 — **Flow Control** (except *Low Loss*)
- 4 — **Breathing Oxygen Supply** (except *Low Loss*)

Figure 3: Stationary Unit Components

See Figure 4 for the following items:

- *Fill Connector* - is in the center of the top of the Stationary unit. When filling the Portable unit, it connects to the coupling in the bottom of the Portable unit.
- *Fill Connector Cover (Optional)* - protects the top fill connector when not in use. Always replace the fill connector cover after filling your Portable unit.

WARNING

Extreme Cold Hazard. Liquid oxygen can be released from fill connector. Do not depress or disturb fill connector poppet.

- *Portable Release Button* - is used to remove bottom-fill type Portable units from a Stationary unit when a fill is completed.


Figure 4: Stationary Unit Components

USING THE STATIONARY UNIT

When the Stationary unit is delivered to you by your liquid oxygen supplier, it is filled with oxygen. The Stationary unit will be refilled as required by your individual usage.

The Stationary unit may be used to fill Portable oxygen units. For filling procedures, see the '*Filling the Portable Unit*' section of this manual.

The Companion Stationary can be used directly as a source of breathing oxygen (except Companion *Low Loss* units). For this application, a humidifier (optional) and a nasal cannula or other breathing device is directly connected to the Stationary unit's oxygen supply outlet as follows:

1. If a *Humidifier* is used, fill it with distilled water to the level indicated on the humidifier's instructions. Attach the *Humidifier* to the *Breathing Oxygen Supply* fitting on the Stationary unit. (Figure 5).


Figure 5: Preparing Stationary Unit for Breathing

2. Attach the *Oxygen Tubing* from the cannula to the *Breathing Oxygen Supply* connection or to the humidifier connection, as appropriate. Adjust the cannula to your face (Figure 6).


Figure 6: Attaching a Nasal Cannula

3. Turn the *Flow Control* on the top of the Stationary unit to the setting prescribed by your physician.

NOTE: *Use only marked settings - oxygen will not flow if the control is set between flow settings. Do not set the flow above the flow rate prescribed by your physician.*

If a humidifier is used, check for a steady stream of bubbles. This indicates that oxygen is flowing.

4. The Stationary unit's *Contents Indicator* will show how much liquid oxygen remains in the unit. This indicator should be checked periodically to ensure the supply is adequate.
5. During use, moisture will condense on the internal parts of the unit. A Moisture Collection Container is provided.

COMPANION PORTABLE UNITS

Companion 1000/1000T Portable Unit Components

Become familiar with the various parts of the unit described below and shown in Figure 7.

- *Vent Valve* - is used to fill the unit. You must lift the vent valve lever to fill the unit when it is coupled to a Stationary unit. The vent valve lever is returned to its original position to stop the fill.
- *Fill Connector* - for connecting to a Stationary unit, and is found on the bottom of the unit.


Figure 7: 1000/1000T Portable Unit

- *Contents Indicator* - is found on the top of the Portable unit. A pointer tells you how much oxygen remains in the unit by unsnapping the plastic connector and suspending the unit from the end of the carrying strap that is closest to the indicator.
- *Oxygen Outlet* - allows connection of the cannula that delivers oxygen for breathing.
- *Moisture Cup and Pad (1000T Only)* - helps contain condensation formed when liquid oxygen is warmed to a gaseous state.
- *Carrying Strap* — This is an adjustable strap that enables the user to wear the unit over the shoulder. The plastic connector in the

carrying strap may be engaged to create a short carrying handle. The strap is also used to measure the contents of the Companion 1000/1000T (see *Contents Indicator* on page 17).


Figure 8: 1000/1000T Portable Unit

COMPANION 500/550 UNIT COMPONENTS

Become familiar with the various parts of the unit described below and shown in Figure 9.

- *Vent Valve* — This lever operated valve enables the user to fill the unit. Upon connecting the Companion 500/550 onto the Companion Stationary unit, the user begins the filling sequence by opening the vent valve and terminates the fill by closing the valve.
- *Flow Control* — This adjustable rotary valve controls the rate of oxygen flow from the unit and should be set for the user's prescribed oxygen flowrate.
- *Oxygen Outlet Connector* — Oxygen flow from the unit is delivered through this connector. The oxygen supply tube is connected here.

- *Contents Indicator* — The amount of liquid oxygen contained in the unit is measured by this internal weight scale mechanism. Liquid oxygen contents is indicated by unsnapping the plastic connector and suspending the unit from the end of the carrying strap that is closest to the indicator.
- *Fill Connector* — This connector, when coupled to a mating connector on a Companion stationary unit, provides a means for transferring liquid oxygen from the stationary unit to the Companion 500/550.
- *Carrying Strap* — This is an adjustable strap that enables the user to wear the unit over the shoulder. The plastic connector in the carrying strap may be engaged to create a short carrying handle. The strap is also used to measure the contents of the Companion 500/550 (see *Contents Indicator*).

C550 Portable Unit Components ONLY:

- *Continuous / Demand Flow Switch* — This toggle switch sets the oxygen delivery mode of the unit for either continuous flow or flow only during the inspiratory part of each breath.
- *Inspiration Sensor Connector* — The user's inspiratory effort is sensed through this connector. Either tube of the dual-lumen cannula is connected here.
- *Dual-Lumen Cannula* — This is a special cannula with two connecting tubes for use with the Companion 550. One cannula tube is connected to the unit's oxygen outlet connector and routes oxygen flow to the user. The second tube is connected to the unit's inspiration sensor connector and routes the user's inspiratory signal to the unit. The cannula tubes are interchangeable.


Figure 9: 550 Portable Unit

Filling the Portable Unit

The Companion Portable unit is filled from the Stationary unit as follows.

CAUTION

The fill connectors must be clean and dry in order to avoid malfunction due to freezing.

1. Check the contents indicator on the stationary unit to be sure that there is an ample supply of liquid oxygen for filling purposes.

NOTE: *Do not attempt to fill the Companion Portable unit if the stationary unit contents indicator is near the empty area.*

NOTE: For Companion 1000T only. *Remove Moisture Cup by rotating thumb nut $\frac{1}{4}$ turn counter-clockwise. Wring out Moisture Pad, held with Velcro®.*

2. Using a clean, dry, lint-free cloth, dry the male and female fill connectors.
3. Hold the Companion Portable unit with both hands and position the contoured case over the matching recessed area in the Stationary unit cover (Figure 10).
4. Lower the Companion Portable unit carefully into place taking care to assure proper engagement of fill connectors.


Figure 10: Filling the Portable Unit

5. Place one hand on top of the Companion Portable unit directly over the fill connector and press straight down. This will lower the Companion Portable unit approximately 10mm, and assure proper engagement of the fill connectors.
6. While holding the Portable unit in the fill position, move the vent valve lever to the open position (90° from normal off position). This will result in a rather loud hissing noise (Figure 11).

NOTE: *Note the time at the start of the fill .*

1 — Vent Valve

Figure 11: Filling the Portable Unit

WARNING

Do not leave the Companion Portable unit unattended during filling operation.

NOTE: *One hand should maintain a slight downward pressure on the unit during filling to assure stability and proper filling position. Approximately 20-30 seconds into the filling procedure it is advisable to close and reopen the vent valve one or more times. This will break up any ice that may begin to form around the valve stem and serve to avoid any problems with the vent valve freezing open.*

7. When there is a noticeable change in the sound of venting gas, followed by the emission of a dense, white vapor around the cover of the Stationary unit, close the vent valve. Fill time may vary according to the temperature of the container being filled. Maximum fill time is approximately 1 ½ minutes.

NOTE: *Should for any reason the vent valve fail to close and the hissing continue, remove the Companion Portable unit by depressing the portable release button on the Stationary unit. The Companion Portable unit will stop venting in a few minutes. Allow the unit to warm until the vent valve can be closed. A period of 2-3 hours at no flow may be required to allow the Companion Portable unit to restore adequate pressure for accurate oxygen flow.*

8. Disengage the Companion Portable unit from the Stationary unit by holding the carrying strap above the unit and depressing the release button (). Always hold the Companion Portable unit with one hand when attempting to disengage it. Should the units not disengage easily, they may have become frozen. DO NOT USE FORCE. Simply allow a few moments for the frozen parts to warm and disengage when the ice has melted.

WARNING

Should a minor liquid oxygen leakage occur when the Companion Portable unit is disengaged, proceed to re-engage and disengage the unit. This will help dislodge any ice or other obstruction. If liquid leakage is still present, engage unit again and notify your liquid oxygen supplier.


Should a major liquid oxygen leakage occur when the Companion Portable unit is disengaged, stay away from the unit and notify your liquid oxygen supplier.

9. Check the liquid oxygen contents indicator (Figure 13). The amount of liquid contained in the Companion Portable unit is measured by an internal scale which is built into the unit.

The scale is actuated by simply lifting the Companion Portable unit by the strap nearest the indicator. The liquid contents is indicated on the colour-coded gauge.


NOTE: *The plastic connector in the carrying strap must be disengaged when checking liquid oxygen contents.*

1 — Release Button

Figure 12: Disengaging the Portable Unit

1 — Contents Indicator

Figure 13: Checking Contents Indicator

NOTE: *For shorter planned use time, you can partially fill the Companion Portable unit by closing the vent valve sooner than you normally would. In this way, your Portable unit will be even lighter and you will waste less oxygen.*

For Companion 1000T only

Attach Moisture Cup by inserting fastener and rotating thumb nut clockwise until it snaps into position.

10. Place the oxygen tube on the Companion Portable oxygen outlet connector and adjust the cannula or other breathing device to the face to receive oxygen comfortably.
11. Adjust the cannula to comfortably fit your face (Figure 14). When the cannula is in place, turn the *Flow Control* to the prescribed rate.


Figure 14: Adjusting the Cannula

For Companion 550 only

Place one of the tubes from the dual-lumen oxygen cannula on the Companion 550 oxygen outlet connector (upper connector). Place the other cannula tube on the unit's inspiration sensor connector (lower connector). Adjust the cannula to the face to receive oxygen comfortably.


NOTE: *The dual-lumen cannula tubes may be interchanged on either Companion 550 connector.*

12. Turn the Companion Portable unit flow control to the prescribed rate (Figure 15).

NOTE: *Do not set the flow above the flow rate prescribed by your physician. Do not attempt to set the flow control between settings. Oxygen will flow from the unit only when the flow control is set to one of its marked increments.*


Flow Control (set for 2 L/min)

Figure 15 Flow Control Knob


For Companion 550 only

Set the Continuous / Demand Flow Switch for the desired oxygen delivery mode — continuous oxygen flow or flow only during the inspiratory part of each breath (Figure 16).

Demand Flow

(flows from 1 to 6 L/min only)

Continuous Flow

Figure 16 Continuous/Demand Flow Switch

NOTE: *The demand flow mode should be used at flow settings of 1 to 6 L/min only. Use of the demand flow mode below 1 L/min may shorten planned use time as compared to continuous flow. Slightly elevated blood oxygen levels may result if the demand flow mode is selected at flows equal to or less than 1 L/min.*

NOTE: *With the Companion 550 in the demand flow mode, you may experience continuous flow from the oxygen outlet when using the unit immediately after filling. If this is the case, normal demand flow operation should be evident within approximately 10 minutes.*

Using the Portable Unit

The Companion Portable unit may be worn on either side of the body. The unit may be worn over the shoulder or in an across the body fashion. An adjustable shoulder strap and pad are provided for maximum comfort. The plastic connector in the shoulder strap may be engaged to create a short carrying handle.

For Companion 550 only

When using the Companion 550 in the demand flow mode, oxygen is delivered when the unit senses the user's inspiratory effort. An initial “puff” of oxygen may be noticed followed by oxygen flow during the remainder of the inspiration. No oxygen flow occurs during the exhalation phase of breathing. When the unit is set for the continuous flow mode, oxygen flows during the entire breathing cycle.

At flow settings of 1 to 6 L/min, the demand flow mode can significantly increase the planned use time while reducing the drying effect that continuous flow oxygen has on the nasal passages.

CLEANING AND MAINTENANCE

Do not use alcohol, solvents, polishes, or any oily substance on oxygen equipment. If cleaning is necessary, use only warm water and a mild dish washing detergent. Dampen a cloth in the detergent and water solution and wipe the outside surfaces of the equipment until clean.

WARNING

Do not allow water into any of the controls, the fill connector, or the oxygen supply connector.

WARNING

Never attempt to repair or disassemble this equipment. You could create a hazardous condition or cause equipment failure. If you have problems, questions, or are unsure if equipment is operating properly, call your liquid oxygen system supplier.

CAUTION

Never apply cleaning agents or disinfectants to Contents Indicator. These agents attack the plastic which can lead to leaking or cracking of the indicator. Use only a cloth dampened with water to clean the indicator as needed.

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Companion is a trademark of CAIRE Inc.
Velcro is a trademark of Crown Packaging.
©2016 CAIRE Inc. All rights reserved.
Companion es una marca comercial de CAIRE Inc.
©2016 CAIRE Inc. Reservados todos los derechos.
Companion é uma marca de fábrica da CAIRE Inc.