



Innovation. Experience. Performance.®

SERVICE BULLETIN

PN: 21136744

DATE: 2/22/2017

MODEL: Liberators 30 & 60 Reservoirs

ISSUE: Empty and full weight changes

NOTES: Effective 6/10/17 the thickness of the outer shell will be increased on all new manufactured CAIRE Liberators 30's & 60's Reservoirs. This was implemented to make them more robust and to increase the longevity of the Reservoirs. Empty weights will increase from 45 lbs to 47.5 lbs for the Lib 30 and from 66 lbs to 74.5 lbs for the Lib 60. Also, full weights will increase from 120 lbs. to 122.5 lbs for Lib 30 and from 210 lbs to 219.4 lbs for the Lib 60. Below are the new specifications of the Liberators 30 and 60 Reservoirs:

Models	Weight Empty		Weight Full	
	Pounds	Kg	Pounds	Kg
Liberator 30	47.5	21.54	122.5	55.55
Liberator 60	75.4	34.19	219.4	99.5