



# *SERVICE BULLETIN*

PN: 21330667


**RELEASE DATE:** October 10, 2018

**REVISION DATE(S):** March 5, 2019 (rev B)

**MODEL:** CAIRE Liberator Reservoirs

**ISSUE:** Hazard Warning Label

**NOTES:** As of October 2018, CAIRE will no longer ship its Liberator products with the Hazard warning label affixed to the reservoir. This item will be placed in the accessory bag and customers should ensure that this label, or another suitable label, is affixed prior to transportation if the equipment is filled with liquid oxygen.



**Figure 1:** Label Delivered in Accessory Bag



**Figure 2:** Label Affixed on Reservoir

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

[www.cairemedical.com](http://www.cairemedical.com)

Europe: +44 (0) 118 936 7060

[www.cairemedical.com](http://www.cairemedical.com)

For ordering information or general inquires, contact Customer Service:

USA: 800-482-2473

[www.cairemedical.com](http://www.cairemedical.com)

Europe: +44 (0) 118 936 7060

[www.cairemedical.com](http://www.cairemedical.com)

*Copyright © 2019 CAIRE Inc. CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.*