



# *SERVICE BULLETIN*

PN: 21356419

**DATE:** 12/20/2018

**MODEL:** Stroller and Hi-Flow Stroller Portable Liquid Oxygen Devices

**ISSUE:** Symbolic Label Change

**NOTES:** Effective **Jan 1<sup>st</sup>, 2019** all newly manufactured Stroller and Hi-Flow Stroller portable liquid oxygen devices will be equipped with new symbolic warning labels. Prior to this change, CAIRE maintained a large variety of LOX labels in stock due to the number of different languages throughout our markets. The use of standard symbolic labels on LOX portables will greatly streamline and simplify the ordering process for customers, improve repair turn-around times, and simplify inventory management for both CAIRE and our customers.

All newly manufactured portables will now be equipped with auto-close vent valve going forward.


The new symbols on the warning label are defined in the accompanying user manual included with each unit.

Many of the country specific portable will be replaced with one, global item number that can be sold into multiple regions and countries. This global model will include a multi-lingual manual. Table 1 shows the new global portable models that will be available going forward (as well as unchanged final models). See Table 2 for list of replacement item numbers. Contact customer service for pricing and availability.

Table 3 shows models that are obsolete and will no longer be available.


For replacement labels for service use reference service bulletin **21174579**. Language specific labels will also no longer be available, refer to service bulletin **21174579** for the correct replacement label kit. Replacement labels required for service use will depend on manufacture date and the region of the device.

Table 4 shows language specific patient manuals that will be replaced by a multi-lingual version – MN236-C4.



**Fig 1(a). Location of New Labels on Stroller soft case and High Flow Stroller**

**Fig 1(b). Location of New Labels on Stroller**



**Fig 2. Portable Symbolic Warning Label**

**Table 1. New Global Stroller & Hi-Flow Stroller Part Numbers**

<b>Previous Item Number</b>	<b>Description</b>	<b>Description 2</b>	<b>Change</b>	<b>New Part Number</b>
10845200	FNL STROLLER T SF CA M 6LPM	20PSI METER GERMAN	New Global PN	10845200
11933113	FNL STROLLER T SF CA M 2LPM	20PSI METER GERMAN	New Global PN	11933113
10734067	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL GERMAN	New Global PN	10734067
14483752	FNL STROLLER T SF CA M 2LPM	20PSI METER SOL ITALIAN	New Global PN	14483752
14909784	FNL STROLLER T SF CA M 6LPM	20PSI SCALE GERMAN	New Global PN	14909784
10956680	FNL STROLLER T TF PB F 6LPM	20PSI METER GERMAN	New Global PN	10956680
11730957BZ	FNL STROLLER T TF PB F 6LPM	20PSI METER PORTUGUESE BRAZIL	New Global PN	11730957BZ
13259433	FNL STROLLER T TF PB F 2LPM	20PSI METER GERMAN	New Global PN	13259433
14306748	FNL STROLLER T TF PB F 2LPM	20PSI METER GERMAN SC GRY BLK	New Global PN	14306748
14306756	FNL STROLLER T TF PB F 2LPM	20PSI SCALE GERMAN SC GRY	New Global PN	14306756
12993791	FNL STROLLER T TF PB F 6LPM	20PSI SCALE GERMAN	New Global PN	12993791
11908161	FNL STROLLER T TF PB F 6LPM	20PSI SCALE SAPIO ITALIAN	New Global PN	11908161
10769147	FNL STROLLER T TF PEN F 6LPM	20PSI METER ENGLISH	New Global PN	10769147
10935723	FNL STROLLER T TF PEN F 6LPM	20PSI METER SOL ITALIAN	New Global PN	10935723
11899808BZ	FNL STROLLER T TF TAE F 6LPM	20PSI METER PORTUGUESE BRAZIL	New Global PN	11899808BZ
14779060	FNL STROLLER T SF CA M 6LPM	20PSI METER ENGLISH SC BLK	New Global PN	14779060
14768264	FNL STROLLER T TF PB F 6LPM	20PSI METER ENGLISH SC BLK	New Global PN	14768264
14765549	FNL STROLLER T TF PB F 6LPM	20PSI SCALE ENGLISH SC BLK	New Global PN	14765549
14308591	FNL STROLLER T TF PB F 6LPM	20PSI SCALE ENGLISH SC GRY	New Global PN	14308591
14778770	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL GERMAN BLK	New Global PN	14778770

Previous Item Number	Description	Description 2	Change	New Part Number
14089341	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL GERMAN GRY	New Global PN	14089341
14768299	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER GERMAN BLK	New Global PN	14768299
14081315	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER GERMAN GRY	New Global PN	14081315
14081500	FNL HI FLOW T TF PEN F 15L G4	20 PSI METER ISRAELI GRY	New Global PN	14081500
20663674	FNL HI FLOW T TF PB F 15LPM G4	20 PSI SCALE GERMAN BLK	New Global PN	20663674
11490156	FNL STROLLER C SF CA M 6LPM	20PSI METER KOIKE JAPANESE	Keep the PN	No changes to labels / Manuals

**Table 2. Replacement Stroller & Hi-Flow Stroller Part Numbers**

Previous Item Number	Description	Description 2	Change	New Part Number
14303934	FNL STROLLER C SF CA F 6LPM	50PSI METER KOIKE JAPANESE SC	Obsolete replaced by	11490156
14304152	FNL STROLLER C SF CA M 6LPM	20PSI METER KOIKE JAPANESE SC	Obsolete replaced by	11490156
10678464	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL ITALIAN	Obsolete replaced by	10734067
20705127	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL SPANISH	Obsolete replaced by	10734067
10733283	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL FRENCH	Obsolete replaced by	10734067
11202636	FNL STROLLER T SF CA M 6LPM	20PSI SCALE SPANISH	Obsolete replaced by	14909784
10998054	FNL STROLLER T SF CA M 6LPM	20PSI SCALE ENGLISH	Obsolete replaced by	14909784
14259427	FNL STROLLER T SF CA M 6LPM	20PSI METER ITALIAN	Obsolete replaced by	14909784
10564109	FNL STROLLER T TF PB F 6LPM	20PSI METER ENGLISH	Obsolete replaced by	10956680
10875556	FNL STROLLER T TF PB F 6LPM	20PSI METER FRENCH	Obsolete replaced by	10956680
10734593	FNL STROLLER T TF PB F 6LPM	20PSI METER SPANISH	Obsolete replaced by	10956680
10924469	FNL STROLLER T TF PB F 6LPM	20PSI METER SPANISH	Obsolete replaced by	10956680
11488267	FNL STROLLER T TF PB F 6LPM	20PSI METER ITALIAN	Obsolete replaced by	10956680

Previous Item Number	Description	Description 2	Change	New Part Number
10815861	FNL STROLLER T TF PB F 6LPM	20PSI METER DUTCH	Obsolete replaced by	10956680
10736548	FNL STROLLER T TF PB F 6LPM	20PSI METER SWEDISH	Obsolete replaced by	10956680
11730957	FNL STROLLER T TF PB F 6LPM	20PSI METER PORTUGUESE	Obsolete replaced by	10956680
11873771	FNL STROLLER T TF PB F 6LPM	20PSI METER HUNGARIAN	Obsolete replaced by	10956680
11873704	FNL STROLLER T TF PB F 6LPM	20PSI METER GREEK	Obsolete replaced by	10956680
11046696	FNL STROLLER T TF PB F 6LPM	20PSI METER FINNISH	Obsolete replaced by	10956680
10776873	FNL STROLLER T TF PB F 6LPM	20PSI METER DANISH	Obsolete replaced by	10956680
11928402	FNL STROLLER T TF PB F 2LPM	20PSI METER ENGLISH	Obsolete replaced by	13259433
14306730	FNL STROLLER T TF PB F 2LPM	20PSI METER ENGLISH SC GRY	Obsolete replaced by	13259433
10951563	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL GERMAN	Obsolete replaced by	10956680
11735336	FNL STROLLER T TF PB F 6LPM	20PSI SCALE ITALIAN	Obsolete replaced by	12993791
11063816	FNL STROLLER T TF PB F 6LPM	20PSI SCALE SPANISH	Obsolete replaced by	12993791
10997991	FNL STROLLER T TF PB F 6LPM	20PSI SCALE ENGLISH	Obsolete replaced by	12993791
10564125	FNL STROLLER T SF CA M 6LPM	20PSI METER ENGLISH	Obsolete replaced by	10845200
20892874	FNL STROLLER T SF CA M 6LPM	20PSI METER SWEDISH	Obsolete replaced by	10845200
10724985	FNL STROLLER T SF CA M 6LPM	20PSI METER SPANISH	Obsolete replaced by	10845200
10734340	FNL STROLLER T SF CA M 6LPM	20PSI METER DANISH	Obsolete replaced by	10845200
10733216	FNL STROLLER T SF CA M 6LPM	20PSI METER FRENCH	Obsolete replaced by	10845200
10734817	FNL STROLLER T SF CA M 6LPM	20PSI METER DUTCH	Obsolete replaced by	10845200
14767456	FNL STROLLER T TF PB F 6LPM	20PSI M MEDICAIR ITAL SC BLK	Obsolete replaced by	14768264
14767510	FNL STROLLER T TF PB F 6LPM	20PSI METER BOC ENGLISH SC BLK	Obsolete replaced by	14768264
14767536	FNL STROLLER T TF PB F 6LPM	20PSI METER NORWEGIAN SC BLK	Obsolete replaced by	14768264

Previous Item Number	Description	Description 2	Change	New Part Number
14767552	FNL STROLLER T TF PB F 6LPM	20PSI METER FRENCH SC BLK	Obsolete replaced by	14768264
14769814	FNL STROLLER T TF PB F 6LPM	20PSI METER SPANISH SC BLK	Obsolete replaced by	14768264
14769822	FNL STROLLER T TF PB F 6LPM	20 PSI METER SWEDISH SC BLK	Obsolete replaced by	14768264
14769831	FNL STROLLER T TF PB F 6LPM	20PSI METER DANISH SC BLK	Obsolete replaced by	14768264
14769865	FNL STROLLER T TF PB F 6LPM	20PSI METER DUTCH SC BLK	Obsolete replaced by	14768264
14769881	FNL STROLLER T TF PB F 6LPM	20PSI METER SPANISH SC BLK	Obsolete replaced by	14768264
14769890	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL GERMAN SC BLK	Obsolete replaced by	14768264
14769902	FNL STROLLER T TF PB F 6LPM	20PSI METER GERMAN SC BLK	Obsolete replaced by	14768264
14769911	FNL STROLLER T TF PB F 6LPM	20PSI METER DVLBS ENG SC BLK	Obsolete replaced by	14768264
14769945	FNL STROLLER T TF PB F 6LPM	20PSI METER FINNISH SC BLK	Obsolete replaced by	14768264
14769961	FNL STROLLER T TF PB F 6LPM	20PSI METER ITALIAN SC BLK	Obsolete replaced by	14768264
14769970	FNL STROLLER T TF PB F 6LPM	20PSI METER PORTUGUESE SC BLK	Obsolete replaced by	14768264
14769988	FNL STROLLER T TF PB F 6LPM	20PSI METER MEDIGAS ITA SC BLK	Obsolete replaced by	14768264
14769996	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL DUTCH SC BLK	Obsolete replaced by	14768264
14770516	FNL STROLLER T TF PB F 6LPM	20PSI METER GREEK SC BLK	Obsolete replaced by	14768264
14777486	FNL STROLLER T TF PB F 6LPM	20PSI METER HUNGARIAN SC BLK	Obsolete replaced by	14768264
14777494	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL FRENCH SC BLK	Obsolete replaced by	14768264
14777523	FNL STROLLER T TF PB F 6LPM	20PSI METER ISRAELI SC BLK	Obsolete replaced by	14768264
14777540	FNL STROLLER T TF PB F 6LPM	20PSI METER TURKISH SC BLK	Obsolete replaced by	14768264
14168322	FNL STROLLER T SF CA M 6LPM	SOL GREEK	Obsolete replaced by	14779060

Previous Item Number	Description	Description 2	Change	New Part Number
14304179	FNL STROLLER T SF CA M 6LPM	20PSI METER ENGLISH SC GRY	Obsolete replaced by	14779060
14304195	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL ITALIAN SC GRY	Obsolete replaced by	14779060
14304208	FNL STROLLER T SF CA M 6LPM	20PSI METER CARBUROS ES SC GRY	Obsolete replaced by	14779060
14304216	FNL STROLLER T SF CA M 6LPM	20PSI METER SPANISH SC GRY	Obsolete replaced by	14779060
10733697	NL STROLLER T SF CA M 6LPM	20PSI METER SOL DUTCH	Obsolete replaced by	10734067
11676225	FNL STROLLER T TF PEN F 6LPM	20PSI METER ISRAELI	Obsolete replaced by	10769147
13328937	FNL STROLLER T SF CA M 6LPM	20PSI SCALE ITALIAN	Obsolete replaced by	10845200
10716619	FNL STROLLER T SF CA M 6LPM	20PSI METER CARBUROS ES	Obsolete replaced by	10845200
10733873	FNL STROLLER T SF CA M 6LPM	20PSI METER SWEDISH	Obsolete replaced by	10845200
10734278	FNL STROLLER T SF CA M 6LPM	20PSI METER ENGLISH	Obsolete replaced by	10845200
10948727	FNL STROLLER T SF CA M 6LPM	20PSI METER GREEK	Obsolete replaced by	10845200
11015814	FNL STROLLER T SF CA M 6LPM	20PSI METER PORTUGUESE	Obsolete replaced by	10845200
10733398	FNL STROLLER T TF PB F 6LPM	20PSI METER MEDICAIR ITALIAN	Obsolete replaced by	10956680
10733451	FNL SPRINT T TF PB F 6LPM	20PSI METER BOC ENGLISH	Obsolete replaced by	10956680
10733460	FNL STROLLER T TF PB F 6LPM	20PSI METER BOC ENGLISH	Obsolete replaced by	10956680
10733718	FNL STROLLER T TF PB F 6LPM	20PSI METER NORWEGIAN	Obsolete replaced by	10956680
13783551	FNL STROLLER T TF PB F 6LPM	20PSI METER ISRAELI	Obsolete replaced by	10956680
13957397	FNL STROLLER T TF PB F 6LPM	20PSI METER TURKISH	Obsolete replaced by	10956680
14807374	FNL STROLLER T TF PB F 6LPM	20PSI METER FRENCH W/A-C	Obsolete replaced by	10956680
11008104	FNL STROLLER T TF PB F 6LPM	20PSI SCALE DVLBS ENGLISH	Obsolete replaced by	11908161
14516841	FNL STROLLER T TF PB F 6LPM	PSI METER SAPIO ITALIAN	Obsolete replaced by	11908161

Previous Item Number	Description	Description 2	Change	New Part Number
14779318	FNL STROLLER T SF CA M 2LPM	20PSI METER SOL GERMAN SC BLK	Obsolete replaced by	11933113
11707781	FNL STROLLER T TF PB F 6LPM	20PSI SCALE MEDICAIR ITALIAN	Obsolete replaced by	12993791
12965820	FNL STROLLER T TF PB F 6LPM	20PSI SCALE PORTUGUESE	Obsolete replaced by	12993791
20598672	FNL STROLLER T TF PB F 2LPM	20PSI METER ITALIAN	Obsolete replaced by	13259433
14765522	FNL STROLLER T TF PB F 2LPM	20PSI SCALE GERMAN SC BLK	Obsolete replaced by	14306748
14777718	FNL STROLLER T TF PB F 2LPM	20PSI METER ENGLISH SC BLK	Obsolete replaced by	14306748
14777726	FNL STROLLER T TF PB F 2LPM	20PSI METER GERMAN SC BLK	Obsolete replaced by	14306748
14765531	FNL STROLLER T TF PB F 6LPM	20PSI SCALE DVLBS ENG SC BLK	Obsolete replaced by	14765549
14765637	FNL STROLLER T TF PB F 6LPM	20PSI SCALE SPANISH SC BLK	Obsolete replaced by	14765549
14766841	FNL STROLLER T TF PB F 6LPM	20PSI MEDICAIR ITAL SC BLK	Obsolete replaced by	14765549
14766859	FNL STROLLER T TF PB F 6LPM	20PSI SCALE ITALIAN SC BLK	Obsolete replaced by	14765549
14766867	FNL STROLLER T TF PB F 6LPM	20PSI SCALE SAPIO ITAL SC BLK	Obsolete replaced by	14765549
14766875	FNL STROLLER T TF PB F 6LPM	20PSI SCALE PORTUGUESE SC BLK	Obsolete replaced by	14765549
14766883	FNL STROLLER T TF PB F 6LPM	20PSI SCALE GERMAN SC BLK	Obsolete replaced by	14765549
14779326	FNL STROLLER T SF CA M 2LPM	20PSI METER GERMAN SC BLK	Obsolete replaced by	14765549
14779078	NL STROLLER T SF CA M 6LPM	20PSI METER SOL ITALIAN SC BLK	Obsolete replaced by	14779060
14779086	FNL STROLLER T SF CA M 6LPM	20PSI METER CARBUROS ES SC BLK	Obsolete replaced by	14779060
14779107	FNL STROLLER T SF CA M 6LPM	20PSI METER SPANISH SC BLK	Obsolete replaced by	14779060
14779115	FNL STROLLER T SF CA M 6LPM	20PSI METER FRENCH SC BLK	Obsolete replaced by	14779060



Previous Item Number	Description	Description 2	Change	New Part Number
14779123	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL FRENCH SC BLK	Obsolete replaced by	14779060
14779158	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL DUTCH SC BLK	Obsolete replaced by	14779060
14779166	FNL STROLLER T SF CA M 6LPM	20PSI METER SWEDISH SC BLK	Obsolete replaced by	14779060
14779174	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL GERMAN SC BLK	Obsolete replaced by	14779060
14779203	FNL STROLLER T SF CA M 6LPM	20PSI METER DANISH SC BLK	Obsolete replaced by	14779060
14779211	FNL STROLLER T SF CA M 6LPM	20PSI METER DUTCH SC BLK	Obsolete replaced by	14779060
14779220	FNL STROLLER T SF CA M 6LPM	20PSI METER GERMAN SC BLK	Obsolete replaced by	14779060
14779246	FNL STROLLER T SF CA M 6LPM	20PSI METER GREEK SC BLK	Obsolete replaced by	14779060
14779254	FNL STROLLER T SF CA M 6LPM	20PSI METER DVLBS ENG SC BLK	Obsolete replaced by	14779060
14779262	FNL STROLLER T SF CA M 6LPM	20PSI METER PORTUGUESE SC BLK	Obsolete replaced by	14779060
14779271	FNL STROLLER T SF CA M 6LPM	SOL GREEK SC BLK	Obsolete replaced by	14779060
14779289	FNL STROLLER T SF CA M 6LPM	20PSI METER ITALIAN SC BLK	Obsolete replaced by	14779060
14779342	FNL STROLLER T SF CA M 6LPM	SOL GREEK SC BLK	Obsolete replaced by	14779060
11913981	FNL STROLLER T SF CA M 6LPM	20PSI SCALE CARBUROS SPANISH	Obsolete replaced by	14909784
11563503	FNL STROLLER T TF PB F 6LPM	PSI METER SOL DUTCH	Obsolete replaced by	10956680
13733772	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL FRENCH	Obsolete replaced by	10956680
20705125	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL SPANISH	Obsolete replaced by	10956680
14308698	FNL STROLLER T TF PB F 6LPM	20PSI SCALE DVLBS ENG SC GRY	Obsolete replaced by	14308591
14308743	FNL STROLLER T TF PB F 6LPM	20PSI SCALE SPANISH SC GRY	Obsolete replaced by	14308591
14308815	FNL STROLLER T TF PB F 6LPM	20PSI SCALE MEDIGAS ITA SC GRY	Obsolete replaced by	14308591

Previous Item Number	Description	Description 2	Change	New Part Number
14308874	FNL STROLLER T TF PB F 6LPM	20PSI S MEDICAIR ITAL SC GRY	Obsolete replaced by	14308591
14308891	FNL STROLLER T TF PB F 6LPM	I SCALE ITALIAN SC GRY	Obsolete replaced by	14308591
14308946	FNL STROLLER T TF PB F 6LPM	20PSI SCALE SAPIO ITAL SC GRY	Obsolete replaced by	14308591
14308954	FNL STROLLER T TF PB F 6LPM	20PSI SCALE PORTUGUESE SC GRY	Obsolete replaced by	14308591
14308962	FNL STROLLER T TF PB F 6LPM	20PSI SCALE GERMAN SC GRY	Obsolete replaced by	14308591
14304161	FNL STROLLER T SF CA M 2LPM	20PSI METER SOL GERMAN SC GRY	Obsolete replaced by	14306756
14304187	FNL STROLLER T SF CA M 2LPM	20PSI METER GERMAN SC GRY	Obsolete replaced by	14306756
11925113	FNL STROLLER T SF CA M 2LPM	20PSI METER SOL GERMAN	Obsolete replaced by	14483752
14778948	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER GERMAN BLK	Obsolete replaced by	14778770
14779019	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER ENGLISH BLK	Obsolete replaced by	14778770
14777849	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL ITALIAN BLK	Obsolete replaced by	14778770
14777822	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SPANISH BLK	Obsolete replaced by	14778770
14777857	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER FRENCH BLK	Obsolete replaced by	14778770
20636855	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER NORWEGIAN BLK	Obsolete replaced by	14778770
14778825	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER DANISH BLK	Obsolete replaced by	14778770
14777873	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL FRENCH BLK	Obsolete replaced by	14778770
14190651	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER ENGLISH GRY	Obsolete replaced by	14089341
14089384	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER GERMAN GRY	Obsolete replaced by	14089341
14089309	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER FRENCH GRY	Obsolete replaced by	14089341
14089296	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL ITALIAN GRY	Obsolete replaced by	14089341

Previous Item Number	Description	Description 2	Change	New Part Number
20679129	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL ITALIAN GRY	Obsolete replaced by	14089341
14767421	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER ENGLISH BLK	Obsolete replaced by	14768299
14768379	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER ITALIAN BLK	Obsolete replaced by	14768299
14767499	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER NORWEGIAN BLK	Obsolete replaced by	14768299
14768256	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SOL GERMAN BLK	Obsolete replaced by	14768299
14767464	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SPANISH BLK	Obsolete replaced by	14768299
14768352	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER FINNISH BLK	Obsolete replaced by	14768299
14768416	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER HUNGARIAN BLK	Obsolete replaced by	14768299
14081227	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER ENGLISH GRY	Obsolete replaced by	14081315
14081260	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SPANISH GRY	Obsolete replaced by	14081315
14081323	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER FRENCH GRY	Obsolete replaced by	14081315
14081382	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER ITALIAN GRY	Obsolete replaced by	14081315
14081391	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER PORTUGUESE GRY	Obsolete replaced by	14081315
14081294	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER DUTCH GRY	Obsolete replaced by	14081315
14818111	FNL HI FLOW T TF PB F 15LPM G4	20PSI SCALE ENGLISH BLACK	Obsolete replaced by	20663674
14081471	FNL HI FLOW T TF PEN F 15L G4	20 PSI METER ENGLISH GRY	Obsolete replaced by	14081500
14081489	FNL HI FLOW T TF PEN F 15L G4	20 PSI METER SOL ITALIAN GRY	Obsolete replaced by	14081500
14081243	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER BOC ENGLISH GRY	Obsolete replaced by	14081315
14081251	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER NORWEGIAN GRY	Obsolete replaced by	14081315
14081278	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SWEDISH GRY	Obsolete replaced by	14081315
14081286	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER DANISH GRY	Obsolete replaced by	14081315

Previous Item Number	Description	Description 2	Change	New Part Number
14081307	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SOL GERMAN GRY	Obsolete replaced by	14081315
14081366	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER FINNISH GRY	Obsolete replaced by	14081315
14081411	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SOL DUTCH GRY	Obsolete replaced by	14081315
14081420	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER GREEK GRY	Obsolete replaced by	14081315
14081438	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER HUGARIAN GRY	Obsolete replaced by	14081315
14081446	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SOL FRENCH GRY	Obsolete replaced by	14081315
14081454	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER ISRAELI GRY	Obsolete replaced by	14081315
14081462	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER TURKISH GRY	Obsolete replaced by	14081315
14089261	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SPANISH GRY	Obsolete replaced by	14089341
14089270	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER CARBUROS ES GRY	Obsolete replaced by	14089341
14089288	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER ITALIAN GRY	Obsolete replaced by	14089341
14089317	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL FRENCH GRY	Obsolete replaced by	14089341
14089325	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL DUTCH GRY	Obsolete replaced by	14089341
14089333	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SWEDISH GRY	Obsolete replaced by	14089341
14089368	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER DANISH GRY	Obsolete replaced by	14089341
14089376	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER DUTCH GRY	Obsolete replaced by	14089341
14089392	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER GREEK GRY	Obsolete replaced by	14089341
14089413	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER PORTUGESE GRY	Obsolete replaced by	14089341
14767448	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER BOC ENGLISH BLK	Obsolete replaced by	14768299
14767561	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SWEDISH BLK	Obsolete replaced by	14768299
14767659	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER DANISH BLK	Obsolete replaced by	14768299

Previous Item Number	Description	Description 2	Change	New Part Number
14768205	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER DUTCH BLK	Obsolete replaced by	14768299
14768281	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER FRENCH BLK	Obsolete replaced by	14768299
14768387	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER PORTUGUESE BLK	Obsolete replaced by	14768299
14768395	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SOL DUTCH BLK	Obsolete replaced by	14768299
14768408	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER GREEK BLK	Obsolete replaced by	14768299
14768441	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER SOL FRENCH BLK	Obsolete replaced by	14768299
14768459	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER ISRAELI BLK	Obsolete replaced by	14768299
14768467	FNL HI FLOW T TF PB F 15LPM G4	20 PSI METER TURKISH BLK	Obsolete replaced by	14768299
14777831	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER ITALIAN BLK	Obsolete replaced by	14778770
14777890	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SOL DUTCH BLK	Obsolete replaced by	14778770
14777902	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER SWEDISH BLK	Obsolete replaced by	14778770
14778833	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER DUTCH BLK	Obsolete replaced by	14778770
14778913	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER CARBUROS ES BLK	Obsolete replaced by	14778770
14778981	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER GREEK BLK	Obsolete replaced by	14778770
14779001	FNL HI FLOW T SF CA M 15LPM G4	20 PSI METER PORTUGUESE BLK	Obsolete replaced by	14778770

**Table 3. Obsolete Stroller & Hi-Flow Stroller Part Numbers**

Previous Item Number	Description	Description 2	Change	New Part Number
10768671	FNL STROLLER T TF CR2 M 6LPM	20PSI METER ENGLISH	Obsolete no replacement	No R/B
10987048	FNL STROLLER T TF CR2 M 6LPM	20PSI METER DVLBS ENGLISH	Obsolete no replacement	No R/B
10987056	FNL STROLLER T TF PEN F 6LPM	20PSI METER DVLBS ENGLISH	Obsolete no replacement	No R/B

Previous Item Number	Description	Description 2	Change	New Part Number
10998003	FNL STROLLER T TF CR2 M 6LPM	20PSI SCALE ENGLISH	Obsolete no replacement	No R/B
11008091	FNL STROLLER T TF CR2 M 6LPM	20PSI SCALE DVLBS ENGLISH	Obsolete no replacement	No R/B
11899808	FNL STROLLER T TF TAE F 6LPM	20PSI METER PORTUGUESE	Obsolete no replacement	No R/B
11899816	FNL STROLLER T TF TAE F 6LPM	20PSI METER SPANISH	Obsolete no replacement	No R/B
11935725	FNL STROLLER T TF TAE F 2LPM	20PSI METER GERMAN	Obsolete no replacement	No R/B
13997962	FNL STROLLER T TF USP 93%	20PSI SCALE INSPIRED TECH	Obsolete no replacement	No R/B
14304224	FNL STROLLER T SF CA M 6LPM	20PSI METER FRENCH SC GRY	Obsolete no replacement	No R/B
14304232	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL FRENCH SC GRY	Obsolete no replacement	No R/B
14304241	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL DUTCH SC GRY	Obsolete no replacement	No R/B
14304275	FNL STROLLER T SF CA M 6LPM	20PSI METER SWEDISH SC GRY	Obsolete no replacement	No R/B
14304419	FNL STROLLER T SF CA M 6LPM	20PSI METER SOL GERMAN SC GRY	Obsolete no replacement	No R/B
14304443	FNL STROLLER T SF CA M 6LPM	20PSI METER DANISH SC GRY	Obsolete no replacement	No R/B
14304451	FNL STROLLER T SF CA M 6LPM	20PSI METER DUTCH SC GRY	Obsolete no replacement	No R/B
14304460	FNL STROLLER T SF CA M 6LPM	20PSI METER GERMAN SC GRY	Obsolete no replacement	No R/B
14304478	FNL STROLLER T SF CA M 6LPM	20PSI METER GREEK SC GRY	Obsolete no replacement	No R/B
14304486	FNL STROLLER T SF CA M 6LPM	20PSI METER DVBLBS ENG SC GRY	Obsolete no replacement	No R/B
14304515	FNL STROLLER T SF CA M 6LPM	20PSI SCALE ENGLISH SC GRY	Obsolete no replacement	No R/B
14304540	FNL STROLLER T SF CA M 6LPM	20PSI SCALE DVLBS ENG SC GRY	Obsolete no replacement	No R/B
14304558	FNL STROLLER T SF CA M 6LPM	20PSI METER PORTUGUESE SC GRY	Obsolete no replacement	No R/B
14304574	FNL STROLLER T SF CA M 6LPM	20PSI SCALE SOL ITALIAN SC GRY	Obsolete no replacement	No R/B

Previous Item Number	Description	Description 2	Change	New Part Number
14304582	FNL STROLLER T SF CA M 6LPM	20PSI SCALE SPANISH SC GRY	Obsolete no replacement	No R/B
14306588	FNL STROLLER T SF CA M 6LPM	20PSI SCALE CARBUROS ES SC GRY	Obsolete no replacement	No R/B
14306596	FNL STROLLER T SF CA M 6LPM	20PSI SCALE ITALIAN SC GRY	Obsolete no replacement	No R/B
14306617	FNL STROLLER T SF CA M 6LPM	SOL GREEK SC GRY	Obsolete no replacement	No R/B
14306641	FNL STROLLER T SF CA M 6LPM	20PSI METER ITALIAN SC GRY	Obsolete no replacement	No R/B
14306772	FNL STROLLER T TF PB F 6LPM	20PSI METER ENGLISH SC GRY	Obsolete no replacement	No R/B
14306781	FNL STROLLER T TF PB F 6LPM	SI M MEDICAIR ITAL SC GRY	Obsolete no replacement	No R/B
14306799	FNL STROLLER T TF PB F 6LPM	20PSI METER BOC ENGLISH SC GRY	Obsolete no replacement	No R/B
14306801	FNL STROLLER T TF PB F 6LPM	20PSI METER NORWEGIAN SC GRY	Obsolete no replacement	No R/B
14306810	FNL STROLLER T TF PB F 6LPM	20PSI METER SPANISH SC GRY	Obsolete no replacement	No R/B
14306828	NL STROLLER T TF PB F 6LPM	20PSI METER SWEDISH SC GRY	Obsolete no replacement	No R/B
14306836	FNL STROLLER T TF PB F 6LPM	20PSI METER DANISH SC GRY	Obsolete no replacement	No R/B
14306844	FNL STROLLER T TF PB F 6LPM	20PSI METER DUTCH SC GRY	Obsolete no replacement	No R/B
14306852	FNL STROLLER T TF PB F 6LPM	20PSI METER FRENCH SC GRY	Obsolete no replacement	No R/B
14306887	FNL STROLLER T TF PB F 6LPM	20PSI METER SPANISH SC GRY	Obsolete no replacement	No R/B
14306895	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL GERMAN SC GRY	Obsolete no replacement	No R/B
14308567	FNL STROLLER T TF PB F 6LPM	20PSI METER GERMAN SC GRY	Obsolete no replacement	No R/B
14308575	FNL STROLLER T TF PB F 6LPM	20PSI METER DVLBS ENG SC GRY	Obsolete no replacement	No R/B
14308727	FNL STROLLER T TF PB F 6LPM	20PSI METER FINNISH SC GRY	Obsolete no replacement	No R/B

Previous Item Number	Description	Description 2	Change	New Part Number
14308760	FNL STROLLER T TF PB F 6LPM	20PSI METER ITALIAN SC GRY	Obsolete no replacement	No R/B
14308786	FNL STROLLER T TF PB F 6LPM	20PSI METER PORTUGUESE SC GRY	Obsolete no replacement	No R/B
14308794	FNL STROLLER T TF PB F 6LPM	20PSI METER MEDIGAS ITA SC GRY	Obsolete no replacement	No R/B
14308831	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL DUTCH SC GRY	Obsolete no replacement	No R/B
14308911	FNL STROLLER T TF PB F 6LPM	SI METER ENGLISH SC GRY	Obsolete no replacement	No R/B
14308920	FNL STROLLER T TF PB F 6LPM	20PSI METER GREEK SC GRY	Obsolete no replacement	No R/B
14308938	FNL STROLLER T TF PB F 6LPM	20PSI METER HUNGARIAN SC GRY	Obsolete no replacement	No R/B
14309455	FNL STROLLER T TF PB F 6LPM	20PSI METER SOL FRENCH SC GRY	Obsolete no replacement	No R/B
14309471	FNL STROLLER T TF PB F 6LPM	20PSI METER ISRAELI SC GRY	Obsolete no replacement	No R/B
14309480	FNL STROLLER T TF PB F 6LPM	20PSI METER TURKISH SC GRY	Obsolete no replacement	No R/B
14309498	FNL STROLLER T TF PEN F 6LPM	PSI METER ENGLISH SC GRY	Obsolete no replacement	No R/B
14309501	FNL STROLLER T TF PEN F 6LPM	20PSI METER SOL ITALIAN SC GRY	Obsolete no replacement	No R/B
14309519	FNL STROLLER T TF PEN F 6LPM	20PSI METER DVLBS ENG SC GRY	Obsolete no replacement	No R/B
14766824	FNL STROLLER T TF PB F 6LPM	20PSI SCALE MEDIGAS ITA SC BLK	Obsolete no replacement	No R/B
14309586	FNL STROLLER T TF PEN F 6LPM	20PSI SCALE ENGLISH SC GRY	Obsolete no replacement	No R/B
14309594	FNL STROLLER T TF PEN F 6LPM	20PSI SCALE DVLBS ENG SC GRY	Obsolete no replacement	No R/B
14309631	FNL STROLLER T TF PEN F 6LPM	20PSI METER ISRAELI SC GRY	Obsolete no replacement	No R/B
14309658	FNL STROLLER T TF TAE F 2LPM	20PSI METER GERMAN SC GRY	Obsolete no replacement	No R/B
14309682	FNL STROLLER T TF TAE F 6LPM	20PSI METER GERMAN SC GRY	Obsolete no replacement	No R/B



Previous Item Number	Description	Description 2	Change	New Part Number
14309691	NL STROLLER T TF TAE F 6LPM	20PSI METER PORTUGUESE SC GRY	Obsolete no replacement	No R/B
14309703	FNL STROLLER T TF TAE F 6LPM	20PSI METER SPANISH SC SC GRY	Obsolete no replacement	No R/B
14309711	FNL STROLLER T TF USP 93%	20PSI SCALE INSPIRED TECH SC	Obsolete no replacement	No R/B
14448236	FNL STROLLER T TF CR2 M 6LPM	20PSI METER ISRAELI	Obsolete no replacement	No R/B
14766891	FNL STROLLER T TF PEN F 6LPM	20PSI SCALE ENGLISH SC BLK	Obsolete no replacement	No R/B
14767270	FNL STROLLER T TF PEN F 6LPM	20PSI SCALE DVLBS ENG SC BLK	Obsolete no replacement	No R/B
14767309	FNL STROLLER T SF CA M 6LPM	20PSI SCALE ENGLISH SC BLK	Obsolete no replacement	No R/B
14767325	FNL STROLLER T SF CA M 6LPM	20PSI SCALE DVLBS ENG SC BLK	Obsolete no replacement	No R/B
14767333	FNL STROLLER T SF CA M 6LPM	20PSI SCALE SOL ITALIAN SC BLK	Obsolete no replacement	No R/B
14767341	FNL STROLLER T SF CA M 6LPM	20PSI SCALE SPANISH SC BLK	Obsolete no replacement	No R/B
14767384	FNL STROLLER T SF CA M 6LPM	20PSI SCALE CARBUROS ES SC BLK	Obsolete no replacement	No R/B
14767405	FNL STROLLER T SF CA M 6LPM	20PSI SCALE ITALIAN SC BLK	Obsolete no replacement	No R/B
14777558	FNL STROLLER T TF PEN F 6LPM	20PSI METER ENGLISH SC BLK	Obsolete no replacement	No R/B
14777582	FNL STROLLER T TF PEN F 6LPM	20PSI METER SOL ITALIAN SC BLK	Obsolete no replacement	No R/B
14777591	FNL STROLLER T TF PEN F 6LPM	20PSI METER DVLBS ENG SC BLK	Obsolete no replacement	No R/B
14777603	FNL STROLLER T TF PEN F 6LPM	20PSI METER ISRAELI SC BLK	Obsolete no replacement	No R/B
14777620	FNL STROLLER T TF TAE F 6LPM	20PSI METER GERMAN SC BLK	Obsolete no replacement	No R/B
14777654	FNL STROLLER T TF TAE F 6LPM	20PSI METER PORTUGUESE SC BLK	Obsolete no replacement	No R/B
14777662	FNL STROLLER T TF TAE F 6LPM	20PSI METER SPANISH SC SC BLK	Obsolete no replacement	No R/B

Previous Item Number	Description	Description 2	Change	New Part Number
1477751	FNL STROLLER T TF TAE F 2LPM	20PSI METER GERMAN SC BLK	Obsolete no replacement	No R/B
14781557	FNL STROLLER T TF TAE F 6LPM	20PSI METER ITALIAN	Obsolete no replacement	No R/B
14842516	FNL STROLLER T TF TAE F 6LPM	20PSI METER ISRAELI	Obsolete no replacement	No R/B
14942411	FNL STROLLER T SF CA M 6LPM	20PSI SCALE GERMAN SC BLK	Obsolete no replacement	No R/B
14942437	FNL STROLLER T SF CA M 6LPM	20PSI SCALE GERMAN SC GRAY	Obsolete no replacement	No R/B
14306650	FNL STROLLER T TF CR2 M 6LPM	20PSI METER ENGLISH SC	Obsolete no replacement	No R/B
14306676	FNL STROLLER T TF CR2 M 6LPM	20PSI METER DVLBS ENGLISH SC	Obsolete no replacement	No R/B
14306692	FNL STROLLER T TF CR2 M 6LPM	20PSI SCALE ENGLISH SC	Obsolete no replacement	No R/B
14306705	FNL STROLLER T TF CR2 M 6LPM	20PSI SCALE DVLBS ENGLISH SC	Obsolete no replacement	No R/B
14306713	FNL STROLLER T TF CR2 M 6LPM	20PSI SCALE ITALIAN SC	Obsolete no replacement	No R/B
10998011	FNL STROLLER T TF PEN F 6LPM	20PSI SCALE ENGLISH	Obsolete no replacement	No R/B
11008112	FNL STROLLER T TF PEN F 6LPM	20PSI SCALE DVLBS ENGLISH	Obsolete no replacement	No R/B
10987064	FNL STROLLER T SF CA M 6LPM	20PSI METER DVBLS ENGLISH	Obsolete no replacement	No R/B
11008075	FNL STROLLER T SF CA M 6LPM	20PSI SCALE DVLBS ENGLISH	Obsolete no replacement	No R/B
13798189	FNL STROLLER T TF LO 6LPM	20PSI METER ITALIAN	Obsolete no replacement	No R/B
14306721	FNL STROLLER T TF LO 6LPM	20PSI METER ITALIAN SC GRY	Obsolete no replacement	No R/B
11897781	FNL STROLLER T TF TAE F 6LPM	20PSI METER GERMAN	Obsolete no replacement	No R/B
10987030	FNL STROLLER T TF PB F 6LPM	20PSI METER DVLBS ENGLISH	Obsolete no replacement	No R/B
14769806	FNL STROLLER T TF LO 6LPM	20PSI METER ITALIAN SC BLK	Obsolete no replacement	No R/B
11193063	FNL STROLLER T SF CA M 6LPM	20PSI SCALE SOL ITALIAN	Obsolete no replacement	No R/B

Previous Item Number	Description	Description 2	Change	New Part Number
13733684	FNL STROLLER T TF PB F 2LPM	20PSI SCALE GERMAN	Obsolete no replacement	No R/B
11544919	FNL STROLLER T TF PB F 6LPM	PSI METER MEDIGAS ITALIAN	Obsolete no replacement	No R/B
11547095	FNL STROLLER T TF PB F 6LPM	0PSI SCALE MEDIGAS ITALIAN	Obsolete no replacement	No R/B
11902894	FNL STROLLER T TF CR2 M 6LPM	20PSI SCALE ITALIAN	Obsolete no replacement	No R/B
10743941	FNL HI FLOW T SF CA M 15LPM	20PSI METER SOL DUTCH G2.5	Obsolete no replacement	No R/B
14769654	FNL HI FLOW T TF PEN F 15L G4	20 PSI METER ENGLISH BLK	Obsolete no replacement	No R/B
14769662	FNL HI FLOW T TF PEN F 15L G4	20 PSI METER SOL ITALIAN BLK	Obsolete no replacement	No R/B
14769671	FNL HI FLOW T TF PEN F 15L G4	20 PSI METER ISRAELI BLK	Obsolete no replacement	No R/B
14081518	FNL HI FLOW T TF TAE F 15L G4	20 PSI METER GERMAN GRY	Obsolete no replacement	No R/B
14081526	FNL HI FLOW T TF TAE F 15L G4	20 PSI METER PORTUGESE GRY	Obsolete no replacement	No R/B
14081534	FNL HI FLOW T TF TAE F 15L G4	20 PSI METER SPANISH GRY	Obsolete no replacement	No R/B
14769726	FNL HI FLOW T TF TAE F 15L G4	20 PSI METER GERMAN BLK	Obsolete no replacement	No R/B
14769734	FNL HI FLOW T TF TAE F 15L G4	20 PSI METER PORTUGUESE BLK	Obsolete no replacement	No R/B
14769742	FNL HI FLOW T TF TAE F 15L G4	20 PSI METER SPANISH BLK	Obsolete no replacement	No R/B
20651171	FNL SPIRIT 300T SF CA M 5LPM	20PSI SCALE DANISH	Obsolete no replacement	No R/B
14768512	FNL HI FLOW T TF CR2 15LPM G4	20 PSI METER ENGLISH BLK	Obsolete no replacement	No R/B
14081198	FNL HI FLOW T TF CR2 15LPM G4	20 PSI METER ENGLISH GRY	Obsolete no replacement	No R/B
14081219	FNL HI FLOW T TF LO 15LPM G4	20 PSI METER ITALIAN GRY	Obsolete no replacement	No R/B
14769241	FNL HI FLOW T TF LO 15LPM G4	20 PSI METER ITALIAN BLK	Obsolete no replacement	No R/B
20551921	FNL HI FLOW T TF PB F 15LPM G4	20PSI SCALE ENGLISH GRAY	Obsolete no replacement	No R/B

**Table 4 (a). Replacement Patient Operating Instructions**

<b>Previous Item Number</b>	<b>Description</b>	<b>Change</b>	<b>New Part Number</b>
13138009	POI PORT SERIES ENGLISH	Obsolete replaced by	MN236-C4
13224185	POI PORT SERIES ITALIAN	Obsolete replaced by	MN236-C4
13224193	POI PORT SERIES GERMAN	Obsolete replaced by	MN236-C4
13224206	POI PORT SERIES SPANISH	Obsolete replaced by	MN236-C4
13224214	POI PORT SERIES FRENCH	Obsolete replaced by	MN236-C4
13224222	POI PORT SERIES DUTCH	Obsolete replaced by	MN236-C4
13224231	POI PORT SERIES NORWEGIAN	Obsolete replaced by	MN236-C4
13224249	POI PORT SERIES DANISH	Obsolete replaced by	MN236-C4
13224257	POI PORT SERIES GREEK	Obsolete replaced by	MN236-C4
13224265	POI PORT SERIES PORTUGUESE	Obsolete replaced by	MN236-C4
13224273	POI PORT SERIES SWEDISH	Obsolete replaced by	MN236-C4
13224281	POI PORT SERIES HEBREW	Obsolete replaced by	MN236-C4
13224290	POI PORT SERIES HUNGARIAN	Obsolete replaced by	MN236-C4
13224311	POI PORT SERIES FINNISH	Obsolete replaced by	MN236-C4
13224450	POI PORT SERIES ITALIAN SOL	Obsolete replaced by	MN236-C4
13224468	POI PORT SERIES GERMAN SOL	Obsolete replaced by	MN236-C4
13224476	POI PORT SERIES FRENCH SOL	Obsolete replaced by	MN236-C4
13224484	POI PORT SERIES DUTCH SOL	Obsolete replaced by	MN236-C4
13856473	POI PORT SERIES TURKISH	Obsolete replaced by	MN236-C4
20596377	POI PORT SERIES ENGLISH SOL	Obsolete replaced by	MN236-C4
20705204	POI PORT SERIES SPANISH SOL	Obsolete replaced by	MN236-C4
13138009-12	POI PORT SERIES ARABIC	Keep	13138009-12

**Table 4 (a). Replacement Patient Operating Instructions**

<b>Previous Item Number</b>	<b>Description</b>	<b>Change</b>	<b>New Part Number</b>
14123888	POI HI FLOW G4 ENGLISH	Obsolete replaced by	MN236-C4
14141955	POI HI FLOW G4 DANISH	Obsolete replaced by	MN236-C4
14141963	POI HI FLOW G4 DUTCH	Obsolete replaced by	MN236-C4
14141971	POI HI FLOW G4 DUTCH SOL	Obsolete replaced by	MN236-C4
14141980	POI HI FLOW G4 FINNISH	Obsolete replaced by	MN236-C4
14141998	POI HI FLOW G4 FRENCH	Obsolete replaced by	MN236-C4
14142018	POI HI FLOW G4 FRENCH SOL	Obsolete replaced by	MN236-C4
14142026	POI HI FLOW G4 GERMAN	Obsolete replaced by	MN236-C4
14142042	POI HI FLOW G4 GERMAN SOL	Obsolete replaced by	MN236-C4
14142051	POI HI FLOW G4 GREEK	Obsolete replaced by	MN236-C4
14142771	POI HI FLOW G4 HEBREW	Obsolete replaced by	MN236-C4
14142798	POI HI FLOW G4 HUNGARIAN	Obsolete replaced by	MN236-C4
14142819	POI HI FLOW G4 ITALIAN	Obsolete replaced by	MN236-C4
14142827	POI HI FLOW G4 ITALIAN SOL	Obsolete replaced by	MN236-C4
14142835	POI HI FLOW G4 NORWEGIAN	Obsolete replaced by	MN236-C4
14142894	POI HI FLOW G4 PORTUGUESE	Obsolete replaced by	MN236-C4
14142923	POI HI FLOW G4 SPANISH	Obsolete replaced by	MN236-C4
14142991	POI HI FLOW G4 SWEDISH	Obsolete replaced by	MN236-C4
14143002	POI HI FLOW G4 TURKISH	Obsolete replaced by	MN236-C4

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

[www.cairemedical.com](http://www.cairemedical.com)

Europe: +44 (0) 119-936-7060

[www.cairemedical.com](http://www.cairemedical.com)

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

[www.cairemedical.com](http://www.cairemedical.com)

Europe: +44 (0) 119-936-7060

[www.cairemedical.com](http://www.cairemedical.com)

*Copyright © 2018 CAIRE Inc. CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.*

---