

Topaz Ultra

PSA Oxygen Generator Parts List

(Serial Number 83587 and Higher)



AS123-7, AS123-8, AS123-9, AS123-10 (120 V ~ /60 Hz)			
	Part Number	Description	Quantity
12,000 Hour/ Service Parts*	CO012-11	Compressor Kit, 2660 (or order CO305-1)	1
	CO305-1	Compressor, 120 V/60 Hz, 2660CE50 (or order CO012-11)	1
	MO007-1	Mount Sandwich, 30 lb.	4
30,000 Hour/ Service Parts*	VA034-1	Feed and Waste Solenoid Valve Rebuild Kit, 316372 (without Coil)	4
	VA048-1	Valve Rebuild Kit, 318-356, for Equalization Valve VA003-1	1
Associated Items	BE303-1S	3" x 18" Bed Assembly	
	CB100-1	Circuit Board, Oxygen Monitor (Serial Number 2656763-1 and Lower)	
	CB196-1	Circuit Board, Oxygen Monitor (Serial Number 2656763-1 and Higher)	
	CB200-1S	Main Circuit Board	
	CC012-8	Capacitor, 20 UF 370 V, Aluminum Case	
	CO305-1	Compressor, 120 V/60 Hz, 2660CE50 (Requires FN049-2 (638424) Fan Replacement)	
	CR001-1	Circuit Breaker, 8 amp	
	F0025-1	Fitting, Plastic, Adaptor, B-Size, ARP-128, F-SW x 1/4" Barb	
	F0135-2	Compression Nuts, Brass, Fitting, Ferrule, 3/8"	
	F0164-2	1/8" MPT x 1/8" Tube, Push-Lok Fitting	
	F0292-3	1/4" MPT x 5/16" Tube, Push-Lok Tube Fitting, Plastic Elbow Swivel	
	F0693-1	1/4" Tee Push-Lok Tube Fitting	
	F0694-1	1/4" Check Valve Push-Lok Tube Fitting	
	F0695-1	1/4" MPT x 1/4" Tube, Push-Lok Tube Fitting, Elbow	
	F0697-1	1/4" MPT x 1/4" Tube, Push-Lok Tube Fitting, Straight	
	F0698-1	1/8" MPT x 1/4" Tube, Push-Lok Tube Fitting, Straight	
	F0750-1	1/8" FPT x 1/4" Tube, Push-Lok Tube Fitting, Straight	
	F0882-1	10-32 MPT x 1/8" Push-Lok Tube Fitting, Plastic Elbow Swivel	
	FM080-3	10 LPM Flowmeter, Key	
	FN022-1	Fan, 120 V/60 Hz	
	FN049-2	Fan, Blade Black (638424)	
	FO003-1	Muffler, Foam Acoustic	
	FO015-1	4 1/2 x 4 1/2 x 1/2" Foam Filter, 15 ppi	
	GA001-2	Gauge (0 – 100 psig)	
	GA001-5	Gauge (0 – 30 psig)	
	HM009-1	Hourmeter 100 – 270 V/50/60 Hz	
	HX002-2	Heat Exchanger	
	MO007-1	Mount Sandwich, 30 lb.	
MU034-1	1/4" Muffler		
OF001-16	Orifice, .132 OD x .250, Brass, Slug, .059, Barbell (Requires 1; Equalization)		
OF060-1	Orifice, Inline, Nylon, 1/8", Tube, .003"		
RG088-1	Regulator (0 – 100 psig)		
SW004-1	Rocker Switch (Green)		
TA149-1S	Mixing Tank Assembly		

*Whichever service interval occurs first.

Topaz Ultra

PSA Oxygen Generator Parts List

(Serial Number 83587 and Higher)



AS123-7, AS123-8, AS123-9, AS123-10 (120 V ~ /60 Hz)			
	Part Number	Description	Quantity
Associated Items	TU002-6	¼" Polyurethane Tubing (Black)	
	TU003-2	⅝" Nylon Tubing (Black)	
	TU003-4	¾" Nylon Tubing (Black)	
	TU004-8	⅝" x ⅜" Vinyl Tubing	
	TU006-2	⅝" PVC Reinforced Tubing	
	TU256-1	¼" Polyurethane Tubing (Green)	
	TU256-2	¼" Polyurethane Tubing (Clear)	
	TW001-1	Tie Wrap, 4"	
	TW002-1	Tie Wrap, 14.9"	
	VA001-1	¼" Solenoid Valve Rebuild Kit with Coil, 100 VDC (Feed and Waste)	
	VA003-1	Solenoid Valve, Equalization, Floater (HV238-915-2), 100 VDC	
	VA290-1	Block Valve (non-anodized)	
	WH149-1*	Wire Harness	
	Supplemental Equipment	AN005-1	Oxygen Analyzer (Maxtec Max O ₂ Plus)
AN021-1		Oxygen Analyzer (Maxtec Handi)	
AN075-1		Oxygen Analyzer/Sensor (Maxtec Max O ₂ Plus)	

*When replacing circuit boards to CB016-1, you may require a wire harness (WH149-1).



© 2019 AirSep Corporation. All Rights Reserved.
 AirSep reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. Any rights expressly stated herein are reserved by us, as applicable.

Ordering Information:
 260 Creekside Drive
 Buffalo, NY 14228-2075 U.S.A.
 Tel: (716) 691-0202 ■ Fax: (716) 691-1255
www.airsepcpd.com ■ info@airsep.com