

Centrox

PSA Oxygen Concentrator Parts List

(All Serial Numbers)



AS074-1 (120 V ~ /60 Hz)			
	Part Number	Description	Quantity
8,000 Hour/ Service Parts*	FI153-1	Filter Element for FI137-1	1
12,000 Hour/ Service Parts*	CO185-1	Compressor, 2660CE40 (or order CO220-1)	2
	CO220-1	Compressor Kit, 2650 & 2660 (or order CO185-1)	2
30,000 Hour/ Service Parts*	VA034-1	Feed and Waste Rebuild Kit, 316372	6
	VA048-1	Equalization Valve Rebuild Kit without Coil, 316-375	3
	VA495-1	1/8" Product Solenoid Valve, 100 VDC, SC8256A002V	1
Associated Items	BE126-1	Three 4" Diameter Beds Assembly	
	BO014-6	Bolt, Carb Stl., Plated Hex Slot, 1/4"-20 x 2" with 3M2510	
	CB079-1	Circuit Board (use CB057-1)	
	CB100-1	Circuit Board, Oxygen Monitor	
	CO185-1	Compressor, 115 V/ 60 Hz, 2660CE40	
	CR021-1	Circuit Breaker, 10 amp	
	F0135-2	Compression Nuts, Brass, Fitting, Ferrule, 3/8"	
	F0292-3	1/4" MPT x 5/16" Tube, Push-Lok Tube Fitting, Plastic Elbow Swivel	
	F0571-1	1/4" MPT x 3/8" Tube Brass Fitting, Elbow, with Tubenut (Required 1 ea. for HX010-2)	
	FI137-1	Coalescing Element, Housing and Filter Element	
	FM025-2	15 LPM Flowmeter	
	FN001-1	Fan, 115 V/50/60 Hz	
	FO015-1	4 1/2 x 4 1/2 x 1/2" Foam Filter, 15 ppi	
	FO039-1	1 7/8 x 1/2" Foam Filter, 15 ppi	
	FU009-1	15 Amp Fuse	
	FU014-1	1.0 Amp Fuse	
	GA001-1	Gauge (0 – 15 psig)	
	GA001-2	Gauge (0 – 100 psig)	
	HM009-1	Hourmeter 100 – 270 V/50/60 Hz	
	HX010-2	Heat Exchanger, Coil, 3/8" OD (Required 1 ea. for F0571-1 and F0135-2)	
	MI039-1	Misc., Spacer, .750 x 1.35"	
	MI188-1	Misc., Spring Plate, Top, Molded	
	MU033-1	Muffler Assembly	
	OF060-1	Orifice, Inline, Nylon, 1/8", Tube, .003" (for oxygen sample line)	
	RG022-1	Regulator (0 – 10 psig)	
	RL013-2	Relay Singlepole, 50 Amp, 3-32 VDC, 24-265 VAC, 50/60 Hz	
	SP005-2	Spring, 15 lb., 1.15 x .080	
	SW003-3	Pressure Switch	
	SW124-1	Switch, Contact Block, 800 F Series, 22 MM, NO	
	SW131-2	Switch, Selector, Assembly, Metal, AB, 22 MM, 2 POS, Green with 1 NO, 120 V	

*Whichever service interval occurs first.

Centrox

PSA Oxygen Concentrator Parts List

(All Serial Numbers)



AS074-1 (120 V ~ /60 Hz)

	Part Number	Description	Quantity
Associated Items	SW141-1	Switch, Selector, Assembly, Metal, AB, 22 MM, 2 POS, Blue with 1 NC, 120 V	
	TR027-1	Transformer	
	TU003-1U	¼" Nylon Tubing	
	TU003-2U	⅜" Nylon Tubing (Black)	
	TU004-10	¾" x ½" Vinyl Tubing	
	TU004-4	¼" x ⅛" Vinyl Tubing	
	TU121-4	⅜" Aluminum Tubing, Alloy (Feed)	
	TU121-5	⅜" Aluminum Tubing, Alloy (Product)	
	VA001-1	¼" Solenoid Valve Rebuild Kit with Coil, 100 VDC (Feed and Waste)	
	VA003-11	⅛" Solenoid Valve Kit with Coil, 100 VDC, Equalization (HV238-915-004), 100 V, XAU11X	
	VA011-1	⅜" Solenoid Check Valve	
	VA034-1	Feed and Waste Solenoid Valve Rebuild Kit (316372) without Coil	
	VA044-1	¼" Check Valve	
	VA048-1	Equalization Valve Rebuild Kit without Coil (316-375)	
	VA275-2	Block Valve (Product), 3-Bed, 15 LPM	
	VA276-2	Block Valve (Feed/Waste), 3-Bed, 15 LPM	
	VA453-1	Dump Valve, 100 VDC (HV288220001)	
	VA495-1	⅛" Product Solenoid Valve, 100 VDC (SC8256A002V)	
	WA023-1	Washer, Carbon Stl., Pl., Flat, 1½ x ⅜ x .09"	
	WA025-2	Washer, Neoprene, 40 Duro, 1⅝" x ⅝" x ⅛"	
Supplemental Equipment	AN005-1	Oxygen Analyzer (Maxtec Max O ₂ Plus)	
	AN021-1	Oxygen Analyzer (Maxtec Handi)	
	AN075-1	Oxygen Analyzer/Sensor (Maxtec Max O ₂ Plus)	



© 2019 AirSep Corporation. All Rights Reserved.
 AirSep reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. Any rights expressly stated herein are reserved by us, as applicable.

Ordering Information:
 260 Creekside Drive
 Buffalo, NY 14228-2075 U.S.A.
 Tel: (716) 691-0202 ■ Fax: (716) 691-1255
 www.airsepcpd.com ■ info@airsep.com

AS-SPL-44 A