

# Reliant

## PSA Oxygen Concentrator Parts List

(All Serial Numbers)



AS072-3 (220 V ~ /50 Hz)			
	Part Number	Description	Quantity
8,000 Hour/ Service Parts*	FI153-1	Filter Element for FI137-1	1
12,000 Hour/ Service Parts*	CO185-2	Compressor, 2660CHI46 (or order CO220-1)	1
	CO220-1	Compressor Kit, 2650 & 2660 (or order CO185-2)	1
30,000 Hour/ Service Parts*	VA034-1	Feed and Waste Rebuild Kit, 316372	6
	VA048-1	Equalization Valve Rebuild Kit without Coil, 316-375	3
	VA495-1	1/8" Product Solenoid Valve, 100 VDC, SC8256A002V	1
Associated Items	BE108-1	Bed Assembly, 3" x 16"	
	BO014-6	Bolt, Carb Stl., Plated Hex Slot, 1/4"-20 x 2" with 3M2510	
	CB057-1	Circuit Board	
	CB100-1	Circuit Board, Oxygen Monitor	
	CO185-2	Compressor, 220 V/ 50 Hz, 2660CHI46	
	CR022-1	Circuit Breaker, 6 amp	
	FI137-1	Coalescing Element, Housing and Filter Element	
	FM091-1	8 LPM Flowmeter	
	FN004-2	Fan, 230 V/50/60 Hz	
	FO038-1	4 3/4 x 7 x 1/2" Foam Filter, 15 ppi	
	FO039-1	1 7/8 x 1/2" Foam Filter, 15 ppi	
	FU004-1	3.0 Amp Fuse	
	GA001-1	Gauge (0 – 15 psig)	
	GA001-2	Gauge (0 – 100 psig)	
	HM009-1	Hourmeter 100 – 270 V/50/60 Hz	
	HX010-1	Heat Exchanger, Coil, 3/8" OD (Required 1 ea. for F0135-2)	
	MI039-1	Misc., Spacer, .750 x 1.35"	
	MI188-1	Misc., Spring Plate, Top, Molded	
	MU033-1	Muffler Assembly	
	OF060-1	Orifice, Inline, Nylon, 1/8", Tube, .003" (for oxygen sample line)	
	RG016-1	1/4" Regulator Wilkerson	
	RG022-1	Regulator (0 – 10 psi)	
	RL013-2	Relay Singlepole, 50 Amp, 3-32 VDC, 24-265 VAC, 50/60 Hz	
	SP005-2	Spring, 15 lb., 1.15 x .080	
	SW003-3	Pressure Switch	
	SW005-1	Rocker Switch (Amber)	
	SW102-2	Rocker Switch (Green)	
	TR027-1	Transformer	
	TU003-1U	1/4" Nylon Tubing	

\*Whichever service interval occurs first.

# Reliant

## PSA Oxygen Concentrator Parts List

(All Serial Numbers)



### AS072-3 (220 V ~ /50 Hz)

	Part Number	Description	Quantity
	TU003-2U	3/8" Nylon Tubing (Black)	
	TU004-10U	3/4" x 1/2" Vinyl Tubing	
	TU004-4U	1/4" x 1/8" Vinyl Tubing	
	TU121-1	3/8" Aluminum Tubing, Alloy (Intake)	
	TU121-2	3/8" Aluminum Tubing, Alloy (Feed)	
	TU121-3	1/4" Aluminum Tubing, Alloy (Product)	
	VA001-1	1/4" Solenoid Valve Rebuild Kit with Coil, 100 VDC (Feed and Waste)	
	VA003-11	1/8" Solenoid Valve Kit with Coil, 100 VDC, Equalization (HV238-915-004), 100 V, XAU1IX	
<b>Associated Items</b>	VA034-1	Feed and Waste Solenoid Valve Rebuild Kit (316372) without Coil	
	VA044-1	1/4" Check Valve	
	VA048-1	Equalization Valve Rebuild Kit without Coil (316-375)	
	VA275-1	Block Valve (Product), 3-Bed	
	VA276-1	Block Valve (Feed/Waste), 3-Bed	
	VA495-1	1/8" Product Solenoid Valve, 100 VDC (SC8256A002V)	
	VA566-1	3/8" Check Valve	
	WA023-1	Washer, Carbon Stl., Pl., Flat, 1 1/2 x 3/32 x .09"	
	WA025-2	Washer, Neoprene, 40 Duro, 1 5/8" x 3/8" x 1/8"	
<b>Supplemental Equipment</b>	AN005-1	Oxygen Analyzer (Maxtec Max O <sub>2</sub> Plus)	
	AN021-1	Oxygen Analyzer (Maxtec Handi)	
	AN075-1	Oxygen Analyzer/Sensor (Maxtec Max O <sub>2</sub> Plus)	



© 2019 AirSep Corporation. All Rights Reserved.  
 AirSep reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. Any rights expressly stated herein are reserved by us, as applicable.

**Ordering Information:**  
 260 Creekside Drive  
 Buffalo, NY 14228-2075 U.S.A.  
 Tel: (716) 691-0202 ■ Fax: (716) 691-1255  
[www.airsepcpd.com](http://www.airsepcpd.com) ■ [info@airsep.com](mailto:info@airsep.com)

AS-SPL-47 A