

Onyx Plus

PSA Oxygen Generator Parts List



AS017-1 (120 V ~ /60 Hz)			
	Part Number	Description	Quantity
12,000 Hour/ Service Parts*	CO012-11	Compressor Kit, 2660 (or order CO240-1)	1
	CO240-9	Compressor, 120 V/60 Hz, 2660CE36 (or order CO012-11)	1
30,000 Hour/ Service Parts*	VA003-9	Valve for Equalization, Plastic, 100 VDC	
	VA034-1	Feed and Waste Solenoid Valve Rebuild Kit, 316372 (without Coil)	
Associated Items	BE157-1S	8L Bed Assembly	
	BO014-6	Bolt, Carb Stl., Plated Hex Slot, 1/4"-20 x 2" with 3M2510	
	CA002-11	Back Panel Cabinet	
	CA005-11	Front Panel Cabinet	
	CA006-11	Side Panel Cabinet	
	N/A	Main Circuit Board (NOT AVAILABLE)	
	CC012-3	Capacitor, 15UF 370 V	
	CO240-9	Compressor, 120 V/60 Hz, 2660CE36	
	CR001-1	Circuit Breaker, 8 amp	
	CS013-1	Caster, Stem	
	CS014-1	Caster, Wheel	
	F0025-1	Fitting, Plastic, Adaptor, B-Size, ARP-128, F-SW x 1/4" Barb	
	F0067-1	Fitting, Brass, Compression Ferrule, Nut & Sleeve, 1/4" OD	
	F0134-1	Fitting, Brass, Compression Ferrule, Nut & Sleeve, 3/8" FPT for TU198-064	
	FI002-1	5 1/4 x 3 7/8 x 7/8" Foam Filter, 15 ppi	
	FM080-2	8 LPM Flowmeter, Key	
	FN022-1	Fan, 120 V	
	FO008-1	Muffler, Foam Acoustic, Front Panel	
	FO101-1	Muffler, Foam Acoustic, Side Panel	
	HM009-2	Hourmeter 100 – 270 V/50/60 Hz, with Bracket	
	HX002-2	Heat Exchanger	
	MI039-1	Misc., Spacer, .750 x 1.35"	
	MI161-3	Misc., Felt, 1.9 x .7", F10, 1/2" LG	
	MI188-1	Misc., Spring Plate, Top, Molded	
	MU043-2S	Muffler, Half-Felt, Assembly	
	OF001-2	Orifice, .132 OD x .250, Brass, Slug, .049, Barbell (Requires 1 Equalization)	
	OF049-1	Orifice, .132 OD x .039, Barbell (Requires 2 Purge - Beds)	
	PA001-1	Packaging, Poly Bag, Clear, 25 x 36"	
	PA002-2	Packaging, Kraft Carton	
	PA054-2	Packaging, Polyam Foam Block for Flowmeter	
	PA274-1	Packaging, Bottom Insert	
	PA275-1	Packaging, Top Left Insert	
PA276-1	Packaging, Top Right Insert		
RG088-1	Regulator, 1/4" FPT, 100 psig		
RN023-1	Standard Resonator without Tube		

*Whichever service interval occurs first.

Onyx Plus

PSA Oxygen Generator Parts List



AS017-1 (120 V ~ /60 Hz)				
	Part Number	Description	Quantity	
Associated Items	SP005-1	Spring, 10 lb., 1.14 x .99, Mounts		
	SW002-3	Temperature Switch, 165°F		
	SW114-1	Rocker Switch, 2-Pole		
	TA132-1	Mixing Tank Assembly		
	TU054-1	¼" OD Vinyl Tubing (Black)		
	TU055-1	¼" OD Vinyl Tubing (Green)		
	TU198-064	⅝" OD Fep Tubing (Clear), .04" Wall, 6.5" LG (Bottom of Beds) (Requires two (2) F0134-1 ferrule compression nuts)		
	TW001-4	Tie Wrap, 4"		
	TW002-1	Tie Wrap, 14.9" (Beds)		
	TW004-1	Tie Wrap, 14.5" (Equalization Valve/Capacitor/Resonator)		
	VA001-1	¼" Solenoid Valve Rebuild Kit with Coil, 100 VDC (Feed and Waste)		
	VA003-9	Valve for Equalization, Plastic, 100 VDC		
	VA006-1	¼" Check Valve		
	VA290-3	Block Valve (non-anodized)		
	WA023-1	Washer, Carbon Stl., Pl., Flat, 1½ x ⅜ x .09"		
	WA025-2	Washer, Neoprene, 40 Duro, 1⅝" x ⅝" x ⅛"		
	Supplemental Equipment	AN005-1	Oxygen Analyzer (Maxtec Max O ₂ Plus)	
		AN021-1	Oxygen Analyzer (Maxtec Handi)	
AN075-1		Oxygen Analyzer/Sensor (Maxtec Max O ₂ Plus)		



© 2019 AirSep Corporation. All Rights Reserved.
 AirSep reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. Any rights expressly stated herein are reserved by us, as applicable.

Ordering Information:
 260 Creekside Drive
 Buffalo, NY 14228-2075 U.S.A.
 Tel: (716) 691-0202 ■ Fax: (716) 691-1255
www.airsep.com ■ cpd@caireinc.com

AS-SPL-51 A