

Portable and Stationary Concentrator Filters

CAIRE® Companion 5™

Intake Filter

Item Number 15062925



AirSep® VisionAire™

Handle Filter (optional)

Item Number FI202-1



AirSep® NewLife® Family

Gross Particle Filter

Item Number FI002-1



Felt Filter

Item Numbers

Elite 110V MI161-1
 Elite 220V MI161-2
 Intensity 8 MI161-3
 Intensity 10 MI161-2



SeQual® Eclipse®

Preventative Maintenance Kit Item Number 5022-SEQ

Includes:

2 Air Intake Filters (7028-SEQ)
 Compressor Intake Filter (8069-SEQ)
 HEPA Filter, External (6986-SEQ)*
 HEPA Filter, Internal (4770-SEQ)*
 9V Battery (8098-SEQ)
 O-Rings for Internal HEPA Filter
 (9251-SEQ & 4382-SEQ)
 10 Zip Ties

**Contact technical service if you are unsure of which HEPA filter you need.*



SeQual® eQuinox®

Air Inlet Filter

Item Number 4845-SEQ



Oxygen Outlet Filter

Item Number 4770-SEQ



Compressor Inlet Filter

Item Number 9759-SEQ



AirSep® FreeStyle™ 3

Intake Filter

Item Number FI144-1



*Compatible with Internal and user-replaceable battery versions

AirSep® FreeStyle™ 5

Intake Filter

Item Number FI194-1



*Compatible with Internal and user-replaceable battery versions