

NEW...

AirSep External Power Cartridge and Charger

Compatible with AirSep Focus™, FreeStyle™, and FreeStyle™ 5
 Portable Oxygen Concentrators

The new AirSep External Power Cartridge provides extended battery life for all AirSep POCs.

- **NEW LIGHTWEIGHT** alternative to the AirSep AirBelt
- External charger can charge 2 batteries simultaneously
- Built-in belt clip for easy wearability
- Fully charges battery in 5 hours or less
- Easy-access battery gauge displays remaining charge

Battery Durations:

Focus

Setting	External Power Cartridge Duration
2	3.5 Hours

Focus adapter cable (CD035-1) must be purchased separately to use the external power cartridge with the Focus.

FreeStyle

Setting	Duration of FreeStyle Internal Battery + External Power Cartridge
1	10 Hours
2	6 Hours
3	5 Hours

FreeStyle 5

Setting	Duration of FreeStyle 5 Internal Battery + External Power Cartridge
1	7 Hours
2	4 Hours
3	3 Hours
4	2.5 Hours
5	2 Hours



Charger contains an easy-to-read light system:

Flashing Green Light: Charging

Solid Green Light: Fully Charged

Flashing Red Light: Error

Solid Red Light: No battery inserted

ORDERING INFORMATION

- BT024-1** AirSep External Power Cartridge
- PW029-1** Charger (Dock) and Power Supply
- BT026-1** External Power Cartridge & Charger Combination (US)*
- BT026-2** External Power Cartridge & Charger Combination (EUROPEAN)*
- BT026-3** External Power Cartridge & Charger Combination (AUSTRALIAN)*
- BT026-4** External Power Cartridge & Charger Combination (UK)*

Additional batteries may be purchased independently.

*Combinations include one battery, one charger, and one power supply.