

# FREESTYLE<sup>®</sup>

Portable Oxygen Concentrator

## Discover your independence

The lightweight and wearable FreeStyle<sup>®</sup> offers you the ultimate freedom to go anywhere. Easy to use while active in your everyday activities and FAA-approved for use on a commercial flight — it is perfect for everywhere you go.



### A LIGHTWEIGHT SOLUTION

Replace your tanks with the lightweight FreeStyle<sup>®</sup> at under 5 lbs. This device can be easily carried by handle, worn by a comfortable shoulder strap, or attached to a backpack. You really can go anywhere and enhance your activity level.



### PEACE OF MIND

Enjoy the convenience of a portable oxygen concentrator with sensitive breath detection — UltraSense<sup>™</sup> Technology — which ensures the device detects each of your breaths quickly to give you the oxygen you need each time you breathe in.



### STAY ACTIVE

Quickly swap the battery packs on the go so you never have to worry about running out of oxygen. The convenient pack will re-charge anytime wall (AC) or motor vehicle (DC) power is connected and is easily accessible from the front of the device.



#### CAIRE Contact Information:

Providers Toll Free  
1.800.482.2473  
Patients Toll Free  
1.877.704.0878

[www.cairemedical.com](http://www.cairemedical.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

# FREESTYLE®

Portable Oxygen Concentrator

## Includes

- Portable Oxygen Concentrator
- Carry Bag
- AC/DC Universal Power Supply and Cords
- Carry-All Accessory Bag
- Shoulder Strap

## Product Features

- Under 5 pounds with a small compact frame
- Replaceable battery pack with up to 3½ hour battery life at setting 2
- Easy to carry with a handle, or wear with a shoulder strap or backpack
- FAA-approved for commercial airline use
- Fully functional on wall (AC) or motor vehicle (DC) power



## Specifications

**FLOW RATE**  
Pulse 1-3

**DIMENSIONS**  
8.6 x 6.1 x 3.6 in  
(21.8 x 15.5 x 9.1 cm)

**WEIGHT**  
4.9 lbs (2.2 kg)

**BATTERY DURATION**  
Setting 1 - 6.5 hours (6:30)  
Setting 2 - 3.5 hours (3:30)  
Setting 3 - 3.0 hours (3:00)

**BATTERY RECHARGE**  
4 – 5 hours

**WARRANTY**  
3 years FreeStyle  
1 year battery and power supplies

## Available Accessories



**CARRY BAG**



**AC/DC POWER SUPPLY**  
(CORDS AVAILABLE SEPARATELY)



**CARRY-ALL ACCESSORY BAG**



**BACKPACK HARNESS**



**EXTRA BATTERY PACK**



**DESKTOP CHARGER**



**SPARE INTAKE FILTER**

### CAIRE Contact Information:

Providers Toll Free  
1.800.482.2473  
Patients Toll Free  
1.877.704.0878

[www.cairemedical.com](http://www.cairemedical.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

ML-CONC0122 D