

autoSAT® Technology

autoSAT™ Technology adjusts the Eclipse or eQuinox to keep up with your changing breath rates. Your breath rate changes as you go through tasks of every-day life, such as sitting, standing, or walking. To receive the best oxygen therapy on pulse mode, you need a device that will react and adapt to these changes in your breath rate. autoSAT recognizes your changing breath rates and speeds up the oxygen production to keep up with your activity level. Much like cruise control adjusts a car's power to ensure a constant speed, autoSAT adjusts oxygen production to guarantee that you are delivered your prescribed amount of oxygen with each breath. The unparalleled performance of autoSAT ensures you receive the oxygen you need, with every breath, no matter how active you are.



UltraSense™ Technology



FreeStyle, FreeStyle 5, and Focus are equipped with UltraSense™ technology, which provides a more sensitive pulse dose trigger than any other device on the market. By detecting your breath with less effort, UltraSense helps to ensure that a pulse of oxygen is delivered with each breath you take.

The increased sensitivity provided by UltraSense also ensures the pulse is delivered in the first stage of your breath, when it is most critical to receive oxygen. UltraSense increases your comfort level and allows a wider variety of people to use pulse dose oxygen effectively. With UltraSense, you always receive the oxygen you need, when you need it most.

CAIRE's POCs are available for online purchase, as well as through short- or long-term rentals by your local oxygen provider. All CAIRE respiratory products are covered under a variety of payment or insurance options with full or partial reimbursement. To receive a free informational packet or more information on the purchase or rental of CAIRE products, please contact your local provider.

Regain independence and feel the difference a lighter oxygen system can make for ambulation.

Accessories

SeQual eQuinox

- Accessory Bag
- eQuinox Cart
- AC Power Supply with Power Cord
- DC Power Cord with no transformer brick
- 12-Cell Battery
- Optional Accessories:**
- 24-cell Battery
- Desktop Charger
- Wheelchair Strap
- Fingertip Pulse Oximeter

SeQual Eclipse 5

- AC and DC Power Supply
- Accessory Bag
- Power Cartridge (Battery)
- Eclipse Cart
- Optional Accessories:**
- Desktop Charger
- Eclipse Protective Cover
- Eclipse Wheelchair Pack
- Wheelchair Strap
- Fingertip Pulse Oximeter
- Humidifier Adapter Kit
- Travel Accessory Kit
- Travel Case
- Universal Power Plug Set

AirSep FreeStyle 5

- AC Power Supply with Power Cord
- Carry-All-Bag for Accessories
- Carrying Case with Handle and Shoulder Strap
- DC Power Supply
- Optional Accessories:**
- External Rechargeable Battery
- AirBelt
- Backpack Harness
- Wheelchair Strap

AirSep FreeStyle

- Carry-All-Bag for Accessories
- Carrying Case and Shoulder Strap
- Universal AC/DC Power Supply with Power Cords
- Optional Accessories:**
- External Rechargeable Battery
- AirBelt
- Backpack Harness
- Wheelchair Strap

AirSep Focus

- Carry-All-Bag for Accessories
- Carrying Bag with Handle and Two Battery Pockets
- Carrying Case
- Coil Cord to Charge Battery
- Coil Cord to Power Focus from Battery
- 2 Batteries
- Battery Case with Belt Clip
- Universal AC/DC Power Supply with Power Cords
- Optional Accessories:**
- External Rechargeable Battery
- AirBelt and Battery Belt Interface Cable

All POCs are FAA approved for in-flight use.

CAIRE, Inc.
2200 Airport Industrial Way,
Suite 500 Ball Ground, GA 30107
800.482.2473 Telephone
888-932-2473 Fax

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham,
Berkshire, RG41 2PL
Ph +44(0) 1189 367060
Fax +44(0) 1189 799245

www.CAIREmedical.com

Copyright © 2014 Chart Industries
Chart Industries reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or, processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

PN ML-CONC0002 Rev. G

CAIRE
A Chart Industries Company
AirSep® / SeQual®

CHART

Discover **your**
perfect **fit**
from CAIRE.

Choose from palm-sized portables and continuous flow transportable oxygen concentrators to meet your needs.



**Eclipse 5™ & eQuinox™ with autoSAT™
FreeStyle™, FreeStyle™ 5 & Focus™ with UltraSense™**

Find the perfect fit for YOUR lifestyle.

We understand that finding an oxygen device that accommodates both lifestyle and oxygen needs can be a difficult process for both you and your loved ones. CAIRE offers the largest variety of oxygen equipment on the market, providing you with the option to pick the right device for both your lifestyle and your needs. Whether you need an ultra lightweight, wearable portable oxygen concentrator (POC), or a robust POC that offers pulse and continuous flow modes, CAIRE has the perfect match for you.

You have an oxygen prescription and ...

want a 24/7, continuous and pulse flow mode transportable POC that ...

need a wearable, pulse-only device and ...

need up to a 5-pulse setting.

AirSep FreeStyle™ 5 with UltraSense™

The FreeStyle™ 5 is the lightest POC offering five pulse setting options, providing the highest oxygen production of any wearable POC. Diverse power choices allow you to adapt the FreeStyle 5 by using the rechargeable internal battery, the AC or the DC power supplies, and the optional AirBelt or External Battery Cartridge for extended time away from a power source.

- Weighs 6.2 lb (2.8 kg)
- Pulse flow setting up to 5
- Internal battery duration of up to 3 hours
- AirBelt extends battery life by up to 1 hour at setting 5
- External rechargeable battery provides an additional 2 hours at setting 5



is robust and reliable with a newly improved package.

SeQual Eclipse™ 5 with autoSAT®

The Eclipse™ is a proven, reliable POC providing continuous and pulse flow options for use at home, on the go, and in the car. New improvements make the Eclipse 5 more robust and reliable, as well as make it fully functional on DC power at flow rates up to 3.0 LPM and pulse settings 1-9. The battery cartridge will re-charge in the car while running at 2.0 LPM or less.

- Weighs 18.4 lbs (8.3 kg) with battery installed
- Continuous flow of 0.5 - 3.0 LPM
- 9 pulse settings up to 192mL—the highest pulse dose output of any POC
- Battery duration: 2 hours (at 2LPM continuous flow); 5.1 hours (pulse at a 2 setting)



is lightweight and has cutting-edge technology.

SeQual eQuinox™ with autoSAT®

The eQuinox™ features an advanced voice module interface as well as the highest pulse and continuous flow options available in an easy-to-manuever frame. The eQuinox is 35% smaller than the Eclipse 5, making it the smallest transportable concentrator on the market offering 3.0 LPM flow. The eQuinox has a variety of new features to make using a POC easier than ever. A multi-language voice interface audibly alerts you of setting changes, alarms, and battery time remaining. The screen intermittently displays the remaining battery time, providing an extra layer of comfort and peace of mind. The eQuinox is fully functional on AC, DC, and battery power with pulse settings 1-9 and continuous flow up to 3.0 LPM. All this, and it is only 45 dB at 2 continuous flow—making it world class quiet.

- Weighs 14 lbs (6.4 kg) with standard battery
- Continuous flow of 0.5 - 3.0 LPM
- 9 pulse settings up to 192 mL—the highest pulse dose output of any POC
- Standard (12 cell) Battery Duration: 1.2 hours at 2.0LPM continuous flow, 2.7 hours at pulse setting 2
- Optional (24 cell) Battery Duration: 2.75 hours at 2.0LPM continuous flow, 5.9 hours at pulse setting 2



require moderate flow.

AirSep FreeStyle™ with UltraSense™

The FreeStyle™ offers a uniquely efficient, ultra-compact and flexible pulse-flow delivery system with three settings. At only 4.4 lbs and in an attractive package, the FreeStyle is easily carried or worn with a backpack or shoulder strap. Running on AC, DC, or internal battery power, the FreeStyle provides the oxygen you need without hampering your mobility. The AirBelt and external power cartridge provide extended durations away from a power source.

- Weighs 4.4 lbs (2.0 kg)
- Pulse flow setting up to 3
- Internal battery duration of up to 3.5 hours
- AirBelt extends battery life by up to 3 hours at setting 3
- External rechargeable battery provides an additional 3 hours at setting 3



want the lightest POC available.

AirSep Focus™ with UltraSense™

Focus™ is the lightest POC on the market, small enough to fit in the palm of your hand and incredibly wearable using a belt clip or shoulder strap. Simply plug the Focus in to AC/DC power supply or connect to one of the battery options for unlimited oxygen on the go.

- Weighs 1.75 lbs (0.8 kg)
- Pulse flow setting of 2
- Battery duration of 3 hours (1.5 hours per battery; 2 included)
- AirBelt extends battery life by 4 hours
- External rechargeable battery provides an additional 3.5 hours of use at setting 2

