

Does your POC have **UltraSense™?**

Exclusively available for AirSep POCs!

AirSep FreeStyle™

AirSep FreeStyle™ 5


AirSep Focus™

Contact **Customer Service** for more information
customerservice.usa@chartindustries.com
or 800-482-2473
www.CAIREmedical.com



UltraSense™: Sensitivity Does Matter

Sensitivity values have been arranged on the graph below to demonstrate the amount of work, or effort, required by the patient to trigger a pulse dosage from each of the devices included.



Least Sensitive

Most Sensitive

Average Trigger Sensitivity Values

FreeStyle 5	-0.054
FreeStyle	-0.084
Focus	-0.114
eQuinox	-0.135
Eclipse 5	-0.135
Inogen One	-0.167
Simply Go	-0.195
Simply Go - sleep mode	-0.196
OxLife	-0.203
iGo	-0.238
EasyPulse	-0.360
ActivOx	-0.361
SOLO2	-0.389

*All values are given as negative cm of H₂O (negative pressure)

*Values for each device's trigger sensitivity are an average of the trigger sensitivity at each setting at 20 BPM.