FOCUS[™] Portable Oxygen Concentrator

Our smallest and lightest

Less than 2 pounds (0.9 kg) and small enough to fit in the palm of your hand, the Focus™ is easily worn with a shoulder strap or belt clip. It is also FAA-approved for commercial air travel.





STAY ACTIVE

Enjoy your freedom and stay active thanks to this lightweight device equipped with 3 hours of standard battery life. The batteries store conveniently in the pouches located on either side of the carrying bag.

\¢ |

PEACE OF MIND

Enjoy convenient, portable oxygen with sensitive breath detection — UltraSense™ Technology — which ensures the quick delivery of pulse oxygen in the early stages of your breath and reduces the risk of skipped breaths.

LIGHTWEIGHT AND DISCREET

Equipped with pulse setting 2, the Focus[™] is designed to meet the needs of oxygen users without adding the extra weight of higher flow units. Discreet in size, it is perfect for helping you maintain your active lifestyle whether it is playing a round of golf or going out for a night of dinner and dancing.



CAIRE Contact Information: Providers Toll Free 1.800.482.2473 Patients Toll Free 1.877.704.0878 www.cairemedical.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

FOCUS™

Portable Oxygen Concentrator

Includes

- Portable Concentrator
- 2 Batteries
- AC/DC Universal Power Supply and Cords
- Carry Bag w/ Battery Pockets

Product Features

- Only 1.75 lbs (0.8 kg)
- Simple to use no settings to adjust and one touch operation
- No patient maintenance
- FAA-approved for commercial airline use

Specifications

FLOW RATE Pulse 2

SOUND LEVEL 45 dB(A) @ Setting 2

BATTERY DURATION

3 hrs (1.5 per battery -

DIMENSIONS 6.4 x 4.8 x 2.5 in (16.4 x 12.2 x 6.1 cm)

WEIGHT 1.75 lbs (0.8 kg) 2 included) BATTERY RECHARGE 4 hours

Available Accessories



EXTRA BATTERIES





CARRYING BAG*

Battery Charging and Operation

Carry-All Accessory Bag

• Wearable via belt clip, carry handle, or shoulder strap

• Fully functional on wall (AC)

or motor vehicle (DC) power

WARRANTY 3 vears Focus™

1 year battery and

power supplies

• Can swap and exchange batteries for long durations

Coiled Cords

Shoulder Strap

Battery Case



BAG W/BATTERY POCKETS*

*Items sold together as a bag kit.



AC/DC POWER SUPPLY (CORDS AVAILABLE SEPARATELY)



BATTERY CHARGING CABLE

CAIRE Contact Information: Providers Toll Free 1.800.482.2473 Patients Toll Free 1.877.704.0878

www.cairemedical.com



BATTERY-TO-UNIT CABLE



CARRY-ALL ACCESSORY BAG

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

