

# NEWLIFE® ELITE®

5 LPM At-Home Oxygen Concentrator

## THE ONE THAT SET THE STANDARDS

Performing in the harshest of environments, this 5 LPM concentrator is simple to service with time-proven quality and reliability.



### ROBUST AND EASY TO SERVICE

The robust Elite is highly-rated for efficiency and performance in harsh environments, extreme temperatures, tropical climates, and developing countries where service capabilities are limited. It is also a perfect choice for long term care facilities.

The simplistic internal design makes service and maintenance simple for the Elite, along with the entire NewLife series line of concentrators, making it a favorite of technicians everywhere.



### COMFORT AND EASE OF USE

Enjoy peace of mind and reliable oxygen therapy delivery with the Elite — your at-home oxygen solution. A simple user interface and an easy-to-use humidification attachment take the guesswork away and allow you to enjoy the little things.

Rest easy, stream your favorite TV show, cheer on your sports team, or spend time with loved ones, and never worry about your oxygen delivery.



CAIRE Contact Information:  
1.800.482.2473  
[www.caireinc.com](http://www.caireinc.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.  
© 2021 CAIRE Inc. All Rights Reserved.

# NEWLIFE® ELITE®

5 LPM At-Home Oxygen Concentrator

## Includes

- Oxygen Concentrator
- Humidifier Adaptor Tube

## Product Features

- Highly rated in harsh environments
- Lockable flow meter
- Time-proven NewLife® ruggedness and reliability
- Easy to service and maintain

## Specifications

### FLOW RATE

1.0 - 5.0 LPM

### DIMENSIONS

28.5 x 15.7 x 14.5 in  
(72.4 x 40.0 x 36.8 cm)

### WEIGHT

54 lbs (24.5 kg)

### SOUND LEVEL

53 dB(A)

### O<sub>2</sub> CONCENTRATION

90 % (+5.5 / -3 %)

### OPERATING

#### ENVIRONMENT

41° to 104°F (5° to 40°C)  
15-90% humidity

### STORAGE

#### ENVIRONMENT

-13°F to 158°F (-25°C to  
70°C)  
0-90% humidity

### MAXIMUM TUBING

50 ft (15.2 m)

### MAINTENANCE

#### SCHEDULE

Intake Filter – Clean  
Weekly

Felt Filter (excluding models  
AS005-4, AS005-241, and  
AS005-263) –

1 year replacement

### WARRANTY

3 Years



## Available Accessories



### HUMIDIFIER ADAPTOR & BOTTLE

Humidifier Adaptor (9 in/23 cm)

Ask your doctor or therapist if  
humidification is right for you.



### NASAL CANNULAS

Various lengths are available.

(7 ft and 25 ft)

Talk to our Customer Service team  
to learn more.



### OXYGEN TUBING

Comes in: 7 ft (2.1 m), 25 ft (7.6 m),  
and 50 ft (15.2 m).



CAIRE Contact Information:

1.800.482.2473

[www.caireinc.com](http://www.caireinc.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.  
© 2021 CAIRE Inc. All Rights Reserved.

ML-CONC0085 H