

CONCENTRATOR FILTERS

AIRSEP® FREESTYLE® 3

Intake Filter


Item Number FI144-1



AIRSEP® FREESTYLE® 5

Intake Filter

Item Number FI194-1



CAIRE® COMPANION 5™

Intake Filter

Item Number 15062925



AIRSEP® VISIONAIRE™

Handle Filter (Optional)

Item Number FI202-1




SEQUAL® ECLIPSE® PREVENTATIVE MAINTENANCE KIT

Item Number 5022-SEQ

Includes:

- 2 Air Intake Filters (7028-SEQ)
- Compressor Intake Filter (8069-SEQ)
- HEPA Filter, Internal (21360366)
- 9V Battery (BT002-1)
- O-Rings for Internal HEPA Filter (9251-SEQ AND 4382-SEQ)
- 10 Zip Ties



AIRSEP® NEWLIFE® FAMILY

Gross Particle Filter

Item Number FI002-1



Felt Filter

Item Numbers

- Elite 110V MI161-1
- Elite 220V MI161-2
- Intensity 8 MI161-3
- Intensity 10 MI161-2



SEQUAL® EQUINOX™ PREVENTATIVE MAINTENANCE KIT

Item Number 21079406

Includes:

- Air Intake Filter (4845-SEQ)
- Oxygen Outlet Filter (21360366)
- O-Rings for HEPA Filter (9251-SEQ and 9250-SEQ)
- Compressor Inlet Filter (9759-SEQ)
- 9V Battery (BT002-1)
- 10 Zip Ties



CAIRE Contact Information:

Providers Toll Free
1.800.482.2473
Patients Toll Free
1.877.704.0878
www.caireinc.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.