

# CAIRE O2 Hospital Discharge Program

**CAIRE**<sup>®</sup>  
O<sub>2</sub> is what we do.

## OFFERING YOUR PATIENTS THE CONVENIENCE OF PORTABLE OXYGEN THERAPY

The CAIRE oxygen hospital discharge program creates an efficient discharge process with a clinically-advanced transportable oxygen concentrator, providing the best service and care for your oxygen therapy patients.

### ✓ CLINICALLY-EFFECTIVE OXYGEN DELIVERY

With 24/7 capabilities, the award-winning Eclipse 5<sup>®</sup> transportable oxygen concentrator offers smart technology that tailors the delivery of oxygen to the individual patient and their disease state to ensure saturation

- ▶ autoSAT<sup>®</sup> Technology automatically adjusts the delivery of oxygen to the patient's respiratory rate
- ▶ UltraSense<sup>®</sup> with adjustable sensitivity breath detection and adjustable rise time ensures consistent timing of oxygen delivery
- ▶ Offers continuous flow up to 3 LPM and pulse flow up to a 192 mL bolus, meeting a wide variety of patient needs

### ✓ QUICK AND EFFICIENT DISCHARGE PROCESS

- ▶ Eclipse 5 is easy to train on and operate, affording unlimited ambulation; a rechargeable battery and motor vehicle and electrical power options encourage the patient to use the device as prescribed
- ▶ Reduce delays for patient discharges and conveniently schedule low-contact in-home setup
- ▶ Facilitate improved service during unique discharge situations, including large rural service areas, inclement weather or times of high volume oxygen setups

### ✓ IMPROVED COMPLIANCE AND CONTINUOUS CARE

- ▶ Optimize patient care by pairing the FreeStyle<sup>®</sup> Comfort<sup>®</sup> with one of CAIRE's high-quality stationary concentrators, covering a broad range of prescriptions and at-home needs, while encouraging mobility
- ▶ Secure visibility to device performance and ability to remotely troubleshoot with myCAIRE telehealth solution
- ▶ CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds



Find out more about CAIRE's Hospital Discharge Program and accessory products at

[CAIREinc.com](http://CAIREinc.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2021 CAIRE Inc. All Rights Reserved. ML-CONC0112 C