# Eclipse 5<sup>®</sup>

# Designed with your active patients in mind



### Key features and benefits.



#### **Unmatched Versatility**

Setting the standard in transportable oxygen with the broadest oxygen delivery capabilities in continuous and pulse dose. In a compact, transportable design, your patients can go anywhere using the easyto-attach traveling cart with a telescoping handle and durable wheels. The extended battery life, AC/DC power options and broad portfolio of accessories help to keep your patients active.



#### **Clinically Effective O2**

Designed to adapt to your patients needs with unique clinical features including autoSAT®, adjustable rise time, UltraSense<sup>TM</sup> and myCAIRE<sup>TM</sup> to help ensure your patients will get the oxygen they need, when they need it.



#### **Proper O2 Saturation**

With 6 settings up to 3 LPM in continuous flow and 9 settings up to a 192 mL bolus in pulse mode, the Eclipse 5 serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



#### **Proven Performance**

The Eclipse 5 utilizes
Vacuum Pressure Swing
Adsorption (VPSA)
powerful air separation
technologythat efficiently
eliminates residual
nitrogen while delivering
the highest purity of
oxygen. The robust
design features allow for
minimum maintenance
and long life, providing a
device you can rely on.

## Specifications.

#### **FLOW RATE**

Continuous Flow: 6 settings up to 3 LPM

Pulse Flow: 9 settings up to 192 mL bolus

#### **DIMENSIONS**

19.3 x 12.3 x 7.1 in (49.0 x 31.2 x 18.0 cm)

#### **WEIGHT**

15 lbs (without battery) 18.4 lbs (with battery)

#### **SOUND LEVEL**

40.0 dB(A) @ Setting 3 Pulse Mode

#### TRIGGER SENSITIVITY

-0.135-0.37 cm/H2O (adjustable)

#### **BATTERY DURATION**

up to 5.4 hours in pulse mode

up to 4.4 hours in continuous flow mode

#### **POWER**

AC Power: 100-240 VAC, 50-60Hz DC Power: 12 VDC minimum (all settings available in AC and DC power supply) Battery: Lithium Ion

#### **WARRANTY**

3 years (Includes sieve beds) I year on battery and power supplies



# Eclipse 5<sup>®</sup>

## Takes you anywhere you want to go





#### **Item Number**

US: 20852326S EU: 20852325S UK: 20852324S



DC POWER SUPPLY (with vehicle accessory outlet plug)

**Item Number** 5942-SEQ



#### POWER CARTRIDGE (BATTERY)

**Item Number** 7082-SEQ



## DESKTOP BATTERY CHARGER (Battery Not Included)

**Item Number** 7112-SEQ



#### TRAVEL ADAPTER

**Item Number** 5591-SEQ



See https://www.iec.ch/world-plugs for details on plug types.

Find out more about Eclipse 5 and accessory products at

UNIVERSAL CART (with telescoping handle)

**Item Number** 5991-SEO



HUMIDIFIER ADAPTER KIT
(Humidifier bottle not included)

**Item Number** 7116-SEQ



ACCESSORY BAG

**Item Number** 7104-SEO



TRAVEL CASE

**Item Number** 5010-SEQ.



TRAVEL ACCESSORY KIT (Power Cartridge, Desktop

Charger, Travel Case)

Item Number 5093-SEQ

**PROTECTIVE COVER** 

Item Number 5052-SEQ



WHEELCHAIR PACK

Item Number 5220-SEQ



SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.