FreeStyle® Comfort® The perfect fit



Key features and benefits



Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.



Smart O2 delivery

Unique clinical features including UltraSense®, autoDOSE, and myCAIRE™ help to ensure your patients will get the oxygen they need, when they need it.



Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds.

Specifications



FLOW RATE

Pulse Dose 1-5
O2 Output Settings:
Setting 1-210 mL/min
Setting 2-420 mL/min
Setting 3-630 mL/min
Setting 4-840 mL/min
Setting 5-1050 mL/min

DIMENSIONS

10.0 x 7.3 x 3.1 in (25.4 x 18.5 x 7.9 cm) 8-cell battery pack 11.0 x 7.3 x 3.1 in (27.9 x 18.5 x 7.9 cm) 16-cell battery pack

WEIGHT

5 lbs (2.3 kg) 8-cell battery pack 6 lbs (2.7 kg) 16-cell battery pack

SOUND LEVEL

39.9 dB(A) @ Setting 2

BATTERY DURATION

8-cell battery up to 8 hours

16-cell battery up to 16 hours

POWER

AC Power: 100-240 VAC, 50-60Hz DC Power: 11-18 VDC (10 max amp) Battery: Lithium Ion

WARRANTY

3 years (Includes sieve beds) 1 year on battery and power supplies

FreeStyle Comfort

Open the door to a bigger, more enjoyable life





Item Number MI396-1



BATTERY PACK

8-Cell Item Number

BT034-1S



16-Cell Item Number

BT035-1S



POWER SUPPLY WITH AC INPUT CABLE

Item Numbers

PW036-1S (Type A Plug)*

PW036-2S (Type C Plug)*



AC INPUT CABLE (Regional Options Available)

Item Numbers

4997-SEQ (Type A Plug)* 4998-SEQ (Type C Plug)*

DESKTOP BATTERY CHARGER (Battery Not Included)

Item Numbers BT036-1

(Type B Plug)*



VEHICLE POWER CORD (12 VDC Outlet)

Item Number CD041-1

POWER SUPPLY (120W 15V)

Item Number PW036-1



CARRY-ALL ACCESSORY BAG (Holds Unit and Accessories)

Item Number

MI372-2

Replacement Strap Item Number MI397-1



BACKPACK

Item Number MI459-1



TRAVEL ADAPTER

Item Number 5591-SEQ



GROSS PARTICLE FILTERS

Right Filter Item Number

FI226-1

Left Filter Item Number FI227-1





Find out more about FreeStyle Comfort and accessory products at

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

^{*} See https://www.iec.ch/worldplugs/map.htm# for details on plug types.

^{**} Purchase 5591-SEQ along with this item to properly adapt the plug for IEC type I power sockets.