

QUICK START GUIDE

myCAIRE™ PATIENTS

1 Download myCAIRE app from the Google Play or Apple app store, or myCAIRE.com.



2 Turn on Bluetooth from Settings.


Confirm Bluetooth is on.



3 Locate myCAIRE app onto your device.



4 Open the myCAIRE app and work with your home care provider to create an account and pair your device.



© Copyright 2021 CAIRE Inc. All Rights Reserved. SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information. CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable. myCAIRE is a telehealth solution intended to share specified device data from the enabled CAIRE oxygen concentrator. myCAIRE does not affect oxygen concentrator settings, performance, or oxygen therapy, and is not intended for patient monitoring, diagnosis or treatment. Please seek a doctor's advice before making any medical decisions.