

EDAT Software

Expanded data access to optimize fleet management

Eclipse Data Acquisition Tool (EDAT) is a computer-based data retrieval tool for use with the Eclipse 5®. EDAT can be utilized in investigating, identifying, and resolving device-related issues.



EDAT MAXIMIZES TECHNICAL RESOURCES IN MAINTAINING YOUR ECLIPSE FLEET

- Provides real-time diagnostic data to efficiently identify system issues and user errors, view active alarms and check battery functionality
- Offers simulated pulse mode for testing and calibration to ensure optimal operation
- Saves and prints reports for improved record keeping
- Allows installation of firmware upgrades and calibrates specific PCB replacement boards after installation



EDAT OFFERS A GLOBAL SOLUTION TO EXPAND SERVICE AND SUPPORT NEEDS

- Field service reps can transmit data and troubleshoot with in-house tech support personnel to quickly address equipment needs from the patients' home
- Optimize operations with improved fleet management by reducing service time and improving equipment utilization



Find out more about EDAT, Eclipse 5, and accessory products at

eclipseoxygen.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2021 CAIRE Inc. All Rights Reserved. ML-CONC0220 A