

HELIOS[™] U36 & U46

Liquid Oxygen Reservoirs

One Reservoir Fills all Top-Fill Portables

CAIRE Inc.'s versatile liquid oxygen reservoir works with our popular HELIOS Personal Oxygen System as well as Companion, Stroller, Hi-Flow, and Spirit portable liquid oxygen units. Available in two sizes, the flexible HELIOS Universal Reservoir can help you make optimal use of inventory to reduce your operating costs. The reservoir keeps delivery costs to a minimum because it is designed with a low-loss evaporative control system. The HELIOS Universal Reservoir can also provide oxygen directly to patients with an optional external flow control valve..

Features

- Improve asset utilization by using one reservoir for top-fill portable units
- Oxygen-saving features help reduce refill deliveries and operating costs
- Electronic contents indicator eliminates guesswork



Call 1.800.482.2473 to order!
www.CAIREmedical.com

HELIOS™ U36 & U46

Liquid Oxygen Reservoirs



Product Specifications

	U36	U46
LOX capacity	85 lb (38,6 kg)	110 lb (49,9 kg)
Gaseous equivalent capacity	29069 L	37916 L
Weight, empty	53 lb (24,0 kg)	60 lb (27,2 kg)
Weight, filled	138 lb (62,6 kg)	170 lb (77,1 kg)
Height	33.5 in. (851 mm)	37.5 in. (952 mm)
Diameter	15.4 in. (391 mm)	15.4 in. (391 mm)
Typical use time at 2 LPM	10.9 days	13.9 days
Operating pressure	22 psi (152 kPa)	22 psi (152 kPa)
Normal evaporation rate	1.2 lb/day (0,54 kg/day)	1.2 lb/day (0,54 kg/day)
Standard flow control range	Off, .5, .75, 1, 1.5, 2, 3, 4, 5, 6, 8, 10 LPM	Off, .5, .75, 1, 1.5, 2, 3, 4, 5, 6, 8, 10 LPM

Ordering Information

Description	Catalog Number
HELIOS U46 Universal Reservoir	B-702099-00
HELIOS U36 Universal Reservoir	B-702100-00
Accessories	
External flow control 22 psi, 1-10 LPM	B-701655-FLO
HELIOS roller base assembly	14880803
HELIOS moisture container	B-701546-00

CAIRE, Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107
Ph 770-721-7759 • Toll-Free 1-800-482-2473
Fax 888-932-2473

Chart Industries reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or, processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham RG41 2PL, United Kingdom
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245

Copyright © 2014 Chart Industries ML-LOX0006 F