PSA Oxygen Generator Reliant





The Reliant is an oxygen generator ideal for use in many applications. The Reliant's ability to deliver up to 95% oxygen concentration at flows of up to 17 SCFH (0.44 Nm³/hr or 8 LPM) make it a great fit for medical clinics requiring up to 50 psig (345 kPa or 3.45 barg).

Features


- Self-contained (internal air compressor)
- Microprocessor controlled
- Low operating cost
- Automatic and unattended operation
- Transportable design
- Easy to install and maintain
- Oxygen concentration monitor and alarm
- Optional auxiliary oxygen receiver (free standing)

Typical Applications

- Medical Oxygen Supply
- Veterinary Medicine

Note: The 220 V ~ \pm 10%, 50 Hz configuration is for medical use outside of the USA.

Specifications		
Product Flow ¹	17 SCFH (0.44 Nm ³ /hr or 8 LPM)	
Product Pressure	50 psig (345 kPa or 3.45 barg)	
Product Concentration	Up to 95%	
Product Dew Point	-100°F (-73°C)	
Dimensions (W x D x H) (Nominal)	22.5 x 24.5 x 36.25 in (57.2 x 62.3 x 92.1 cm)	
Weight	175 lb (79 kg)	
Physical Connections Product Gas Outlet	B size oxygen demand valve	
Sound Level	58 dBA @ 1 meter, open field conditions	
Ambient Operating Conditions	Locate the oxygen concentrator in a well-ventilated area that is protected from weather elements and remains between 40°F (4°C) and 112°F (44°C)	
Control Power Requirements (Single Phase)	220 V ±10%, 50 Hz, 3.0 A	
Typical Power Consumption	680 W	
60 Gallon Oxygen Receiver (Optional)		
Dimensions (Dia. x H)	26 x 52 in (65 x 131 cm)	
Weight	159 lb (72 kg)	



Ordering Information				
Model	Part Number	Description		
Reliant (FOR EXPORT ONLY)	AS072-3	220 V ~ \pm 10%, 50 Hz with oxygen concentration monitor (LED indicator)		
	AS072-5	220 V ~ ±10%, 50 Hz with oxygen concentration monitor (LED indicator), Russian Language		
Required Accessory	KI046-1	Accessory Kit (10 ft Stainless Braided Oxygen Hose & Primary/Secondary Oxygen Ball Valve Assembly)		
Optional Accessories	KI048-1	Accessory Kit (10 ft Stainless Braided Oxygen Hose & Oxygen Receiver Purge Ball Valve)		
	KI168-1	Accessory Starter Kit		
	TA150-1	60 Gallon Oxygen Receiver		
	AN021-1	Oxygen Analyzer (Maxtec Handi)		
	AN005-1	Oxygen Analyzer (Maxtec Max O ₂ Plus)		
	AN075-1	Oxygen Analyzer/Sensor (Maxtec Max O ₂ Plus)		
Shipping Information		Reliant	60 Gallon Oxygen Receiver	
Class		92.5	70	
Commodity Classification Number		8421.39.8040	7311.00.0000	
Dimensions (W x D x H)		28 x 30 x 45 in (71 x 76 x 114 cm) with Accessory Kit	28 x 30 x 56 in (71 x 76 x 142 cm)	
Gross Weight		233 lb (106 kg) with Accessory Kit	177 lb (80 kg)	

Warranty: 1 Year Parts and Factory Labor*

* An unprotected or inadequately ventilated environment, or improper control power may cause damage to the oxygen generator

not covered under warranty.

All performance ratings based on an ambient temperature up to 100°F (38°C), up to 1,000 feet elevation, and 80% relative humidity.





260 Creekside Drive Buffalo, NY 14228-2075 U.S.A. Tel: (716) 691-0202
Fax: (716) 691-1255 www.airsep.com
info@airsep.com ML-IND0057 A



© 2020 AirSep Corporation. All Rights Reserved. AirSep reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. Any rights expressly stated herein are reserved by us, as applicable.