


FREESTYLE[®] 5

Portable Oxygen Concentrator

Unleash your potential

The high capacity and wearable FreeStyle[®] 5 offers you the ultimate freedom to go anywhere. Easy to use while active in your everyday activities and FAA-approved for use on a commercial flight — it is perfect for everywhere you go.




A LIGHTWEIGHT, HIGH CAPACITY SOLUTION

Replace your tanks with the lightweight, wearable FreeStyle[®] 5. This device can be easily carried by handle, worn by a comfortable shoulder strap, or attached to a backpack. You really can go anywhere and enhance your activity level.



PEACE OF MIND

Enjoy the convenience of a portable oxygen concentrator with sensitive breath detection — UltraSense[™] Technology — which ensures the device detects each of your breaths quickly to give you the oxygen you need each time you breathe in.



STAY ACTIVE

Quickly swap the battery packs on the go so you never have to worry about running out of oxygen. The convenient pack will re-charge anytime wall (AC) or motor vehicle (DC) power is connected and is easily accessible from the front of the device.



CAIRE Contact Information:

Providers Toll Free
1.800.482.2473
Patients Toll Free
1.877.704.0878

www.cairemedical.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

FREESTYLE[®] 5

Portable Oxygen Concentrator

Includes

- Portable Concentrator
- Carry Bag
- Battery Pack
- AC Power Supply
- DC Power Supply
- Carry-All Accessory Bag
- Shoulder Strap

Product Features

- User-replaceable battery pack with up to 2½ hour battery life at setting 2
- Easy to carry with a handle, or wear with a shoulder strap, or backpack
- FAA-approved for commercial airline use
- Fully functional on wall (AC) or motor vehicle (DC) power.
- Simple, one-touch operation

Specifications

FLOW RATE

Pulse 1-5

DIMENSIONS

10.7 x 6.6 x 4.4 in
(27.2 x 16.8 x 11.2 cm)

WEIGHT

6.7 lbs (3.0 kg)

BATTERY DURATION

Setting 1 - 4.25 hours (4:15)
Setting 2 - 2.5 hours (2:30)
Setting 3 - 2.0 hours (2:00)
Setting 4 - 1.75 hours (1:45)
Setting 5 - 1.25 hours (1:15)

BATTERY RECHARGE

4 – 5 hours

WARRANTY

3 years FreeStyle[®] 5
1 year battery and
power supplies



Available Accessories



CARRY BAG



AC POWER SUPPLY
(CORDS AVAILABLE
SEPARATELY)



DC POWER SUPPLY



BACKPACK HARNESS



EXTRA BATTERY PACK



DESKTOP CHARGER



SPARE INTAKE FILTER



CARRY-ALL
ACCESSORY BAG

CAIRE Contact Information:

Providers Toll Free
1.800.482.2473
Patients Toll Free
1.877.704.0878

www.cairemedical.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.
Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

ML-CONC0123 D