

Hi, just a reminder that you're receiving this email because you have expressed an interest in CAIRE Inc. Don't forget to add taylor.erwin@chartindustries.com to your address book so we'll be sure to land in your inbox.

You may [unsubscribe](#) if you no longer wish to receive our emails.



CAIRE®

A Chart Industries Company

AirSep®  **SeQual®**

www.CAIREmedical.com

Med Tips

In This Issue

[Helpful Hints / FAQs](#)
[Service Schools &
Tradeshows](#)
[Contact Us](#)

PDF Version

To view Med Tips in PDF,
[click here.](#)

Join Our List

[Join Our Mailing List!](#)

April 2014

PRODUCT INFORMATION

CAIRE Launches SeQual eQuinox™ Portable Oxygen Concentrator

CAIRE is pleased to announce the release of the SeQual eQuinox portable oxygen concentrator (POC).

Weighing in at only 14 pounds, the newly designed, easy-to-maneuver frame makes the SeQual eQuinox the lightest weight POC able to offer 3.0 LPM continuous flow. This small package is powerful enough to give patients the freedom to enjoy life on the go. The eQuinox is fully functional on AC Power, DC Power, and battery power. Wherever patients go, the eQuinox goes with them and keeps up with the various activities of daily living. With the option of a longer lasting 24-cell battery pack, the eQuinox can achieve 2.75 hours of operation at 2.0 LPM.

The Multi-Language Voice Interface offers a new layer of comfort to users by providing verbal confirmation of flow rate settings, battery times and alarms. The user-friendly interface is adjustable to support a variety of patient

languages.

Please note: The eQuinox has only been released in the US at this time. The eQuinox is expected to be released outside of the US in the third quarter of this year.

Features

- **NEW: Multi-Language Voice Interface** makes using a POC simpler than ever
- **NEW: Only 14 lbs.** with 12-cell battery pack installed
- **NEW: Up to 2.75 hours** of battery life at 2.0 LPM with the optional 24-cell battery pack
- **NEW: Only 45 dB** at 2 continuous flow
- **NEW: Service Cabinet** makes filters and the 9V battery accessible for service technicians without removing the case - performing maintenance and upkeep is easier than ever
- **NEW: DC Power supply** consists of a lightweight 12V cable - no transformer brick
- **NEW: Screen** intermittently displays battery charge remaining in minutes
- Approved for Rx of 0.5 - 3.0 LPM Continuous Flow
- 9 Pulse Settings with bolus sizes from 16mL to 192mL
- autoSAT® Technology helps maintain a consistent FiO2 by servo-controlling the device to meet changing respiratory rates
- Three available rise time options: Fast, Medium, Slow
- Three adjustable Pulse Dose Sensitivity settings with a maximum sensitivity of -0.13 cm H2O
- FAA Approved for use on commercial aircraft
- Full functionality on DC power
- Ability to recharge battery cartridge on DC power when excess power is available

Providers can contact their local sales representative to determine availability in their area.





AirSep SureFlow Oxygen Flow Station

AirSep SureFlow Oxygen Flow Station is a device that allows oxygen flow to be split from a single source and distributed to multiple patients. The SureFlow unit enables up to 50 feet of pediatric or adult cannula and safety tubing to be used at the end of its flow meter oxygen outlet ports. Two or more SureFlow oxygen flow stations can be connected together to serve even more patients, thereby greatly reducing the cost of providing therapeutic oxygen treatments. SureFlow works with most oxygen concentrators and oxygen sources from 7 to 50 psig (48-345 kPa). SureFlow can be wall-mounted or can sit on a desktop.

When ordering the SureFlow, use part number **FM069-1** for the five 0 - 1 LPM (1/10 LPM increments) flow meter option or part number **FM069-2** for the five 0 - 2 LPM (1/8 LPM increments) flow meter option. Either SureFlow package will include 2 DISS Adaptors (1/4 in. Barb), 1 - 7ft. length of cannula tubing (1/4 in. I.D.), 1 - 7ft. length of pneumatic hose with DISS fittings on each end, and a User's Manual.





HELPFUL HINTS / FAQs

Liquid Oxygen Time Use Chart

Q: How long will my liquid oxygen portable and reservoir units last at various flow rates?

A: CAIRE has a published liquid oxygen Time Use Chart for all of its liquid oxygen products which can be downloaded, printed, or viewed at the link below:

<http://www.cairemedical.com/getattachment/9c10376d-210e-4059-8658-ecb1fbe49b7f/.aspx>

The chart also includes specifications such as liquid/gaseous oxygen capacity, typical normal evaporation rates (NER), heights, and weights when empty or full.

ACCESSORIES

Universal Power Adapter Plug Set

The Universal Power Adapter includes adapters to connect into standard wall outlets in Europe, UK, New Zealand, North/ South America, China, Japan, and Taiwan and is compatible with all portable and stationary oxygen concentrators. The Universal Power Adapter Plug Set can be ordered using **PN 5591-SEQ**.





Universal Liquid Oxygen Portable Pull-Behind Cart

CAIRE now offers a universal LOX portable cart for all of our portables which can be ordered using **PN 20773616**. This will replace previous LOX portable cart part numbers **B-775295-00, 97404539, and 12917730**. The Universal LOX Portable Cart has new dimensions and a flexible velcro strap that makes it compatible with all of CAIRE's LOX portables: HELiOS H300 and H850, Companion 1000 and 1000T, Sprint, Stroller, Spirit family, and Hi Flow Stroller.



SERVICE SCHOOLS & TRADESHOWS

Mark Your Calendars for our Upcoming 2014 Service Schools!

USA

2014 Service School dates have been released!

For additional details and registration information about these trainings, please visit the ["Events & Service Schools"](#) tab of our website.

Concentrator Training: June 9, 10, 11 in Ball Ground, GA

LOX Training: June 12 & 13 in Ball Ground, GA

**Ball Ground trainings include a tour of our CAIRE and SeQual manufacturing facility.*

Concentrator Training: August 18, 19, 20 in Buffalo, NY

LOX Training: August 21 & 22 in Buffalo, NY

**Buffalo trainings include a tour of our AirSep manufacturing facility.*

Concentrator Training: November 10, 11, 12 in Ball Ground, GA

LOX Training: November 13 & 14 in Ball Ground, GA

**Ball Ground trainings include a tour of our CAIRE and SeQual manufacturing facility.*

Europe

2014 Service School dates have been released!

For additional details and registration information about these trainings, please visit the ["Events & Service Schools"](#) tab of our website.

English LOX Training: April 15 & 16 in Wuppertal, Germany

English Concentrator Training: May 20, 21, 22 in Wuppertal Germany

French LOX Training: June 2 & 3 in France

French Concentrator Training: June 4, 5, 6 in France

Italian Concentrator Training: June 10, 11, 12 in Padova, Italy

Spanish LOX Training: June 16 & 17 in Barcelona, Spain

Spanish Concentrator Training: June 18, 19, 20 in Barcelona, Spain

German LOX Training: September 9 & 10 in Wuppertal, Germany

German Concentrator Training: October 14, 15, 16 in Wuppertal, Germany

English LOX Training: October 28 & 29 in Wokingham, UK

Italian LOX Training: October 29 & 30 in Padova, Italy

English Concentrator Training: November 18, 19, 20 in Wokingham, UK

Tradeshows

We will be attending the following trade shows this month. We would love for you to stop by our booth!

NYMEP

March 31-April 1
White Plains, NY
CAIRE

FAHCS

April 2nd
Orlando, FL
CAIRE

CO-SRC Annual Meeting

April 30-May 2
Vail, CO
CAIRE

CONTACT US

For ordering information, contact Customer Service:

United Kingdom	+44 (0) 1189 367060
France	+33 (0) 561 429 411
Germany	+49 (0) 202 739 55420
Italy	+39 049 879 9601
	customerservice.europe@chartindustries.com
Australia/New Zealand	+61 2 9749 4333
Asia/Pac Rim	customerservice.australia@chartindustries.com
	770.721.7759
	csasia@chartindustries.com
The Americas	770.721.7759
US Toll-Free	800.482.2473
	customerservice.usa@chartindustries.com

For technical information, contact Technical Service:

United Kingdom	+44 (0) 1189 367060
France	+33 (0) 561 429 411
Germany	+49 (0) 202 739 55420
Italy	+39 049 879 9160
	techservice.europe@chartindustries.com
Asia/Pac Rim	+61 2 9749 4333
	techservice.asia@chartindustries.com

The Americas	770.721.7759 techservice.usa@chartindustries.com
US Toll-Free	800.482.2473 techservice.usa@chartindustries.com

[Forward email](#)



Try it FREE today.

This email was sent to snystrom@sequal.com by taylor.erwin@chartindustries.com | [Update Profile/Email Address](#) | Instant removal with [SafeUnsubscribe™](#) | [Privacy Policy](#).

Chart Industries | 2200 Airport Industrial Dr. | Ste. 500 | Ball Ground | GA | 30107