

Product Updates

Liberator Tube Color Change

Going forward customers may receive part number 14224654 as either the current black color or a white color. This is the length of tube over the breathing coil near the relief/economizer assembly. The two colors are interchangeable in terms of service needs.



Discontinued ATF Spare Compressors

Due to change of provider and product, the compressors listed below are obsolete and are no longer available for pairing with ATF modules. Contact CAIRE Technical Service for assistance in selecting a compressor for your ATF oxygen generator systems.

Part Number	Description	Replaced By
7495-SEQ	Packaged Compressor, 115 V/60 Hz	Contact CAIRE Technical Service
1313-SEQ	Packaged Compressor, 120 V/60 Hz	Contact CAIRE Technical Service
1243-SEQ	Packaged Compressor, 115 V/60 Hz	Contact CAIRE Technical Service
1244-SEQ	Packaged Compressor, 220-230 V/50-60 Hz	Contact CAIRE Technical Service
1318-SEQ	Packaged Compressor, 115/230 V	Contact CAIRE Technical Service

Helpful Hints & FAQs

Eclipse 5[®] Preventative Maintenance Kit Part Numbers

Part number 5022-SEQ includes all the replacement parts needed to perform preventative maintenance on Eclipse[®] units. Previously, the PM kit included two different product filters, which allowed it to be compatible with every model Eclipse[®] units produced with the HEPA filter internal to the product tank and units with the disk type HEPA filter external to the product tank. We will no longer supply HEPA Filter 6986-SEQ since the previous style product tank was discontinued 2011. If you require the External type HEPA filter you will need to purchase the new product tank assembly, PN SP4378-SEQ.



FAA Guidelines

Q.) Why can't I find my portable oxygen concentrator listed on the FAA site?

A.) On May 24, 2016 the Federal Aviation Administration (FAA) issued a final rule to replace the existing process by which the Agency approves portable oxygen concentrators (POC) for use on board aircraft in air carrier operations. This final rule replaced the existing FAA approval process in Special Federal Aviation Regulation (SFAR) No. 106 to currently allow passengers to use a POC on board an aircraft if the POC satisfies the acceptance criteria and labeling requirement stated in the final rule. The labeling requirement only affects POCs intended for use on board aircraft that were not previously approved prior

to May 24th ruling for use on aircraft by the FAA2. The FreeStyle Comfort is the only CAIRE product that has been released after the May 24, 2016 ruling and is therefore not included on the SFAR 106 listing. The FreeStyle Comfort is labelled according to FAA labelling requirements.

The link to the Final Rule for the Acceptance Criteria for Portable Oxygen Concentrators Used On-Board Aircraft can be found here:

<https://www.federalregister.gov/documents/2016/07/05/2016-15770/acceptance-criteria-for-portable-oxygen-concentrators-used-on-board-aircraft-correction>

Link to a CAIRE letter referencing self-certification of portable oxygen concentrators:

http://files.caireinc.com/CAIRE_FAA_LETTER_Updated_4_26.pdf

Q.) What CAIRE portable oxygen concentrators meet FAA guidelines?

A.) Here are the CAIRE portable oxygen concentrators that meet FAA guidelines:

AirSep® Focus™

FreeStyle® Comfort®

AirSep® FreeStyle®

AirSep® FreeStyle® 5

AirSep® LifeStyle™

SeQual® Eclipse® (All Generations)

SeQual® eQuinox™ Oxygen System
(model 4000)

SeQual® Oxywell™ Oxygen System
(model 4000)

SeQual® SAROS™

Demonstration Videos

Q.) Where can I find any videos available that shows end users about how to use our concentrators?

A.) Visit our CAIRE YouTube page at youtube.com/c/caireinc

New Website

Q.) What is the new web address for CAIRE Medical?

A.) The new CAIRE website is www.caireinc.com

The Respiratory Healthcare Knowledge Center is now located at:

<https://www.caireinc.com/providers/business-support-programs/knowledge-center/>

The manufacturer of this POC has determined this device conforms to all applicable FAA acceptance criteria for POC carriage and use on board aircraft.



LA438-1 B

**CONTACT:**

For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

techservice.global@caireinc.com

Europe: +44 (0) 118 936 7060

techservice.global@caireinc.com

For ordering information or general inquires, contact Customer Service:

USA: 800-482-2473

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com

© Copyright 2019 CAIRE Inc. All Rights Reserved. CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.