



Pulse Bolus Size Chart

for each *BREATH*
your patient takes

2 is not 2 is not 2




CAIRE
A Chart Industries Company

Call 1.800.482.2473 to order!
www.chartindustries.com/Respiratory-Healthcare

MILD

2 is not 2 is not 2

CAIRE's solutions give providers the ability to help both the most clinically difficult patient and the most difficult to please—all with one portable product portfolio.


MODERATE

2 is not 2 is not 2

2 is not 2 is not 2

SEVERE

CAIRE solutions for LIFE.



Titrate to Saturate



CAIRE solutions for LIFE.

Breath Rate of 10 bpm

		Pulse Volume	Pulse Duration	Pulse Delay	Weight
	Setting	ml	ms	ms	lbs
Focus	1	33.6	393	74	1.75
Freestyle 3	1	23.8	323	70	4.4
	2	41.5	454	69	
	3	50.1	509	71	
Freestyle 5	1	25.9	406	62	6.2
	2	46.4	568	62	
	3	64.8	653	63	
	4	86.7	878	66	
	5	103.3	924	62	
Eclipse 5	1	17.8	126	35	18.4
	2	34.1	225	34	
	3	50.9	310	34	
	4	64.7	377	32	
	5	81.3	432	32	
	6	102	507	32	
Equinox	1	18.0	117	34	14
	2	33.4	217	34	
	3	49.5	301	33	
	4	62.6	368	33	
	5	78.1	428	33	
	6	97.9	489	31	

Breath Rate of 20 bpm

		Pulse Volume	Pulse Duration	Pulse Delay	Weight
	Setting	ml	ms	ms	lbs
Focus	1	17.2	194	73	1.75
Freestyle 3	1	12.6	166	163	4.4
	2	21.4	227	71	
	3	26.0	250	70	
Freestyle 5	1	13.1	198	63	6.2
	2	23.3	270	63	
	3	33.7	327	62	
	4	45.5	448	65	
	5	48.4	424	62	
Eclipse 5	1	17.4	122	35	18.4
	2	34.2	223	34	
	3	51.0	310	34	
	4	65.6	371	32	
	5	81.3	425	32	
	6	102.1	477	34	
Equinox	1	17.1	122	34	14
	2	33.2	218	34	
	3	49.5	304	33	
	4	62.7	366	32	
	5	78.4	420	32	
	6	95.7	473	33	

Breath Rate of 30 bpm

		Pulse Volume	Pulse Duration	Pulse Delay	Weight
	Setting	ml	ms	ms	lbs
Focus	1	11.7	134	72	1.75
Freestyle 3	1	8.6	111	72	4.4
	2	14.6	152	70	
	3	17.6	170	70	
Freestyle 5	1	9.0	133	63	6.2
	2	14.9	167	62	
	3	23.1	224	62	
	4	30.1	298	65	
	5	32.4	282	62	
Eclipse 5	1	17.4	121	34	18.4
	2	34.1	222	34	
	3	50.7	306	34	
	4	65.2	368	33	
	5	81.6	425	32	
	6	101.2	474	34	
Equinox	1	17.1	123	33	14
	2	33.2	218	32	
	3	49.2	301	34	
	4	62.9	365	33	
	5	78.3	419	33	
	6	98.1	469	32	

Chart Industries reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or, processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Copyright © 2015 Chart Industries
ML-CONC0075 A



Call 1.800.482.2473 to order!