

SeQual® eQuinox™

PORTABLE OXYGEN CONCENTRATOR

THE PREMIER PORTABLE CONCENTRATOR

This quiet and advanced portable oxygen concentrator is requested by patients because of its small size and countless features. The eQuinox™ will help drive your referrals and can treat many of your oxygen users at various activity levels.

- **The Most O₂ to Treat the Most Users**

Producing up to 3.0 LPM continuous flow and with pulse settings 1-9 with a bolus size up to 192 mL the highest of any concentrator, the eQuinox™ can cover a wide variety of oxygen users at high levels of ambulation.

- **The Most User-Friendly Features**

Users and providers alike will get added peace of mind with the unique multi-language voice module that verbally announces any alarms, low batteries, or changes in settings. Other features include an easy access front-loading battery, a 24-cell battery for extended durations, and a hour:minute display of remaining battery life.

- **A True 24/7 Device**

Robust enough to run 24/7, it can be used during sleep, at home, on the go, and on the road — an all-in-one oxygen solution.

Continuous flow options eliminate the worry of desaturation during sleep, and the high pulse output allows patients to be active with a portable device and lengthy battery life.

- **Advanced Clinical Features**

Featuring autoSAT® Technology that detects the user's breath rate and automatically adjusts the motor speed to ensure oxygen is always delivered with every breath. This helps to maintain consistent FiO₂, even in high flow patients.

The eQuinox allows oxygen delivery to be tailored to each individual user for the best therapy. Both sensitivity and bolus rise time can be adjusted to ensure proper treatment.




CAIRE®
A Chart Industries Company
www.CAIREMedical.com



** SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.
PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.
WARNING: AS WITH ANY ELECTRICALLY POWERED DEVICE, THE USER MAY EXPERIENCE PERIODS OF NON-OPERATION AS A RESULT OF ELECTRICAL POWER INTERRUPTION, OR THE NEED TO HAVE THE EQUINOX SERVICED BY A QUALIFIED TECHNICIAN. THE EQUINOX IS NOT APPROPRIATE FOR ANY PATIENT WHO WOULD EXPERIENCE ADVERSE HEALTH CONSEQUENCES AS THE RESULT OF SUCH TEMPORARY INTERRUPTION.

© 2017 CAIRE Inc. All Rights Reserved ML-CONC0088 Rev D

PROVIDER FEATURES

- 24/7 capability and reliability
- Service door for easy maintenance without opening cabinet
- Bill both E1390 & E1392 codes
- Drive referrals and retail sales
- Ideal for travel programs and hospital discharge programs
- Diagnostic alarm codes for easier service calls

PATIENT FEATURES

- Run on AC, DC or battery power
- Only 14 lbs (6.4 kg)
- Multi-language voice interface
- Up to 5.9 hour battery life at 2 Pulse
- Front-loading battery
- HH:MM Battery Life Display
- FAA approved
- Robust, adjustable cart

CLINICAL FEATURES

- Proprietary autoSAT® Technology
- Adjustable sensitivity and rise time
- Bolus sizes up to 192 mL
- O₂ monitor standard on all units
- Highest pulse dose of any concentrator
- Ideal for pulmonary rehab patients

SPECIFICATIONS

Continuous Flow	0.5–3.0 LPM
Pulse Flow	Settings 1-6 (16-96 mL) Settings 7-9 (128 – 192 mL)
Weight w/12-Cell Battery	14.0 lbs (6.4 kg)
Dimensions	13.6 x 10.6 x 7.4 inches (34.5 x 26.9 x 18.8 cm)
O₂ Concentration	90% (+5.5 / -3 %)
Sound Level	45 dB(A) @ 2.0 LPM
Operating Temp.	50°F to 104°F (10°C to 40°C)
Operating Altitude	-1254 to 13123 ft (-382 to 4000 m)
12-Cell Battery Weight	1.85 lbs (0.84 kg)
12-Cell Battery Duration	2.8 hours @ Pulse 2 1.2 hours @ 2LPM
24-Cell Battery Weight	3.6 lbs (1.6 kg)
24-Cell Battery Duration	5.9 hours @ Pulse 2 2.8 hours @ 2LPM
Battery Re-Charge	1.5 – 3.5 hours
AC Power	87-264 VAC, 47-63Hz
DC Power	12 V Nominal
Power Draw	110 W @ 3LPM 50 W Battery Charging
Warranty	3 Years eQuinox** 1 Year Battery & Accessories
Maintenance Schedule	Air Inlet, Compressor Intake, HEPA Filters – 1 Yr Replacement

ORDERING INFORMATION

eQuinox™ Includes: Portable Concentrator, Battery, AC Power Supply & Cord, DC Power Cord, Removable Cart, Accessory Bag, Spare Intake Filter

4807-SEQ	USA, 12 cell batt
4837-SEQ	USA, 24 cell batt
4807-08-SEQ	UK, 12 cell batt
4837-08-SEQ	UK, 24 cell batt
4807-C1-SEQ	EU, 12 cell (EN, FR, DU, ES, IT, PT, NL)
4837-C1-SEQ	EU, 24 cell (EN, FR, DU, ES, IT, PT, NL)
4807-C2-SEQ	EU, 12 cell (SV, DA, NO, FI)
4837-C2-SEQ	EU, 14 cell (SV, DA, NO, FI)

other countries' power supplies may be available - contact customer service

ACCESSORIES

20952900	12-Cell Battery
20952906	24-Cell Battery
4920-SEQ	Accessory Bag
4830-SEQ	DC Power Cable
4836-SEQ	AC Supply (no cord)
4997-SEQ	AC Cord (US)
4998-SEQ	AC Cord (EU)
4999-SEQ	AC Cord (UK)
4980-SEQ	Desktop Charger (US) (cord included)




CAIRE®
 A Chart Industries Company
www.CAIREMedical.com



**** SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.**
 PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.
WARNING: AS WITH ANY ELECTRICALLY POWERED DEVICE, THE USER MAY EXPERIENCE PERIODS OF NON-OPERATION AS A RESULT OF ELECTRICAL POWER INTERRUPTION, OR THE NEED TO HAVE THE EQUINOX SERVICED BY A QUALIFIED TECHNICIAN. THE EQUINOX IS NOT APPROPRIATE FOR ANY PATIENT WHO WOULD EXPERIENCE ADVERSE HEALTH CONSEQUENCES AS THE RESULT OF SUCH TEMPORARY INTERRUPTION.

© 2017 CAIRE Inc. All Rights Reserved ML-CONC0088 Rev D