



APRAGAZ



Chée de Vilvorde, 156
 Vilvoordestwg, 156
 B-1120 Bruxelles – Brussel
 Tel. : 32 (0)2 264 03 60
 Fax. : 32 (0)2 268 89 58
 http://www.apragaz.com
 e-mail: info@apragaz.com

NOTIFIED BODY EC TYPE EXAMINATION CERTIFICATE		P16526/009
In accordance with module B of the Pressure Equipment Directive 97/23/EC		Certificate N° 16/US/3920-0 Rev.1
Manufacturer	Name	AIRSEP CORPORATION
	Address	260 Creekside Drive NY 14228 BUFFALO USA
Description :	AS-Series – Oxygen Generator	Equipment : / Assembly of : Vessel, accessories, safety valves Module : B
Serial numbers :	AS-D+, AS-E, AS-G, AS-J, AS-K, AS-L	
Drawing numbers :	AS-097 (sheet 4)-Rev.D, AS-D* (sheet 7)-Rev.E, AS-103 (sheet 4)-Rev.H, AS-E (sheet 7)-Rev.E, AS-104 (sheet 4)-Rev.E, AS-G (sheet 7)-Rev.E, AS-105 (sheets 3 & 5)-Rev.J, AS-J (sheet 7)-Rev.J, AS-106 (sheet 5)-Rev.F, AS-K (sheet 7)-Rev.E, AS-107 (sheets 8 & 9)-Rev.F, AS-L (sheet 10)-Rev.F	
Maximum allowable Pressure (PS) :	0/100 PSI	
Minimum / Maximum allowable temperature (TS) :	40/114°F	
Safety valves :		
Capacity :	ST25 : 94,6 SCFM (AS-D+, AS-E, AS-G) SF50 : 190,3 SCFM (AS-J, AS-K) SB50 : 412,8 SCFM (AS-L)	
Set pressure :	100 PSI	
Hydrostatic test		
Test Pressure (PT) :	9,9 bar (sub-assemblies)	
The undersigned certifies that the design of this pressure equipment conforms with the requirements of the Directive 97/23/EC on the basis of the information supplied. The manufacturer is allowed to affirm the mark <<CE0029>> to approved equipment in the conditions described in article 15 of the Directive, only if the requirements stated in module D are fully complied with.		
Date : May 4, 2016	Name : F. Romans	Position : Technical Department
Notified body identification number :	0029	Signature
Notified body stamp :	APRAGAZ Belgium Inspecting Authority	
The approval is valid until : May 4, 2026		
Notified body reference :	P16526/009	Eng. F. Yernaux Head Technical Department