

Comparison of portable oxygen concentrators using a COPD patient simulation model

Rachel Culbreth, PhD, MPH, RRT; Robert Murray, MS, RRT; Kyle Brandenberger, PhD; Douglas S. Gardenhire, EdD, RRT, RRT-NPS, FAARC
Department of Respiratory Therapy, Byrdine F. Lewis College of Nursing and Health Professions, Georgia State University, Atlanta, GA, USA

Introduction

- Approximately 1.5 million Americans use home oxygen therapy.¹
- Traditional in-home oxygen therapy requires heavy oxygen tanks for patients to transport, limiting mobility and quality of life.¹⁻²
- With the advent of new portable oxygen concentrator (POC) technology, patients are able to travel with their oxygen therapy and have reported a higher quality of life.¹⁻²
- However, studies have shown mixed efficacy in delivering adequate oxygen for various patient scenarios.¹⁻⁵

Objectives

- Determine differences in delivered FiO₂ among POCs at varying respiratory rates (15 breaths per minute (bpm), 20 bpm, 30 bpm, and 40 bpm).
- Examine differences in delivered FiO₂ among POCs and control group oxygen delivery devices using a COPD patient lung simulator.





Methods

- This study utilized IngMar Medical Active Servo Lung 5000 to simulate a COPD patient's lungs.
- POC devices were compared to control devices, including wall oxygen and standalone oxygen concentrators at settings 2, 3, 5, and 6.
- Descriptive statistics and ANOVAs were computed to determine statistically significant differences between POCs and control devices.





Across most respiratory rate scenarios, the wall oxygen (control) obtained a higher FiO₂ compared to the POCs.

Excluding control groups: CAIRE FreeStyle Comfort achieved a higher FiO₂ for the majority of the scenarios.

Inogen G4 obtained a higher FiO₂ on setting 2 for 20 bpm.



Results



Conclusions

- FiO₂ measurements varied across POCs depending on the breathing rate scenario
- The CAIRE FreeStyle Comfort achieved a higher FiO₂ compared to all other POCs in 7 out of the 8 scenarios.
- Clinical providers should account for their patients' respiratory rate demands when recommending specific POCs

Acknowledgement

This study was funded with unrestricted grant from CAIRE and conducted in an independent research lab at Georgia State University by the authors.

References

1. Jacobs SS, Lederer DJ, Garvey CM, Hernandez C, Lindell KO, McLaughlin S, et al. Optimizing Home Oxygen Therapy. An Official American Thoracic Society Workshop Report. *Ann Am Thorac Soc* 2018;**15**:1369–81. <https://doi.org/10.1513/AnnalsATS.201809-627WS>.
2. Martin DC. Contemporary portable oxygen concentrators and diverse breathing behaviours -- a bench comparison. *BMC Pulm Med* 2019;**19**:. <https://doi.org/10.1186/s12890-019-0980-x>.
3. Khor YH, McDonald CF, Hazard A, Symons K, Westall G, Glaspole I, et al. Portable oxygen concentrators versus oxygen cylinder during walking in interstitial lung disease: A randomized crossover trial. *Respirol Carlton Vic* 2017;**22**:1598–603. <https://doi.org/10.1111/resp.13083>.
4. Yáñez AM, Prat JP, Álvarez-Sala JL, Calle M, Díaz Lobato S, García Gonzalez JL, et al. Oxygenation With a Single Portable Pulse-Dose Oxygen-Conserving Device and Combined Stationary and Portable Oxygen Delivery Devices in Subjects With COPD. *Respir Care* 2015;**60**:382–7. <https://doi.org/10.4187/respcare.03470>.
5. Chatburn RL, Williams TJ. Performance comparison of 4 portable oxygen concentrators. *Respir Care* 2010;**55**:433–42.