

Companion 5™

Stationary Oxygen Concentrator

CAIRE®
O₂ is what we do.



TIME-PROVEN DESIGN

PATIENT SATISFACTION

EASE OF SERVICE



The foundation of your patients' oxygen therapy

Designed to Last

The Companion 5 is a state-of-the-art stationary concentrator offering an updated sleek and compact design capable of addressing most of your patients' oxygen delivery needs.

Time-proven design

- **Industry-leading reliability**—greater than 99% reliability rate out of the box
- **Comprehensive warranty**—3-year warranty on the device
- **Proprietary autoFLOW™ feature**—variable cycle time based on flow rates to ensure operation at peak efficiency
- **Proven durability**—engineered with oversized compressor and sieve beds to minimize disruptions in service

Patient satisfaction built-in

- **Effective oxygen delivery**—up to 95.5% concentration from 0.5 to 5.0 LPM, meeting a wide variety of patient prescriptions
- **Simple user interface**—easy to understand alerts
- **Energy efficient**—reduced power consumption to minimize the burden on the patient's energy bill
- **Sleek and compact design**—blends into surroundings

Ease of service

- **Effortless maintenance**—proven quality performance to minimize unscheduled service
- **Modular design**—fewer components reduce time and effort to service the device
- **Enhanced diagnostics**—on board alarm alerts enhance the ability to quickly troubleshoot
- **Telehealth connectivity**—the only stationary concentrator with wireless capability for connection to myCAIRE™ Telehealth Solution



Companion 5 Specifications

FLOW RATE

0.5 – 5.0 LPM

DIMENSIONS

21.5 x 12.5 x 13.5 inches
(54.6 x 31.8 x 34.2 cm)

WEIGHT

36 lbs (16.3 kg)

SOUND LEVEL

50 dB(A) at 2 LPM
MDD Version:
43.74 ± 0.45 dB(A) at 3 LPM*

POWER

285 W @2 LPM,
350 W Maximum

MAXIMUM TUBING

50 ft (15.2 m)

WARRANTY

3 years



CE-certified

Accessories

NASAL CANNULA

Adult

7 ft. (2.1 m) CU002-1
25 ft. (7.6 m) CU002-4

Pediatric

7 ft. (2.1 m) CU002-7



HUMIDIFIER ADAPTOR

(9 in. / 23 cm.)

Item Number 20843882



Intake Filter

Item Number 15062925



HUMIDIFIER BOTTLE

Item Number HU003-1



MASK

Adult, Concentration
(Over the Ear Style)

7 ft. (2.1 m)
CU006-1



Pediatric, Concentration

7 ft. (2.1 m)
CU007-1



OXYGEN TUBING

7 ft. (2.1 m) CU004-2
25 ft. (7.6 m) CU004-3
50 ft. (15.2 m) CU004-1

Serving your patients across the continuum of CAIRE.

Non-Delivery Solution

Reduce operating expenses and time in the field while improving patient care with the CAIRE non-delivery solution—the combination of the industry-leading FreeStyle® Comfort® portable oxygen concentrator with the time-proven Companion 5 stationary oxygen concentrator.



O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy.

CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.



Find out more about Companion 5 and accessory products at

CAIREinc.com


CAIRE®

*per test standard Nr. 14-110/2018 MDS-Hi
SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.
Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2022 CAIRE Inc. All Rights Reserved. ML-CONC0142 G