



Companion 5™ VisionAire™ 5 Stationary Oxygen Concentrators



The durable VisionAire stationary oxygen concentrator is known for its reliability, efficiency, and low maintenance.

VISIONAIRE 5

The robust construction, ease-of-use, and low maintenance make it the perfect choice for facilities, hospitals, clinical settings, and harsh environmental conditions.

- Continuous Flow Settings 0.5 to 5 LPM
- Weight: 30 lbs. (13.6 kg)
- Output pressure: 9 psi (62 KPa)
- Dimensions: 20.8 inches H X 14.1 inches W X 11.5 inches D
- Power draw: 290 watts
- No preventative maintenance required



**VisionAire 5
(AS098-4)**

The
**VisionAire
and Companion 5**
stationary oxygen
concentrators are available
through CAIRE.

Visit our website at
www.caireinc.com.

**COMPACT. POWERFUL.
ENERGY EFFICIENT.**



**Companion 5
(15067005)**

COMPANION 5

One of the most technologically-advanced concentrators available on the market. The simple user interface, and reduced power draw are perfect for patients and care facilities, while the durability, reliability, and ease of maintenance are the perfect combination for equipment providers. Also, patent-pending autoFLOW™ technology optimizes overall performance and reduces wear on the high capacity compressor.

- Continuous Flow Settings 0.5 to 5.0 LPM
- Weight: 36 lbs.
- Dimensions: 21.5 inches H X 12.5 inches W X 13.5 inches D
- Power draw: 285 watts at 2.0 LPM
- Easy-to-use interface
- Wireless connectivity to telehealth app

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR
USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS,
PRECAUTIONS, AND DETAILED SAFETY INFORMATION.
© Copyright 2020 CAIRE Inc. All Rights Reserved
ML-CONC0175 A

NewLife® Elite and Intensity 10 Stationary Oxygen Concentrators

CAIRE'S CONCENTRATORS YOUR SOURCE FOR AT-HOME OXYGEN THERAPY.

The rugged, robust NewLife Intensity stationary oxygen concentrator series is known the world over for its reliability and power to serve high flow patients.

NEWLIFE INTENSITY 10

Ideal for long-term, skilled care facilities where higher oxygen flow is necessary. With 20 psi (138 KPa) output pressure, this unit can safely deliver oxygen to the patient utilizing up to 200 feet of oxygen tubing and can be used in many high pressure applications.

- Continuous Flow Settings 2.0 to 10.0 LPM
- Weight: 58 lbs.
- Output pressure: 20 psi (138 KPa)
- Dimensions: 27.5 inches H X 16.5 inches W X 14.5 inches D
- Power draw: 600 watts



**NewLife Intensity 10
Single Flow
(AS099-101)**

NEWLIFE INTENSITY 10 DUAL FLOW

Two separate flow meters can provide a total of 10 LPM.

NEWLIFE ELITE

The robust construction, ease-of-use, and low maintenance make it the perfect choice for facilities, hospitals, clinical settings, and harsh environmental conditions.

- Continuous Flow Settings 0.5 to 5 LPM
- Weight: 54 lbs. (24.5 kg)
- Output pressure: 20 psi (138 KPa)
- Dimensions: 28.5 inches H X 15.7 inches W X 14.5 inches D
- Power draw: 350 watts



**NewLife Elite
(AS005-4)**

ACCESSORIES



**Humidifier Bottle
(HU003-1)**



**Humidifier Adaptor
(20843882)**



Visit our website at
www.caireinc.com.
Buy a NewLife concentrator today!