


TAKE CHARGE OF YOUR OXYGEN THERAPY.



my
CAIRE™

TELEHEALTH SOLUTION

Access your
information 24/7 to
help ensure your
therapy is on track.






YOUR CARE AT YOUR FINGERTIPS

- Free mobile app for iOS, Android or Surface
- Secure, cloud-based technology
- Collects key data from your oxygen devices
- Access device features and user manual
- Quick set up and easy to use




CONNECTIONS ON YOUR TERMS

- Provide access to family members for additional support
- Request service from provider whenever you need it
- Automated notifications sent in seconds



CRITICAL DATA TO STAY ON TOP OF YOUR THERAPY



NEVER BE OUT OF TOUCH

- > myCAIRE-enabled concentrator passes information to your smart device via wireless connectivity, updating every 11 seconds
- > Smart device transmits the information to your care team via WiFi or mobile data



Connect with the oxygen concentrator that's right for your therapy.

myCAIRE is compatible with our leading portable, transportable and stationary concentrators.

FreeStyle® Comfort®



- > Portable oxygen concentrator
- > Pulse dose flow 1-5
- > Smart O2 delivery
- > Unmatched comfort and convenience
- > Durable design

Eclipse 5®



- > Use as a portable or stationary oxygen concentrator
- > Pulse dose flow 1-9
- > Continuous flow 0.5-3.0 LPM
- > Smart O2 delivery

Companion 5™



- > Stationary oxygen concentrator
- > Continuous flow 0.5-5 LPM
- > Simple user interface
- > Proven performance and reliability

O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy.

Learn more at:

myCAIRE.com



CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.


O₂ is what we do.

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information. myCAIRE is a telehealth solution intended to share specified device data from the enabled CAIRE oxygen concentrator. myCAIRE does not affect oxygen concentrator settings, performance, or oxygen therapy, and is not intended for patient monitoring, diagnosis or treatment. Please seek a doctor's advice before making any medical decisions.

© 2021 CAIRE Inc. All Rights Reserved.

ML-CONC0212 B