





TELEHEALTH SOLUTION



A DIGITAL TRANSFORMATION IN OXYGEN PATIENT MANAGEMENT.

Improve efficiency with 24/7 remote access


CAIRE®
O₂ is what we do.

UPGRADE YOUR O2 OPERATING SYSTEM WITH myCAIRE™

Improve care and streamline your operations with the myCAIRE telehealth solution.

- Rapid, remote data collection from CAIRE oxygen concentrators in the field
- Connect to CAIRE equipment via the patient smartphone or tablet
- Secure, cloud-based technology accessible from desktop or tablet anytime



INCREASE EFFICIENCY WHILE IMPROVING CARE

- Device performance information enables fleet management and troubleshooting
- View of key O2 usage metrics to ensure patient is using equipment as prescribed



IMPROVED PATIENT CONFIDENCE AND SATISFACTION

- Instant access to their own usage data provides ongoing peace of mind
- Increased sense of security that they are always connected to their care team



DIFFERENTIATE FROM YOUR COMPETITORS

- Build credibility with a higher level of patient care and responsiveness
- Empower the patient to take an active role with their therapy




CONVENIENT CONNECTION, 24/7




- myCAIRE-enabled concentrator passes information to the patient's smart device via wireless connectivity, updating every 11 seconds
- Smart device transmits the information to the portal via WiFi or mobile data

EASY ACCESS TO USAGE DATA AND CRITICAL TROUBLESHOOTING


- Multiple provider roles to manage access to data
- Customizable view, filtering and search options for data management
- Comprehensive device management with diagnostic data and GPS location
- Improved ability to troubleshoot remotely
- Visibility to key usage data




EFFICIENT PATIENT MANAGEMENT




AVG FLOW SETTING



USAGE HOURS




O2 STATUS



NOTIFICATIONS/ALARMS

SIMPLE AND INFORMATIVE PATIENT EXPERIENCE

- > Free mobile app for iOS, Android and Surface
- > Quick set up and easy to use
- > View settings and notifications
- > Request service from provider
- > Easy access for patient caregiver





Connect with your patients throughout the continuum of CAIRE

Improve care and streamline your operations with the myCAIRE telehealth solution.



FreeStyle® Comfort®

- > Portable oxygen concentrator
- > Pulse dose flow 1-5
- > Smart O2 delivery
- > Unmatched comfort
- > Durable design



Eclipse 5®

- > Use as a portable or stationary oxygen concentrator
- > Pulse dose flow 1-9
- > Continuous flow 0.5-3.0LPM
- > Smart O2 Delivery



Companion 5™

- > Stationary oxygen concentrator
- > Continuous flow 0.5. to 5 LPM
- > Simple user interface
- > Proven performance and reliability

O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy.

CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.

Learn more at:

myCAIRE.com



SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information. myCAIRE is a telehealth solution intended to share specified device data from the enabled CAIRE oxygen concentrator. myCAIRE does not affect oxygen concentrator settings, performance, or oxygen therapy, and is not intended for patient monitoring, diagnosis or treatment. Please seek a doctor's advice before making any medical decisions.

© 2021 CAIRE Inc. All Rights Reserved.

ML-CONC0213 B

CAIRE®
O₂ is what we do.