

# HOMESTYLE™

Cortina

## Stationary Oxygen Concentrator



## Your O<sub>2</sub> Fleet Solution



Accessible  
Pricing



Affordable  
5-Year  
Warranty  
Option  
Available



Quiet  
Operation



Simple User  
Interface

- Continuous flow 0.5-5 LPM
- Up to 93%±3% oxygen purity
- Sound less than or equal to 40 dB(A)\*
- Power supply: 230V~,50Hz
- Power consumption: 320VA
- 14.5 kg
- 380 mm L x 228 mm W x 592 mm H
- Elegant and compact frame

# Serving your patients across the continuum of CAIRE.

## Non-Delivery Solution

Optimize field efficiency and productivity while improving patient care with the CAIRE non-delivery solution, which combines the industry-leading FreeStyle® Comfort® portable oxygen concentrator with the HOMESTYLE™ Cortina stationary oxygen concentrator.



HOMESTYLE Cortina



FreeStyle Comfort



## The foundation of your patients' oxygen therapy

- A top-tier manufacturer of both oxygen concentrators and liquid oxygen systems to provide the appropriate care for your patients.
- Optimize patient care by pairing the FreeStyle Comfort with the HOMESTYLE Cortina for ambulation and at-home oxygen therapy needs.
- Comprehensive warranty—3-year warranty on the device and affordable 5-year extended warranty option available.<sup>†</sup>
- Effective oxygen delivery—up to 93%±3% concentration from 0.5 to 5.0 LPM, meeting a wide variety of patient prescriptions.
- Simple user interface—easy-to-understand alerts.
- Elegant and compact design—blends into surroundings.

CAIRE®



CAIREinc.com

<sup>\*</sup>Per BEO test report 2025\_183.

<sup>†</sup>See Warranty Statement.

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.  
Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2026 CAIRE Inc. All Rights Reserved.

ML-CONC0291 A