

# Companion® 1000 & 1000T

Liquid Oxygen Portables

## Reliable Oxygen Delivery

The Companion® 1000 Portable Liquid Oxygen Unit and Companion T High-Flow Portable Liquid Oxygen Unit offer a lightweight oxygen therapy option for patients who require continuous oxygen flow. The reliable Companion 1000 unit combines comfort and convenience in a durable, simple-to-use portable liquid oxygen system. It fills quickly and easily from any Companion, Liberator®, or HELiOStm Universal Reservoir. The Companion T unit offers the same design features, plus extended liters per minute settings for higher flow rates.

### Features

- Small, lightweight design allows a more active lifestyle
- Easy portability
- Companion 1000 unit has 11 standard flow settings from 1/4 to 6 LPM
- Companion T unit has built-in control valve with settings from 1/2 to 15 LPM
- Both provide approximately 7.7 hours of use at 2 LPM
- Contents indicator does not require batteries
- Recessed controls and indicators minimize damage and accidental settings changes



  
**CAIRE®**

770.721.7759 | Toll Free 1.800.482.2473

[www.CAIREMedical.com](http://www.CAIREMedical.com)



\*\* SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

**PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR  
USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS,  
PRECAUTIONS, AND DETAILED SAFETY INFORMATION.**

© 2017 CAIRE Inc. All Rights Reserved ML-LOX0011 Rev H

## Product Specifications

Companion 1000		Companion 1000T	
LOX Capacity	3.1 lb (1,41 kg)	LOX Capacity	3.1 lb (1,41 kg)
Gaseous Equivalent Capacity	1058 L	Gaseous Equivalent Capacity	1058 L
Weight, Empty	5.0 lb (2,27 kg)	Weight, Empty	5.6 lb (2,59 kg)
Weight, Filled	8.1 lb (3,67 kg)	Weight, Filled	8.7 lb (3,99 kg)
Height	14.5 in. (368 mm)	Height	14.5 in. (368 mm)
Width	6.5 in (165.1 mm)	Width	6.5 in (165.1 mm)
Depth	6.1 in (154.9 mm)	Depth	6.1 in (154.9 mm)
Typical use time at 2 LPM	7.7 hours	Typical use time at 2 LPM	7.7 hours
Operating Pressure	22 psi (152 kPa)	Operating Pressure	22 psi (152 kPa)
Normal Evaporation Rate	1.5 lb/ day (0,68 kg/ day)	Normal Evaporation Rate	1.5 lb/ day (0,68 kg/ day)
Standard Flow Control Range	Off, .25, .5, .75, 1, 1.5, 2, 2.5, 3, 4, 5, 6 LPM	Standard Flow Control Range	Off, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM



## ORDERING INFORMATION

Description	Catalog Number
Companion 1000 Portable Liquid Oxygen Unit	B-775099-00 / B-700885-00
Companion T High-Flow Liquid Oxygen Unit	B-775532-00 / B-700883-00
Accessories	
LOX Backpack	10018138
Cart, Universal Portable	20765565

  
**CAIRE®**

770.721.7759 | Toll Free 1.800.482.2473

[www.CAIREMedical.com](http://www.CAIREMedical.com)



\*\* SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

**PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR  
USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS,  
PRECAUTIONS, AND DETAILED SAFETY INFORMATION.**

© 2017 CAIRE Inc. All Rights Reserved ML-LOX0011 Rev H