

Troubleshooting Form

for Reliant Oxygen Generator


Customer Name
Company Name
Date
Model Number
Serial Number
Hours

Please fill in the data below to help an AirSep Service team member troubleshoot the oxygen generator. It is important to have a copy of the instruction manual to identify the location of the different components and to ensure the proper operation of the oxygen generator. Contact cpd@caireinc.com for the most up-to-date manual.

The data should not be recorded until the oxygen concentration has stabilized at the generators rated oxygen flow, for at least 15 minutes.

The oxygen generator should be operated in the manual position and at the designed oxygen flowrate of 8 LPM.

Oxygen Receiver Pressure:

Highest Reading (psig)
Lowest Reading (psig)

Oxygen Delivery Pressure:

Highest Reading (psig)
Lowest Reading (psig)

Oxygen Concentration:

Highest Reading (%)
Lowest Reading (%)

Oxygen Flow:

Highest Reading (LPM)
Lowest Reading (LPM)
Comments