

Don't miss out on our upcoming Virtual Service Schools!

CAIRE Service School Training Webinars 2022

CAIRE offers Virtual Service Schools covering both LOX and Concentrator lines.

Each class is a comprehensive program that focuses on the technical and service aspects of the CAIRE family products. Class time is divided between lecture and hands-on workshops. The seminars will help the student develop a better understanding of how liquid oxygen equipment or concentrators work, how to identify the symptoms and causes of potential problems and how to use the technical information that is available in the Technical Manuals. Registration for CAIRE's Service Schools is mandatory. Registration forms must be received two weeks prior to the start of the class in order to guarantee availability and materials.

Table 1 - Details of upcoming training.

Type	Language	Location	Training Dates	Registration by
CONC	English	Ball Ground, GA USA	May 16-17-18	May 2
LOX	English	Ball Ground, GA USA	May 19	May 2
CONC	Italian	Padova, Italy	June 7-8-9	May 26
CONC	English	Ball Ground, GA USA	August 15-16-17	August 1
LOX	English	Ball Ground, GA USA	August 18	August 1
LOX	German	Langenfeld, Germany	September 14	September 1
LOX	LATAM	Toulouse, France	September 20	September 6
CONC	LATAM	Toulouse, France	September 21-22-23	September 6
LOX	Italian	Padova, Italy	October 18	October 4
CONC	German	Langenfeld, Germany	October 25-26-27	October 11
CONC	English	Ball Ground, GA USA	November 7-8-9	October 25
LOX	English	Ball Ground, GA USA	November 10	October 25

*We reserve the right to cancel training seminars at short notice. Customers should ensure that they take necessary precautions to ensure costs are covered in case of cancellation.

** Training will be in the form of a Webinar.

Concentrator Service School Topics include:

- Concentrator hazards and safety precautions
- Principles of pressure, flow, and saturation
- Functions of the major components of a portable/stationary concentrator
- Theory of operation
- Set up and use of test equipment
- Troubleshooting, repair and performance verification/testing procedures



LOX Service School Topics include:

- Liquid oxygen (LOX) hazards and safety precautions
- Principles of pressure, flow and liquid oxygen saturation
- Functions of the major components of a liquid oxygen system
- HELIOS™ system theory of operation
- Reservoir and portable filling procedures
- Set up and use of test equipment

Please submit your registrations forms no later than the dates in Table 1 above. If you are interested in the service schools mentioned above, and it is past the registration date, please reach out to the registration emails below. **Please note that as travel restrictions, and COVID-related guidelines evolve, we intend to restore CAIRE Service School to in person.**

To register for Service School, please email:

Techservice.global@caireinc.com for USA Training or jim.gibson@caireinc.com for European Training.

© Copyright 2022 CAIRE Inc. All Rights Reserved. CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.