



Med Tips



February 2011

SeQual INTEGRATION

As we move forward with our integration of SeQual Technologies into the CAIRE family of products, please note that as of February 1, 2011 we will be handling various SeQual office functions at the CAIRE facility in Georgia. This includes Sales Orders, Returns, Credit, Accounts Receivable, Customer Service, and Technical Support. In accordance, please direct all correspondence as follows:

US Customer Service Contact Information:

Email: customerservice.usa@chart-ind.com
Phone: (800) 482-2473 or 770-721-7759 opt 1
Fax: (888) 932-2473 or 770-721-7758

Europe Customer Service Contact Information:

Email: customerservice.europe@chart-ind.com
Phone: +44 (0) 1344 403100
Fax: +44 (0) 1344 429224

US Technical Support Contact Information:

Email: techservice.usa@chart-ind.com
Phone: (800) 482-2473 or (770) 721-7759 opt 3
Fax: (888) 932-2473 or (770) 721-7758

Europe Technical Support Contact Information:

Email: technicalservice.europe@chart-ind.com
Phone: +44 (0) 1344 403109
Fax: +44 (0) 1344 429224

SEQUAL®



Thank you for your continued support and partnership as we work to fulfill your needs. If you have any questions, please do not hesitate to contact your sales representative.

WIRE HARNESS TEMPORARY COLOR CHANGE FOR PORTABLE UNITS

The color of the feed-through wire (PN 13772093) on all CAIRE portables equipped with an LED contents meter will be temporarily changed from blue to white. The wire is located on the wiring harness (PN 13787877) for the G4 LED meter as shown below.



NEW FLOW SETTINGS & PN's FOR PEDIATRIC FCV's

To better serve the needs of patients requiring low flow rates, the pediatric flow control valve settings for the Liberator and Sprint/Stroller lines have been changed as shown below. Also note that when ordering a pediatric FCV, it is necessary to order the flow rate decal separately. The part numbers to use when ordering are shown below.

Pediatric FCV Flow Settings (LPM)	
Old	New
Off	Off
0.03	0.05
0.06	0.10
0.13	0.20
0.25	0.30
0.38	0.40
0.50	0.50
0.63	0.75
0.75	1.00
1.00	1.25
1.50	1.50
2.00	2.00

Pediatric FCV Part Numbers		
Description	Old Part Number	New Part Number
Liberator Pediatric FCV	11818203T	14700243T
Liberator Pediatric FCV – Decal	13758873	14700438
Sprint/Stroller Portable Pediatric FCV	11923417T	14700235T
Sprint/Stroller Pediatric FCV – Decal	13938639	14700411

NEW LARGER BELT STRAP FOR HELIOS PORTABLE

To better suit the needs of our customers, CAIRE now offers a HELIOS belt extension. This extension attaches to the standard belt, and can be ordered using part number 14951835.

Customer & Technical Info

For ordering information contact

Customer Service:

Europe: +44 1344 403100
Australia/New Zealand: +61 2 9749 4333
US Toll Free: 800.482.2473
Asia/Pac Rim: 770.721.7759
The Americas: 770.721.7759

For technical information contact

Technical Service:

Europe:
Phone: +441344 403100
Email: technicalservice.europe@chart-ind.com
US Toll Free
Phone: 800.482.2473
Email: techservice.usa@chart-ind.com
Asia/Pac Rim
Phone: 770.721.7759
Email: techservice.usa@chart-ind.com

www.cairemedical.com