

## Sprint<sup>™</sup> & Stroller<sup>®</sup>

### Liquid Oxygen Portables

#### Giving you the Freedom to Move

CAIRE, Inc.'s lightweight, sleek, liquid oxygen portables provide ambulatory oxygen therapy patients with a worry-free, and user-friendly system designed to fit easily with the contour of the body.

The Stroller and Sprint are designed to maximize patient convenience and confidence. Patients can easily read the LED meter and quickly perform fill operations themselves. The comfortable easy-grip strap allows for shoulder or hand carrying and is removable for cleaning.

Stroller and Sprint models are available with all popular connectors, to allow for utilization of your existing liquid system inventory.

#### Features

- Ultra-lightweight
  - Sprint: 6 lbs.
  - Stroller: 8 lbs.
- Approx. hours of operation at 2 LPM
  - Sprint: 4 hours
  - Stroller: 8 hours
- LED meter easy to read
- Side Fill and Top Fill configurations
- Increased ease of mobility



**Call 1.800.482.2473 to order!**  
[www.CAIREmedical.com](http://www.CAIREmedical.com)

# Sprint™ & Stroller®

## Liquid Oxygen Portables



### Product Specifications

	<b>Sprint™</b>	<b>Stroller®</b>
LOX capacity	1.5 lb (0,68 kg)	3.0 lb (1,36 kg)
Gaseous equivalent capacity	513 L	1,026 L
Weight, filled	6 lb (2,72 kg)	8 lb (3,63 kg)
Weight, empty	4.5 lb (2,04 kg)	5 lb (2,27 kg)
Height	11.75" (298 mm)	13.5" (343 mm)
Typical use time at 2 LPM	4.3 hours	8 hours
Operating pressure	20psi (137 kPA)	20psi (137 kPA)
Normal evaporation rate	1.3 lb/day (0,57 kg/day)	1.3 lb/day (0,57 kg/day)
Standard flow control range	Off, .25, .5, .75, 1, 1.5, 2, 2.5, 3, 4, 5, 6	Off, .25, .5, .75, 1, 1.5, 2, 2.5, 3, 4, 5, 6

### Ordering Information

<b>Description</b>	<b>Catalog Number</b>
Sprint SF (side fill)	10564117
Sprint TF (top fill)	10564096
Stroller SF (side fill) LED	10564125
Stroller TF (top fill) LED	10564109
Stroller TF (top fill) LED soft case	14768264
Stroller SF (side fill) Scale	10998054
Stroller TF (top fill) Scale	10997991
Stroller TF (top fill) Scale soft case	14765549

\* Catalog Numbers reflect US English Part Numbers, for language-specific part numbers please contact Customer Service

CAIRE, Inc.  
2200 Airport Industrial Drive, Suite 500  
Ball Ground, GA 30107  
Ph 770-721-7759 • Toll-Free 1-800-482-2473  
Fax 888-932-2473

Chart BioMedical, Ltd.  
Unit 6, Ashville Way  
Wokingham RG41 2PL, United Kingdom  
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245