

PN 21500320 A

RELEASE DATE: January 13, 2025

MODEL: CAIRE Liquid Oxygen Models: HELiOS™ H300 and H850

ISSUE: Dual Lumen Cannulas
Part Numbers: 6-778058-00, 6-778057-00, 10035468, 10035467

NOTES: Effective January 13, 2025, CAIRE will be discontinuing some of our offerings of the Dual Lumen Cannulas as they are no longer available from the supplier.

Please refer to the table below for details on the part numbers being discontinued and the recommended replacements or alternatives. Those cannulas are certified and compliant to the EU MDR directives.

Obsolete Part Number	Replaced by
6-778058-00 DUAL LUMEN CANNULA, 5'	20574985 DEMAND NASAL CANNULA 4'
6-778057-00 DUAL LUMEN CANNULA 7'	20574978 DEMAND NASAL CANNULA 7'
10035468 OXYGEN DUAL LUMEN CANNULA, 1.5M	20574985 DEMAND NASAL CANNULA 4'
10035467 OXYGEN DUAL LUMEN CANNULA, 2.1M	20574978 DEMAND NASAL CANNULA 7'

CAIRE will continue to provide and stock the alternative part numbers listed, ensuring the continuity of our product accessory offerings and solutions.

Please contact your local Sales Representative or our Customer Service, or Technical Service teams for assistance with any inquiries or concerns.

CONTACT: For technical questions, contact Technical Service:

North America: 800-482-2473 or +1-770-721-7700

Europe: +44 (0) 118 936 7060

Japan: 050-5050-1275

techservice.global@caireinc.com

For ordering information or general inquiries, contact Customer Service:

North America: 800-482-2473 or +1-770-721-7700

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com

Japan: 050-5050-1275

customerservice.japan@caireinc.com

PN 21500493 A

RELEASE DATE: January 15, 2025

MODEL: CAIRE Liquid Oxygen Models: Stroller[®] and Sprint[™]

ISSUE: Warming Coil Tubing Change

NOTES: Effective January 15, 2025, the convoluted tubing for the vaporizer on the Stroller and Sprint units will have a different form due to a change on the manufacturing process. The new coils are backwards compatible with the current design. The part number for the new design coil is TU303-1. Please see below pictures.



Current Design 10650277



New Design TU303-1

Please contact your local Sales Representative or our Customer Service, or Technical Service teams for assistance with any inquiries or concerns.

CONTACT: For technical questions, contact Technical Service:

North America: 800-482-2473 or +1-770-721-7700

Europe: +44 (0) 118 936 7060

Japan: 050-5050-1275

techservice.global@caireinc.com

For ordering information or general inquiries, contact Customer Service:

North America: 800-482-2473 or +1-770-721-7700

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com

Japan: 050-5050-1275

customerservice.japan@caireinc.com