

**PN 21500897 A**

**RELEASE DATE:** June 12, 2025

**MODEL:** FreeStyle® Comfort® Single Battery

**SUBJECT:** Part Number Change BT034-1 to BT034-2

**DETAILS:** Effective June 11, 2025, the part number for the BT034-1, FreeStyle Comfort 8 cell battery, will change to BT034-2. This change will not affect the form, fit, or function of the battery.

**NOTE:** Customers ordering spare batteries will now need to order part number BT034-2S instead of the previous BT034-1S.

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquires, contact Customer Service:

USA: 800-482-2473

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

**PN 21500717 A**

**RELEASE DATE:** April 22, 2025

**MODEL:** Eclipse 5® Oxygen Concentrator

**SUBJECT:** Power Cartridge Battery Capacity Change

**DETAILS:** CAIRE announces the Federal Aviation Administration (FAA) recently updated its website to include new safety content - PackSafe – Lithium Batteries. This information addresses carrying and use of chargeable and non-rechargeable lithium batteries on commercial air flights. Additional details can be found on the FAA website: <https://www.faa.gov/hazmat/packsafe/lithium-batteries>

The Eclipse® 5 portable oxygen concentrator meets these FAA Guidelines as its power cartridge is currently composed of two (2) isolated 98.9-watt-hour lithium-ion battery packs, each fall under the FAA limit of 160 Wh per battery. The Eclipse 5 does not contain lithium metal batteries.

This change does not affect the charging and usage times of the Eclipse power cartridge. CAIRE recommends that anyone planning to travel with their Eclipse contact their airline directly to see whether use of this device on board is permissible before their trip to ensure a smooth travel experience. Please refer to the *Travel and Safety* section within the CAIRE website, as well as the Eclipse power cartridge label for more information.

### Current Version of Power Cartridge Label



Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Ref 21500717 Rev A

**CONTACT:** For technical questions or concerns, contact Technical Service:

North America: 800-482-2473 or +1-770-721-7700

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquires, contact Customer Service:

North America: 800-482-2473 or +1-770-721-7700

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

**PN 21500973 A**

**RELEASE DATE:** July 21, 2025

**MODEL:** FreeStyle® Comfort®

**SUBJECT:** FreeStyle Comfort Outlet Barb Update

**DETAILS:** To ensure the FreeStyle Comfort portable oxygen concentrator complies with updated requirements to ISO 80601, CAIRE is announcing a change to the oxygen outlet barb, part number F0971-1, to new part number F0992-1 effective 6/12/2025. As part of this update, the F0992-1 barb will have maximum PSI and LPM ratings stamped into the fascia of the outlet barb. See figure 1 below for differences between the previous vs. new part numbers.

**Note: As with the F0971-1, the new F0992-1 barb will only fit the FreeStyle Comfort 2.0 units.**



**F0971-1**



**F0992-1**

**Figure 1: Previous vs. Current Outlet Barb**

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Ref 21500973 Rev A

[www.caireinc.com](http://www.caireinc.com)

Page 1 of 2

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

**PN 21500995 A**

**RELEASE DATE:** July 23, 2025

**MODEL:** VisionAire™ 5 Concentrator

**SUBJECT:** VisionAire 5 Outlet Barb Update

**DETAILS:** To ensure the Visionaire 5 oxygen concentrator complies with updated requirements to ISO 80601, CAIRE is announcing a change to the oxygen outlet barb, part number F0671-2, to new part number 21498726 beginning 7/14/2025. As part of this update, the 21498726 barb will have maximum PSI and LPM ratings stamped into the fascia of the outlet barb. See figure 1 below for differences between the previous vs. new part numbers.



**F0671-2**



**21498726**

**Figure 1: Previous vs. Current Outlet Barb**

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

## PN 21501001 A

**RELEASE DATE:** July 28, 2025

**MODEL:** Eclipse 5™

**SUBJECT:** Eclipse 5 Outlet Barb Update

**DETAILS:** To ensure the Eclipse 5 transportable oxygen concentrator complies with updated requirements to ISO 80601, CAIRE is announcing a change to the oxygen outlet barb, part number 3922-SEQ, to new part number 21498724 beginning 7/09/2025. As part of this update, the 21498724 barb will have maximum PSI and LPM ratings stamped into the fascia of the outlet barb. See figure 1 below for differences between the previous vs. new part numbers.



**3922-SEQ**



**21498724**

**Figure 1: Previous vs. Current Outlet Barb**

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

## PN 21501074 A

**RELEASE DATE:** October 10, 2025

**MODEL:** CAIRE FreeStyle® Comfort® (3-Button version *without* autoSAT®)

**ISSUE:** Discontinuation of FreeStyle Comfort 3-Button Versions

**NOTES:** As part of our continued commitment to innovation, the production and sale of CAIRE’s FreeStyle Comfort portable oxygen concentrator (3-button version *without* autoSAT) will cease **effective December 15, 2025**.

CAIRE will honor the original warranty on all units from the date of purchase, as well as provide support for the product’s service life (5 years) from the date of purchase.

For information about ordering the FreeStyle Comfort version *with* autoSAT (4-Button), please contact your CAIRE sales representative or customer service.

The table below provides details on the part numbers being discontinued.

Discontinued PN	Description	2 <sup>ND</sup> Description
AS200-1	FNL FREESTYLE COMFORT TYPE A	3B, ENG, SINGLE BTRY
AS200-101	FNL FREESTYLE COMFORT TYPE A	3B, ENG, DOUBLE BTRY
AS200-101DTC	FNL FREESTYLE COMFORT TYPE A	3B, ENG, DOUBLE BTRY
AS200-108	FNL FREESTYLE COMFORT TYPE G	ENGLISH, DOUBLE BTRY
AS200-109	FNL FREESTYLE COMFORT	3B,CHINA, DOUBLE BTRY
AS200-112	FNL FREESTYLE COMFORT	3B,ARABIC, DOUBLE BTRY
AS200-12	FNL FREESTYLE COMFORT	3B,ARABIC, SINGLE BTRY
AS200-126	FNL FREESTYLE COMFORT	3B,TAIWAN,DOUBLE BTRY
AS200-1DTC	FNL FREESTYLE COMFORT TYPE A	3B,ENG,SINGLE BTRY
AS200-201	FNL FREESTYLE COMFORT TYPE A	3B, ENG,SINGLE BTRY

Copyright © 2025 CAIRE Inc.  
 CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

AS200-209	FNL FREESTYLE COMFORT	3B,CHINA,2 SINGLE BTRY
AS200-212	FNL FREESTYLE COMFORT	3B,ARABIC,2 SINGLE BTRY
AS200-226	FNL FREESTYLE COMFORT	3B, TAIWAN, 2 SINGLE BTRY
AS200-228	FNL FREESTYLE COMFORT	3B, BRAZIL, 2 SINGLE BATTERY
AS200-26	FNL FREESTYLE COMFORT	3B,TAIWAN,SINGLE BTRY
AS200-8	FNL FREESTYLE COMFORT TYPE G	3B, ENG, SINGLE BTRY
AS200-9	FNL FREESTYLE COMFORT	3B,CHINA,SINGLE BTRY
AS200-1AA	FNL FREESTYLE COMFORT	3B, TAIWAN, ALL ACCESSORY BAG
AS200-12AA	FNL FREESTYLE COMFORT	3B, TAIWAN, ALL ACCESSORY BAG, 2 SINGLE BATTERY
AS200-26AA	FNL FREESTYLE COMFORT	3B, TAIWAN, ALL ACCESSORY, SINGLE BATTERY

**CONTACT:** For technical questions or concerns, contact Technical Service:

North America: 800-482-2473 or +1-770-721-7700

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquires, contact Customer Service:

North America: 800-482-2473 or +1-770-721-7700

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

**PN 21500493 B**


**RELEASE DATE:** January 15, 2025

**REVISION DATE:** April 25, 2025 (Rev B)

**MODEL:** CAIRE Liquid Oxygen Models: HI Flow Stroller<sup>®</sup>, Stroller<sup>®</sup> and Sprint<sup>™</sup>

**ISSUE:** Warming Coil Tubing Change

**NOTES:** Effective January 15, 2025, the convoluted tubing for the vaporizer coil on the Hi Flow, Stroller and Sprint units will have a different form due to a change in the manufacturing process. The new coils are backwards compatible with the current design and do not affect the Vaporizer Coil Kit PN 11033086 for Sprint, PN 11891064 for Stroller and PN 12934695 for HI Flow Stroller that includes the coil, the fitting nuts and the ferrules. The part number for the new design coil included in the above kits is TU303-1 for Sprint and Stroller and TU303-2 for HI Flow Stroller. Please see below pictures.




Current Design 10650277



New Design TU303-1



Current Design 11938918



New Design TU303-2

Note: When running a Hi Flow Stroller at 15 LPM with the newer tube fitted, there may be a slight whistling noise heard for the first 30-60 seconds. This is normal and does not affect the function of the equipment.

Please contact your local Sales Representative or our Customer Service, or Technical Service teams for assistance with any inquiries or concerns.

**CONTACT:** For technical questions, contact Technical Service:

North America: 800-482-2473 or +1-770-721-7700

Europe: +44 (0) 118 936 7060

Japan: 050-5050-1275

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquiries, contact Customer Service:

North America: 800-482-2473 or +1-770-721-7700

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)

Japan: 050-5050-1275

[customerservice.japan@caireinc.com](mailto:customerservice.japan@caireinc.com)

**PN 21501076 A**

**RELEASE DATE:** September 11, 2025

**MODEL:** FreeStyle® Comfort®

**SUBJECT:** Unit Turning Off When Transitioning Power Source

**DETAILS:** CAIRE is reminding our customers if they experience the FreeStyle Comfort unit turning off when transitioning from battery power to AC or DC power, or vice versa, to remove all power sources from the unit prior to reconnecting.

Disconnect AC or DC power source from the unit and remove the battery. Install the battery into unit, connect to AC or DC power source and turn unit on.

Please refer to our maintenance manuals (MN233 for 3-Button Models or MN242 for 4-Button Models) for additional product information.

**CONTACT:** For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

[techservice.global@caireinc.com](mailto:techservice.global@caireinc.com)

For ordering information or general inquires, contact Customer Service:

USA: 800-482-2473

[customerservice.usa@caireinc.com](mailto:customerservice.usa@caireinc.com)

Europe: +44 (0) 118 936 7060

[customerservice.europe@caireinc.com](mailto:customerservice.europe@caireinc.com)