


PN 21500995 A

RELEASE DATE: July 23, 2025

MODEL: VisionAire™ 5 Concentrator

SUBJECT: VisionAire 5 Outlet Barb Update

DETAILS: To ensure the Visionaire 5 oxygen concentrator complies with updated requirements to ISO 80601, CAIRE is announcing a change to the oxygen outlet barb, part number F0671-2, to new part number 21498726 beginning 7/14/2025. As part of this update, the 21498726 barb will have maximum PSI and LPM ratings stamped into the fascia of the outlet barb. See figure 1 below for differences between the previous vs. new part numbers.



F0671-2



21498726

Figure 1: Previous vs. Current Outlet Barb

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Ref 21500995 Rev A

www.caireinc.com

Page 1 of 2

CONTACT: For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

techservice.global@caireinc.com

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com