

PN 21501001 A

RELEASE DATE: July 28, 2025

MODEL: Eclipse 5™

SUBJECT: Eclipse 5 Outlet Barb Update

DETAILS: To ensure the Eclipse 5 transportable oxygen concentrator complies with updated requirements to ISO 80601, CAIRE is announcing a change to the oxygen outlet barb, part number 3922-SEQ, to new part number 21498724 beginning 7/09/2025. As part of this update, the 21498724 barb will have maximum PSI and LPM ratings stamped into the fascia of the outlet barb. See figure 1 below for differences between the previous vs. new part numbers.



3922-SEQ



21498724

Figure 1: Previous vs. Current Outlet Barb

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Ref 21501001 Rev A

www.caireinc.com

Page 1 of 2

CONTACT: For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

techservice.global@caireinc.com

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com