


PN 21501123 A

RELEASE DATE: October 1, 2025

MODEL: AirSep® Topaz™ Series Oxygen Generators (All Variants)

SUBJECT: Tubing Material Change Notification

DETAILS: This bulletin is to inform our customers of a material change to the tubing used in the Topaz Series Oxygen Generators, specifically the (1') section connecting the heat exchanger to the valve block as shown below. This change was implemented in May 2025 and provides improved durability while maintaining the same fit, form, and function as the previous version.



Previous vs. Current Tubing

**Note: Units manufactured with previous version tubing do not require to be updated to the new version.**

CONTACT: For technical questions, ordering information or general inquiries contact:

Within the US: 1-800-874-0202

Outside the US: +1-716-691-0202

[cpd@caireinc.com](mailto:cpd@caireinc.com)

Copyright © 2025 AirSep Corporation

AirSep Corporation reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.