

PN 21500650 A

RELEASE DATE: Revision A, 11 March 2025

MODEL: VisionAire™ 5 Concentrator

ISSUE: VisionAire Muffler Tube.

NOTES: CAIRE® would like to notify our valued customers that commencing with serial number GPB0125010031 we have fitted a tie wrap to the compressor enclosure that secures the muffler tube in position. Figure 1.



Figure 1.

The purpose of this tie wrap is to hold the tubing in position and minimize friction with the compressor braided tube. Figure 2.



Figure 2.

Copyright © 2025 CAIRE Inc.

CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Ref 21500650 Rev A

www.caireinc.com

Page 1 of 2

Note: This is not a compulsory update for previously manufactured units.

For further information please contact the Technical Service Team.

CONTACT:

For technical questions or concerns, contact Technical Service:

USA: 800-482-2473

Europe: +44 (0) 118 936 7060

techservice.global@caireinc.com

For ordering information or general inquiries, contact Customer Service:

USA: 800-482-2473

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com