

PN 21500717 A

RELEASE DATE: April 22, 2025

MODEL: Eclipse 5® Oxygen Concentrator

SUBJECT: Power Cartridge Battery Capacity Change

DETAILS: CAIRE announces the Federal Aviation Administration (FAA) recently updated its website to include new safety content - PackSafe – Lithium Batteries. This information addresses carrying and use of chargeable and non-rechargeable lithium batteries on commercial air flights. Additional details can be found on the FAA website: <https://www.faa.gov/hazmat/packsafe/lithium-batteries>

The Eclipse® 5 portable oxygen concentrator meets these FAA Guidelines as its power cartridge is currently composed of two (2) isolated 98.9-watt-hour lithium-ion battery packs, each fall under the FAA limit of 160 Wh per battery. The Eclipse 5 does not contain lithium metal batteries.

This change does not affect the charging and usage times of the Eclipse power cartridge. CAIRE recommends that anyone planning to travel with their Eclipse contact their airline directly to see whether use of this device on board is permissible before their trip to ensure a smooth travel experience. Please refer to the *Travel and Safety* section within the CAIRE website, as well as the Eclipse power cartridge label for more information.

Current Version of Power Cartridge Label



Service Bulletin

CAIRE®

CONTACT: For technical questions or concerns, contact Technical Service:

North America: 800-482-2473 or +1-770-721-7700

Europe: +44 (0) 118 936 7060

techservice.global@caireinc.com

For ordering information or general inquiries, contact Customer Service:

North America: 800-482-2473 or +1-770-721-7700

customerservice.usa@caireinc.com

Europe: +44 (0) 118 936 7060

customerservice.europe@caireinc.com